



SIGVARIS
GROUP

CATALOGUE
2 0 2 5

Medical Compression Therapy Solutions



sigvaris

mobilis
by SIGVARIS GROUP

CUSTOMER ID



SIGVARIS
GROUP

WE'VE GOT YOU COVERED
FOR **ALL YOUR COMPRESSION**
THERAPY NEEDS!

Our mission

Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our mission for helping people feel their best is why we've been a global market leader for more than 160 years.

Our compression portfolio includes Medical, Consumer Health Care, products for Medically Complex Edema, and Soft Orthopedic Goods that are specifically designed to restore, assist and enhance people's lives. We are passionate about making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern customers' lives.

When heritage, competence, and dedication come together, there's no limit to where we'll go.

Live your life. Be yourself.

For more information, please visit:
sigvaris.com

We produce the best compression garments to enhance our customers' health and quality of life.

We are dedicated to helping legs look and feel their best.



**RECIPROCATED
TOE PRODUCTS**

Crafted with a smooth, flat toe seam and contoured toe box for added comfort.



**THE MOST COMFORTABLE
FABRIC OPTIONS AVAILABLE**

Including the finest cottons, wools, microfibers and double covered yarns.



**WE ARE DEDICATED
TO INNOVATION**

We continue to invest in medical research to help people feel their best. Every day.



**WE OFFER
MORE VARIETY**

More materials, colours, patterns, and styles than any other medical compression brand.



**THE LARGEST SELECTION OF
PREMIUM NATURAL FIBERS**

Including Sea Island Cotton, Supima® Cotton, Linen and Merino Wool.



**THE HIGHEST
QUALITY STANDARDS**

More than 20 people handle our products from start to finish, ensuring every garment is of the highest quality.



**DESIGNED FOR
YOUR BEST FIT**

Using studies on human data measurements, we fit 99.49% of the population.

**WHY
SIGVARIS
GROUP?**

Customer Care

DEALER PORTAL

The Sigvaris Group Dealer Portal is a secure online platform for easily ordering products and managing your account.

24/7 access – order anytime!

View order history
and tracking

Sizing Wizard

Obtain invoice copies

UPC code for a product

Shop the catalogue

Quick order entry

Want to activate your account?
Email: canadaorders@sigvaris.com



ONLINE CHAT

Found on dealer portal

Real time convenience with a live agent
Instant responses
Fast customer inquiry resolution

To Order:

Online Dealer Portal: Sigvaris-online.com

Phone: **1.800.363.4999**

Email: canadaorders@sigvaris.com



Join our Certified Fitter Course

**Elevate your knowledge by enrolling
in our Sigvaris Online Seminar!**

Dive into the world of compression wear and become
a certified expert in the art of comfort and wellness.

Upon completion you will receive:

- ✓ **Certificate**
- ✓ **Fitter pin**
- ✓ **Complimentary pair
of compression socks**

JOIN HERE



		PRODUCT	STYLES	COMPRESSION LEVEL	PAGE
MEDICAL COMPRESSION THERAPY SOLUTIONS	ESSENTIAL	COTTON 230	● calf, calf w/grip-top, thigh w/grip-top ○ calf, thigh w/grip-top	20–30 mmHg 30–40 mmHg	8
		THERMOREGULATING COTTON ETH	● calf	CCL1	10
		OPAQUE 860 (Petite & plus sizes available)	● calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose ○ calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh w/waist attachment	20–30 mmHg 30–40 mmHg	12, 14
	STYLE	SEA ISLAND COTTON 220	● calf	20–30 mmHg	18
		MERINO WOOL 240	● calf	20–30 mmHg	20
		LINEN 250	● calf	20–30 mmHg	22
		MEDIUM SHEER 750 (for women only)	● calf, thigh w/grip-top, pantyhose	20–30 mmHg	24
		SHEER 780 (for women only)	● calf, thigh w/grip-top, pantyhose, ○ calf, thigh w/grip-top, pantyhose	15–20 mmHg 20–30 mmHg 30–40 mmHg	26
		MICROFIBER 820 (for men only)	● calf, calf w/grip-top, thigh w/grip-top (extended foot sizes available)	15–20 mmHg 20–30 mmHg 30–40 mmHg	28
		MICROFIBER PATTERNS 830	● calf	20–30 mmHg	30
		SOFT OPAQUE 840 (for women only)	● calf, calf w/grip-top, thigh w/grip-top, pantyhose, maternity pantyhose ○ calf, thigh w/grip-top, pantyhose	15–20 mmHg 20–30 mmHg 30–40 mmHg	32
	MOTION	CUSHIONED COTTON 360	● calf	20–30 mmHg	36
		HIGH TECH 412C	● calf	20–30 mmHg	38
		THERMOREGULATING WOOL 422	● calf	20–30 mmHg	40
		FLOW TECH® 432	● calf	20–30 mmHg	42
	SPECIALTY	NATURAL RUBBER 500	○ calf, thigh, thigh w/grip-top, thigh w/waist attachment, pantyhose	30–40 mmHg 40–50 mmHg 50–60 mmHg	46
		ADVANCE ARMSLEEVE 910	○ armsleeve w/grip-top w/gauntlet and without gauntlet	15–20 mmHg 20–30 mmHg	48
INELASTIC WRAPS AND MEDICALLY COMPLEX EDEMA WEAR				50	

Closed Toe = ● Open Toe = ○

	PRODUCT	STYLES	COMPRESSION LEVEL	PAGE
CONSUMER HEALTH	SHEER FASHION 120	● calf, thigh w/grip-top, pantyhose ● calf	15–20 mmHg	66
	SOFT SILHOUETTE LEGGINGS 170	● open foot	15–20 mmHg	67
	CASUAL COTTON 146/186	● calf	15–20 mmHg	68
	TRAVENO TRAVEL SOCKS 146/T - 186/T	● calf	15–20 mmHg	68
	SEA ISLAND COTTON 151/191	● calf	15–20 mmHg	68
	ALL-SEASON MERINO WOOL 152/192	● calf	15–20 mmHg	69
	LINEN 153/193	● calf	15–20 mmHg	69
	MICROFIBER SHADES 143/183	● calf	15–20 mmHg	69
	ATHLETIC RECOVERY SOCKS 401	● calf	15–20 mmHg	70
	PERFORMANCE SLEEVES 412V	● open foot	20–30 mmHg	70
	MERINO OUTDOOR SOCKS 421	● calf	15–20 mmHg	71
	EVERSOFT DIABETIC SOCKS 160	● calf	8–15 mmHg	71
	STANDARD & UNISEX SIZING CHARTS			72, 73
	ACCESSORIES			74
RESOURCES	SYMPTOMS & INDICATIONS			82
	GLOSSARIES OF ICONS & TERMS			83
	HOW TO MEASURE & ORDER			85
MOBILIS	ORTHOPEDIC BRACING			87
	ELBOW			92
	WRIST & THUMB			93
	BACK			99
	KNEE			102
	ANKLE			106
	HOW TO ORDER MOBILIS			109

sigvaris ESSENTIAL

Medical Compression Therapy Solutions



Support Yourself

Why choose **Essential**?

- Reduce swelling and discomfort
- Long-distance travel
- Everyday comfort and reliability

ESSENTIAL PRODUCTS

Series	Product	Page
230	Cotton	8
ETH	Thermoregulating Cotton	10
860	Opaque (Women)	12
860	Opaque (Men)	14

BEST SUITED FOR

- Containment of various types of edema, including venous edema, phlebolymphe~~de~~ma, and early stages of lymphedema
- Ideal for dress or casual shoes
- Suited for post-vein treatment

Featured
colour is
Light Beige

COTTON 230

for women & men

Compression levels

20–30 mmHg

30–40 mmHg

Features and Benefits



Sensitive Skin



Cotton



Breathable

- 100% Supima Cotton against the skin
- Ideal for sensitive skin and allergies
- Breathable and durable
- Excellent containment of various types of edema
- Classic ribbed design for everyday wear

Fiber Content

- 50% Nylon, 25% Supima Cotton, 25% Spandex (content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Calf w/ Grip-top
Medium Long in Black

232CMLM99/S



Silicone dotted grip-top available
for calf and thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		NA	
Leg Length: Women & Open Toe	Short	Long	Short	Long	Short	Long		
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")		
Leg Length: Men	Short	Long	Short	Long	Short	Long		
	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")		
Size	SS	SL	MS	ML	LS	LL	XS	XL






Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

Open Toe = ○

				
WHITE 00	NAVY 10	BLACK MIST 14	LIGHT BEIGE 66	BLACK 99

Calf (C)*	●	●	●	●●	●
Calf w/Grip-Top (C/S) Silicone dotted grip-top band				●**	●
Thigh w/Grip-Top (N) Silicone dotted grip-top band				●●	●

* 30–40 mmHg extra large calf only available in black with grip-top

** 20–30 mmHg only

Product Code Example: 232NMLW66

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
23	2	N	ML	W	66

Product Code Example: 232CLSM66/S

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour	Grip-top
23	2	C	LS	M	66	/S

Featured
colour is
Nature

THERMOREGULATING COTTON ETH

for women & men

Compression Class

CCL1

Features and Benefits



Sensitive Skin



Thermoregulating



Opaque

- Fine cotton layer lies directly on the skin and has a thermoregulating effect
- Cotton absorbs excess moisture directly from the skin, stores it inside the fibers and quickly releases it through evaporation
- Ideal for people with sensitive skin
- High quality, reinforced heel
- Opaque, non-transparent look and feel

Fiber Content

- 52% Polyamide, 34% Elastane, 14% Cotton



CCL1
Essential Thermoregulating; CCL1;
Calf; Size: Medium; Length: Long;
Closed Toe Black

ETH1D23199





Quadrinova knobbed grip top
Square nubs feel gentle on the
skin. Features a soft, protective top
seam for comfort.



See product images online
sigvaris.com

SIZING CHART

	Circumference cm							
	XS	XS Plus	S	S Plus	M	M Plus	L	L Plus
cD	27-34	31-38	31-38	34-41	34-41	37-44	37-44	41-48
cC	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49
cB ₁	24-29	24-29	26-32	26-32	29-35	29-35	31-37	31-37
cB	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29
cY	28-33	28-33	29-35	29-35	31-37	31-37	32-38	32-38
cA	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29



Length cm

normal	35 – 40
long	41 – 45

How to Measure:

cB (Ankle): Circumference at the narrowest point just above the ankle bone. This is the reference point for the highest compression and is critical for sizing.

cD (Below Knee): Circumference just below the knee, at the transition from calf to knee.

ℓD (Length): From the floor to just below the knee crease.

Additional measurement points are only necessary for custom made circular knit products

STYLES AND COLOURS



Available in calf with or without **griptop**.
Calf and closed toe only.

HOW TO ORDER

Example Product: Essential Thermoregulating; CCL1; Calf; Size: Medium; Length: Long; Closed Toe Black

Code: ETH1D23199

Segment and Series Name	Compression Level	Style	Plus Size	Ankle Size	Length	Grip Top	Foot	Colour
E = Essential TH = Thermoregulating	CCL1 = 1	D = Calf	If Yes = P If No = blank	0 = x-small 1 = small 2 = medium 3 = large 4 = x-large	2 = Normal 3 = Long	3 = With Grip Top If No = blank	Closed = 1	Black = 99 Nature = 76
ETH	1	D		2	3		1	99



Featured
colour is
Honey

OPAQUE 860

for women

Compression Levels

20–30 mmHg

30–40 mmHg

Features and Benefits



Durable



Largest Size Range



Opaque

- Largest size range, from petite to plus
- Durable construction combined with a smooth opaque look
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema
- Plus size pantyhose are an alternative to maternity pantyhose

Fiber Content

- 64% Nylon, 36% Spandex
(content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Calf Medium Long
in Black

862CMLW99



Silicone dotted grip-top available
for calf and thigh



See product images online
sigvaris.com

SIZING CHART

	Small			Medium			Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Circumference: Calf Plus w/Grip-Top	NA			NA	39.5–53 cm (15.5"–21")		45.5–61 cm (18"–24")		51–71 cm (20"–28")	
Calf Length	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤35.5 cm (14")	≤41 cm (16")	>41 cm (16")	≤35.5 cm (14")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")			48–73 cm (19"–28.5")			52–84 cm (20.5"–33")		64–95 cm (25"–37.5")	
Leg Length	Petite	Short	Long	Petite	Short	Long	Short	Long	Short	Long
	≤65 cm (25.5")	≤72 cm (28.5")	>72 cm (28.5")	≤65 cm (25.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")			87–127 cm (34.5"–50")			98–150.5 cm (38.5"–59.5")		NA	
Plus or Maternity Size Hip Circumference	NA	≤142 cm (56")		NA	≤157.5 cm (62")		≤178 cm (70")		NA	
Size	SP	SS	SL	MP	MS	ML	LS	LL	XS	XL
Calf Plus w/Grip-Top					MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Women's Shoe Size : Petite Leg 4–8, Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

HONEY 33	GOLDEN 36	LIGHT BEIGE 66	BLACK 99
-------------	--------------	----------------------	-------------

Calf (C)	●	●	○	●●
Calf Petite Sizes (C)* Not available with dotted grip-top band	●	●		●
Calf w/Grip-Top (C/S) Silicone dotted grip-top band			●●	●
Calf Plus w/Grip-Top (C/P) Only available with silicone dotted grip-top band			●●	●●*
Thigh w/Grip-Top (N) Silicone dotted grip-top band	●	●	○	●●
Thigh w/Grip-Top (N) XS & XL Silicone dotted grip-top band			●●	○
Thigh Petite Sizes (N)*	●	●		●
Pantyhose (P)**	●	●		●
Pantyhose Petite (P)*	●	●		●
Pantyhose Plus and Maternity (P/P)* **			●	●
Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R)			○	

* 20-30 mmHg only

** Not available in XS or XL sizing

Product Code Example: 863NMSW33

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour
86	3	N	MS	W	33

Product Code Example: 862CLLO66/P

Series	mmHg Level	Style	Size	Gender/ Open Toe	Colour	Plus size with Grip-top
86	2	C	LL	O	66	/P

Featured
colour is
Ribbed Black



OPAQUE 860

for men

Compression Levels


20–30 mmHg

30–40 mmHg

Features and Benefits



Durable



Largest Size Range



Opaque

- Durable construction with a ribbed design on closed toe calf products
- Smooth opaque design on open-toe styles
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema


Fiber Content

- 64% Nylon, 36% Spandex
(content can vary slightly by style)
- Latex free



20–30 mmHg
Open Toe Thigh Medium Long
in Smooth Light Beige

862NMLO66



Silicone dotted grip-top available
for calf and thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Circumference: Calf Plus w/Grip-Top	NA		39.5–53 cm (15.5"–21")		45.5–61 cm (18"–24")		51–71 cm (20"–28")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		64–95 cm (25"–37.5")	
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Size	SS	SL	MS	ML	LS	LL	XS	XL
Calf Plus w/Grip-Top			MS/P	ML/P	LS/P	LL/P	XS/P	XL/P

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	RIBBED LIGHT BEIGE 66	RIBBED BLACK 99	SMOOTH LIGHT BEIGE 66	SMOOTH BLACK 99
Calf (C)	●*	●*	○	○
Calf w/Grip-Top (C/S) Silicone dotted grip-top band			○	
Calf Plus w/Grip-Top (C/P) Only available with silicone dotted grip-top band			○	○*
Open Toe Thigh w/Grip-Top (N) Silicone dotted grip-top band			○	○
Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R)			○	

* 20–30 mmHg only

** Not available in XS or XL with waist attachment

Product Code Example: 862CMLM99

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
86	2	C	ML	M	99

Product Code Example: 862NMLO66

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
86	2	N	ML	O	66

Medical Compression Therapy Solutions



Featured colours
from left to right:
Toasted Almond,
Warm Sand,
Golden, Café.

Express Yourself

Why choose **Style**?

- Fashionable, modern patterns and designs
- Special occasion and formal wear
- Our most luxurious feel

STYLE PRODUCTS

Series	Product	Page
220	Sea Island Cotton	18
240	Merino Wool	20
250	Linen	22
750	Medium Sheer	24
780	Sheer	26
820	Microfiber	28
830	Microfiber Patterns	30
840	Soft Opaque	32

BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- Pre- and post-vein treatment with minimal edema
- Pregnancy
- People who prefer fashionable options
- Ideal for business or dress attire



Featured colour
is Navy

SEA ISLAND COTTON 220

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 68

Features and Benefits



Soft



Cotton



Breathable

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free



20–30 mmHg
Closed Toe Calf Small Long
in Brown

222CSLW11



20–30 mmHg
Closed Toe Calf Medium Long
in Black

222CMLM99



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12
Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●	<div><div>NAVY 10</div><div>BROWN 11</div><div>BLACK 99</div></div>		
	●	●	●

Product Code Example: 222CMLM10

Series	mmHg Level	Style	Size	Gender	Colour
22	2	C	ML	M	10



Featured colour
is Black

MERINO WOOL 240

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 69

Features and Benefits



Wool



Thermoregulating



Stylish

- High quality extra-fine Merino wool makes this sock soft and breathable
- Moisture wicking capabilities keep feet dry
- Naturally odour resistant
- Thermoregulating for year-round wear
- Smooth, flat toe seam

Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free



20–30 mmHg
Closed Toe Calf Medium Long
in Brown

242CMLM11



20–30 mmHg
Closed Toe Calf Small Long
in Navy

242CSLW10



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●	●
-----------------	---	---	---

Product Code Example: 242CMSW99

Series	mmHg Level	Style	Size	Gender	Colour
24	2	C	MS	W	99



Featured colour
is Lavender

LINEN 250

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 69

Features and Benefits



Hypoallergenic
Natural fibers



Anti-bacterial



Moisture
Management

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free



20–30 mmHg
Closed Toe Calf Medium Long
in Light Grey

252CMLM205



20–30 mmHg
Closed Toe Calf Medium Long
in Denim

252CMLM201



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

Gender: Women = **W**
Men = **M**



Calf (C)	●	●	●	●
Gender	W/M	W/M	W/M	W

Product Code Example: 252CMLW207

Series	mmHg Level	Style	Size	Gender	Colour
25	2	C	ML	W	207



Featured colour
is Honey

MEDIUM SHEER 750

for women

Compression Level

20–30 mmHg

Features and Benefits



Medium Sheer



Durable



Breathable

- Fashionable, durable, medium sheer hosiery for everyday wear
- Excellent breathability and moisture management
- Comfortable 2" top band on calf style

Fiber Content

- 68% Nylon, 32% Spandex
(content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Thigh Medium Long
in Black

752NMLW99



Silicone dotted grip-top on thigh



See product images online
sigvaris.com




SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●

		
HONEY 33	NIGHT-SHADE 94	BLACK 99

Calf (C)	●	●	●
Thigh w/Grip-Top (N) Silicone dotted grip-top band	●	●	●
Pantyhose (P)	●	●	●

Product Code Example: 752NMLW33

Series	mmHg Level	Style	Size	Gender	Colour
75	2	N	ML	W	33



Featured colour is Black

SHEER 780

for women

Compression Levels

15–20 mmHg

20–30 mmHg

30–40 mmHg

Features and Benefits



Sheer



Elegant



Breathable

- Fashionable for professional and dress wear
- Consistent sheerness throughout the leg
- Excellent breathability and moisture management

Fiber Content

- 15–20 and 20–30 mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30–40 mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Latex free



20-30 mmHg
Closed Toe Calf Small Short
in Toasted Almond

782CSSW37



SensinnoV® grip-top on thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	DARK NAVY 08	HONEY 33	WARM SAND 34*	GOLDEN 36	TOASTED ALMOND 37*	CAFÉ 73*	MOCHA 85*	NIGHT-SHADE 94	BLACK 99
Calf (C)	●	●○	●○	●○	●○	●○	●○	●	●○
Thigh w/Grip-Top (N) Sensinnov® grip-top band	●	●○	●○	●○	●○	●○	●○	●	●○
Pantyhose (P)	●	●○	●○	●○	●○	●○	●○	●	●○

* 20–30 mmHg only

Product Code Example: 781NMLW99

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
78	1	N	ML	W	99

Product Code Example: 783PSSO99

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
78	3	P	SS	O	99



Featured colour
is Steel Grey

MICROFIBER 820

for men


Compression Levels

15–20 mmHg

20–30 mmHg

30–40 mmHg

Features and Benefits



Microfiber



Soft



Stylish

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking
- Glides on and off with ease
- Comfortable and durable sock for professional and everyday wear
- Available in extended foot sizes
- Smooth, flat toe seam


Fiber Content

- 87% Nylon , 13% Spandex
(content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Thigh Medium Long
in Black

822NMLM99



Silicone dotted grip-top option
available for calf. Sensinnov®
grip-top for thigh.



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")		NA	
Leg Length	Short	Long	Short	Long	Short	Long		
	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")	≤76 cm (30")	>76 cm (30")		
Size Shoe Size	SS 7–11	SL 9.5–14	MS 7–11	ML 9.5–14	LS 7–11	LL 9.5–14	XS 7–11	XL 9.5–14
				MX* 12.5–17.5		LX* 12.5–17.5		XX* 12.5–17.5

Men's shoe sizes: Short Leg 7–11, Long Leg 9.5–14 (For extended foot sizes: long leg 12.5–17.5)

*MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30 mmHg only. Extended foot sizes not available with grip top

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●***	●
Calf w/Grip-Top (C/S)** Silicone dotted grip-top band	●		●
Thigh w/Grip-Top (N)** Sensinnov® grip-top band	●		●

** Not available in 15–20 mmHg

*** Steel Grey only available in 15–20 mmHg and 20–30 mmHg

Product Code Example: 822CMLM93

Series	mmHg Level	Style	Size	Gender	Colour
82	2	C	ML	M	93

Product Code Example: 822CMLM32/S

Series	mmHg Level	Style	Size	Gender	Colour	Grip-top
82	2	C	ML	M	32	/S



Featured colour
is Black Plaid

MICROFIBER PATTERNS 830


for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 69

Features and Benefits



Microfiber



Stylish



Patterns

- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking.
- Fashionable, durable sock for everyday wear
- Glides on and off with ease
- Smooth, flat toe seam

Fiber Content

- 87% Nylon, 13% Spandex
(content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Calf Small Long
in Royal Blue Argyle

832CSLM48



20–30 mmHg
Closed Toe Calf Medium Short
in Wisteria Windows

832CMSW61



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●

Gender:
Women = **W**
Men = **M**

						
Calf (C)	●	●	●	●	●	●
Gender	W	W	W	M	W/M	W/M
						
Calf (C)	●	●	●	●	●	●
Gender	M	W	W/M	W/M	W	W/M
						
Calf (C)	●	●	●	●		
Gender	W/M	W	W	W/M		

Product Code Example: 832CMLM42

Series	mmHg Level	Style	Size	Gender	Colour
83	2	C	ML	M	42



SOFT OPAQUE 840

for women

Compression Levels

15–20 mmHg

20–30 mmHg

30–40 mmHg

Features and Benefits



Opaque



Microfiber




Soft

- Fashionable opaque hosiery for professional and everyday wear
- Long-lasting microfiber softness provides all day comfort and allows product to easily glide on and off
- Durable and breathable
- Maternity pantyhose option available


Fiber Content

- 79% Nylon, 21% Spandex
(content can vary slightly by style)
- Latex free



20–30 mmHg
Closed Toe Thigh with Grip Top
Medium Long in Pecan

842NMLW40



Silicone dotted grip-top for calf
and Sensinnov® grip-top for thigh



See product images online
sigvaris.com

SIZING CHART

	Small		Medium		Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")	
Calf Length	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Thigh Circumference	45.5–66 cm (18"–26")		48–73 cm (19"–28.5")		52–84 cm (20.5"–33")	
Leg Length	Short	Long	Short	Long	Short	Long
	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")	≤72 cm (28.5")	>72 cm (28.5")
Hip Circumference	83–112.5 cm (32.5"–44.5")		87–127 cm (34.5"–50")		98–150.5 cm (38.5"–59.5")	
Maternity hip circumference	≤142 cm (56")		≤157.5 cm (62")		≤178 cm (70")	
Size	SS	SL	MS	ML	LS	LL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

Closed Toe = ●
Open Toe = ○

	MIDNIGHT BLUE 09*	CHAI 35	PECAN 40	MULBERRY 60*	ESPRESSO 89*	GRAPHITE 91*	BLACK 99
Calf (C)	●	●○	●	●	●	●	●○
Calf w/Grip-Top (C/S)* Silicone dotted grip-top							●
Thigh w/Grip-Top (N) Sensinnov® grip-top band	●	●○	●	●	●	●	●○
Pantyhose (P)	●	●○	●	●	●	●	●○
Maternity Pantyhose (M)*							●

* 20–30 mmHg only

Product Code Example: 842CMSW09

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
84	2	C	MS	W	09

Medical Compression Therapy Solutions



Push Yourself

Why choose **Motion**?

- Cushioned sole for added comfort
- Durable thermoregulating wool (422)
- Moisture-wicking yarns

MOTION PRODUCTS

Series	Product	Page
360	Cushioned Cotton	36
412C	High Tech	38
422	Thermoregulating Wool	40
432	Flow Tech®	42

BEST SUITED FOR

- Management of Chronic Venous Disease (CVD)
- People who stand for prolonged periods of time
- Casual and occupational footwear
- Walking and everyday activities

Featured
colour is
White

CUSHIONED COTTON 360

for women & men

Compression Level

20–30 mmHg

Features and Benefits



Cushioned Sole



Cotton



Sensitive Skin

- Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Ideal for sensitive skin and allergies
- Smooth, flat toe seam
- Diabetic friendly
- Ideal for people who stand for prolonged periods of time

Fiber Content

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Latex free



20–30 mmHg
Closed Toe Calf Medium Long
in Black

362CMLW99



Flat, low friction toe seam



See product images online
sigvaris.com

SIZING CHART


	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●



Calf (C)	●	●
-----------------	---	---

Product Code Example: 362CMLM00

Series	mmHg Level	Style	Size	Gender	Colour
36	2	C	ML	M	00



Featured colour
is Steel Blue

HIGH TECH 412C

for women & men

Compression Level

20–30 mmHg

Features and Benefits



**Moisture
Wicking**



**Anatomically
Contoured**



Colours

- High-Tech fibers for moisture wicking and improved thermoregulation
- Cushioned footbed, toe and Achilles tendon zones for optimal comfort
- Anatomically contoured to support the natural shape of the foot and calf (left and right foot specific)
- Smooth, flat toe seam
- Perfect for everyday activities and active individuals


Fiber Content

- 64% Nylon, 17% Polyester, 19% Spandex
- Latex free



20–30 mmHg
Closed Toe Calf Medium Medium
in White

412CMM00



20–30 mmHg
Closed Toe Calf Small Long
in Pink

412CSL56










See product images online
sigvaris.com

UNISEX SIZING CHART

	Small			Medium			Large			X Large		
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")			29–33 cm (11.5"–13")		
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")			40.5–58 cm (16"–23")		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL	XS	XM	XL

STYLES AND COLOURS

Closed Toe = ●

						
WHITE 00	BLUE 50	LIME 54	PINK 56	RED 59	STEEL BLUE 62	BLACK 99

Calf (C)	●	●	●	●	●	●	●
-----------------	---	---	---	---	---	---	---

Product Code Example: 412CML62

Series	mmHg Level	Style	Size	Colour
41	2	C	ML	62



Featured
colour is
Charcoal

THERMO-REGULATING WOOL 422

for women & men

Compression Level

20–30 mmHg

Also available in 15–20 mmHg, see page 71

Features and Benefits



Wool



Cushioned Sole



Thermoregulating

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned soles help prevent blisters
- Soft and comfortable
- Smooth, flat toe seam
- Perfect for outdoor work environments and activities

Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free



20–30 mmHg
Closed Toe Calf Medium Medium
in Charcoal

422CMM12



20–30 mmHg
Closed Toe Calf Medium Short
in Olive

422CMS31



See product images online
sigvaris.com

UNISEX SIZING CHART

	Small			Medium			Large		
Ankle Circumference	18–21.5 cm (7"–8.5")			21.5–25.5 cm (8.5"–10")			25.5–29 cm (10"–11.5")		
Calf Circumference	28–39 cm (11"–15.5")			31.5–45 cm (12.5"–17.5")			35.5–51 cm (14"–20")		
Women's Shoe Size	4.5–7.5	8–10.5	11–13.5	4.5–7.5	8–10.5	11–13.5	8–10.5	11–13.5	>14
Men's Shoe Size	≤6.5	7–9.5	10–12.5	≤6.5	7–9.5	10–12.5	7–9.5	10–12.5	13–15.5
Size	SS	SM	SL	MS	MM	ML	LS	LM	LL

STYLES AND COLOURS

Closed Toe = ●

CHARCOAL
12OLIVE
31

Calf (C)



Product Code Example: 422CSL12

Series	mmHg Level	Style	Size	Colour
42	2	C	SL	12



Featured colour is Orchid

FLOW TECH® 432

for women & men

Compression Level

20–30 mmHg

Features and Benefits



Breathable



Moisture Wicking



Cushioned Sole

- MeshLink™ channels are designed to deliver optimal airflow to the parts of the foot and leg that generate the most heat, linking the bottom of the foot to the top of the calf and allowing heat to escape
- Cushioned footbed and toe with higher Achilles tendon zones for optimal comfort
- SENSIL® Aquarius fibers keep your legs and feet dry with next-generation moisture management
- Mesh flex zone across the instep for donning ease
- Smooth, flat toe seam

Fiber Content

- 62% Nylon, 24% SENSIL Aquarius Nylon, 14% Spandex
- Latex free



20–30 mmHg
Closed toe Calf Medium Long
in Deep Sea



MeshLink™ channels, our unique pattern design for optimal airflow



See product images online
sigvaris.com

432CML212









UNISEX SIZING CHART

	Small		Medium		Large		X Large	
Ankle Circumference	18–21.5 cm (7"–8.5")		21.5–25.5 cm (8.5"–10")		25.5–29 cm (10"–11.5")		29–36 cm (11.5"–14")	
Calf Circumference	28–39 cm (11"–15.5")		31.5–45 cm (12.5"–17.5")		35.5–51 cm (14"–20")		41–61 cm (16"–24")	
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long
	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")	≤41 cm (16")	>41 cm (16")
Size	SS	SL	MS	ML	LS	LL	XS	XL

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = ●				
Calf (C)	●	●	●	●
Closed Toe = ●				
Calf (C)	●	●	●	●

Product Code Example: 432CML210

Series	mmHg Level	Style	Size	Colour
43	2	C	ML	210

sigvaris SPECIALTY

Medical Compression Therapy Solutions



Take care of Yourself

Why choose **Specialty**?

- Specific medical needs
- High Containment
- Natural rubber

SPECIALTY PRODUCTS

Series	Product	Page
500	Natural Rubber	46
910	Advance Armsleeve	48

BEST SUITED FOR

- Severe varicosities
- Chronic edema post DVT
- Severe venous stasis
- Early stages of lymphedema
- Active venous ulcers and prevention of ulcer recurrence
- People with early stages of arm lymphedema
- Post-mastectomy care



NATURAL RUBBER 500

for women & men

Compression Levels

30–40 mmHg

40–50 mmHg

50–60 mmHg

Features and Benefits



Natural Rubber



Durable



Largest Size Range

- Durable, long-term containment solution with minimal fiber fatigue throughout the day
- Extended sizing for a custom-like fit
- **Ideal for patients with more serious indications, including:**
 - Severe varicosities
 - Chronic edema post DVT
 - Prevention and treatment of Active Venous Ulcers
 - Severe venous stasis
 - Lymphedema


Fiber Content

- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh-high styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibers (Pantyhose style)
- **Caution:** contains natural rubber latex that may cause allergic reaction



30–40 mmHg
Open Toe Thigh Medium Average
Long in Beige

503NM2077



Silicone dotted grip-top available
on thigh



See product images online
sigvaris.com

SIZING CHART

	Small				Medium				Large				X Large			
Ankle Circumference	19–23 cm (7.5"–9")				23–30 cm (9"–12")				30–40 cm (12"–15.5")				40–46 cm (15.5"–18")			
Calf Circumference	Average		Full		Average		Full		Average		Full		Average		Full	
	28–38 cm (11"–15")		≤43 cm (17")		32–43 cm (12.5"–17")		≤48 cm (19")		36–48 cm (14.5"–19")		≤57 cm (22.5")		45–57 cm (18"–22.5")		≤68 cm (27")	
Calf Circumference for 505C (Full only)	505C: full calf only		32–43 cm (12.5"–17")		505C: full calf only		36–48 cm (14.5"–19")		505C: full calf only		45–57 cm (18"–22.5")		505C: not available in XL			
Calf Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long
	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")	≤40 cm (≤15.5")	>40 cm (>15.5")
Thigh Circumference	Average		Full		Average		Full		Average		Full		N/A			
	41–54 cm (16"–21")		≤61 cm (24")		48–61 cm (19"–24")		≤71 cm (28")		58–71 cm (23"–28")		≤81 cm (32")					
Thigh Circumference for 505T (Full only)	505T: full thigh only		≤61 cm (24")		505T: full thigh only		≤71 cm (28")		505T: full thigh only		≤81 cm (32")					
Leg Length	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long	Short	Long				
	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")	≤71 cm (28")	>71 cm (28")				
Hip Circumference	<122 cm (48")				<137 cm (54")				<155 cm (61")							
Size	S1	S2	S3	S4	M1	M2	M3	M4	L1	L2	L3	L4	X1	X2	X3	X4

STYLES AND COLOURS

Open Toe = ○



Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications, cautions, and contraindications, please see page 82 for details. Sigvaris accepts no liability for non-observance of contraindications or cautions.

Calf (C)	○	Calf only available in S, M, and L full circumference for 50–60 mmHg
Thigh w/out Grip-Top (T)	○	
Thigh w/Grip-Top (N) Silicone dotted grip-top band	○	Not available in 50-60 mmHg or XL
Pantyhose (P)	○	Pantyhose only available in 30–40 mmHg and only in full thigh circumference
Thigh w/Waist Attachment (W) Indicate Leg LEFT (/L) or RIGHT (/R)	○	Thigh w/ waist attachment is available in average and full thigh circumferences in 30–40 mmHg and in full thigh circumferences only in 40–50 mmHg

Product Code Example: 504CS2O77

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour
50	4	C	S2	O	77

Product Code Example: 503WL1O77/L

Series	mmHg Level	Style	Size	Gender/Open Toe	Colour	Leg
50	3	W	L1	O	77	/L



ADVANCE ARMSLEEVE 910

for women & men

Compression Levels

15–20 mmHg

20–30 mmHg

Features and Benefits



Moisture
Wicking



SensinnoV
Grip-top Band



Breathable

- Ideal for patients with arm lymphedema or post-mastectomy care
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable and soft
- SensinnoV® grip-top band ensures armsleeve stays in place all day long
- Available with or without gauntlet

Fiber Content

- 50% Nylon, 28% Rayon, 22% Spandex
- Latex free



With
Gauntlet



SensinnoV® grip-top band

20–30 mmHg
Armsleeve With Gauntlet
Small Long in Beige

912ASLG77/S

SIZING CHART

	Small		Medium		Large		X Large	
Palm (A)	17–21 cm (6.5"–8")		18–22 cm (7"–8.5")		19–23 cm (7.5"–9")		20–24 cm (8"–9.5")	
Wrist (C)	14–16 cm (5.5"–6.5")		16–18 cm (6.5"–7")		18–20 cm (7"–8")		20–23 cm (8"–9")	
Upper Mid-Arm (F)	Regular		Regular		Plus		Regular	
	26–34 cm (10"–13.5")		29–37 cm (11.5"–14.5")		37–43 cm (14.5"–17")		40–46 cm (15.5"–18")	
Arm Length (C-G)*	Regular	Long	Regular	Long	Regular	Long	Regular	Long
	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")	≤43 cm (17")	>43 cm (17")
Size	SR	SL	MR	ML	MR/P	ML/P	LR	LL

* Should have **slight bend** in arm when measuring length


STYLES AND COLOURS

Open Sleeve = ○



Armsleeve w/Gauntlet (G) Sensinnov® grip-top band	○
Armsleeve w/out Gauntlet (O) Sensinnov® grip-top band	○

Note: all armsleeves come with Grip-top (/S)
15–20 mmHg not available with gauntlet or plus sizing



ARMSLEEVE WITHOUT GAUNTLET

Product Code Example: 912ASLO77/S

Series	mmHg Level	Style	Size	W/out Gauntlet	Colour	Grip-Top
91	2	A	SL	O	77	/S

ARMSLEEVE WITH GAUNTLET

Product Code Example: 912ALRG77/PS

Series	mmHg Level	Style	Size	W/Gauntlet	Colour	Plus Size	Grip-Top
91	2	A	LR	G	77	/P	S

sigvaris INELASTIC WRAPS AND MEDICALLY COMPLEX EDEMA WEAR

Solutions for the active treatment and management of
conditions impacting the lymphatic and venous systems

Move your
life with
WrapAbility



Why Choose Inelastic Wraps and Medically Complex Edema Wear?

- Adjustable pressure
- Large variety of solutions

INELASTIC WRAPS

Product	Page
How Inelastic Wraps Work	52
Coolflex®	53
Compreflex®	55
Compreflex® Options	57
Foot Sizing Charts	58
Foot Compression Garments	60
Liners	60
Accessories	61
Compreshort®	63
Comprevest®	63

The SIGVARIS GROUP line of products for medically complex edema (MCE) are for the treatment and management of conditions impacting the lymphatic and venous circulatory systems.


The product portfolio contains ready-to-wear garments for the active treatment and maintenance of lymphedema and for patients managing chronic venous disorders in the lower extremities. Sigvaris is looking towards a brighter future with new therapies for complex venous and lymphatic diseases. Our inelastic wraps deliver maximum therapeutic results with increased patient compliance.

To view our entire product portfolio, please contact our Customer Care Team or download a digital copy on the Dealer Portal home page.

HOW INELASTIC WRAPS WORK

EFFECTIVE WORKING PRESSURE

- Inelastic compression reduces leg edema and improves venous blood flow.
- Inelastic compression can increase arterial perfusion in patients with mild to moderate peripheral arterial disease.



SAFE, PRECISE RESTING PRESSURE


- Under the supervision of a healthcare provider, inelastic compression can be used safely for patients with mild to moderate peripheral arterial disease.

HOW COMPREFLEX® IS PROVEN TO ACCURATELY SET COMPRESSION LEVELS

Accutab®



Accutab® has been rigorously tested according to established guidelines in the International Compression Club Consensus documents. Testing was conducted using PicoPress, which is recognized for its high accuracy and reproducibility.



HOW COOLFLEX® IS PROVEN TO ACCURATELY SET COMPRESSION LEVELS

EZ-Close™

The integrated **EZ-Close™** Pressure System has visual indicators for the accurate delivery of 20-30 mmHg or 40-50 mmHg on each tab.



COOLFLEX®

20–50 mmHg



ONE-HANDED DONNING

EZ-Close™ Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues

BREATHABLE GARMENT

The Spacer Fabric™ is designed to be highly breathable for cool-wearing comfort

INNOVATIVE DESIGN

Durable, metal hook closures for maximum strength and durability

Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Well-managed Lymphedema
- Non-specific or post-surgical edema
- Patients with limited dexterity, arthritis, and mobility limitations

Features and benefits

- Uses metal hooks, increasing durability and long-term wear
- Spacer Fabric™ is highly breathable for cool-wearing comfort
- Features a stretch panel inside the garment for easy one-handed donning and doffing



Coolflex® Standard

CALF

20–50 mmHg



Ideal for both venous and lymphatic limb disorders, the garments are available in standard and max sizing and can accommodate limbs that feature disproportionately small ankles with a larger calf.

Features

- Lightweight and breathable for cool-wearing comfort
- Easy to put on for patients with limited dexterity
- Uses metal hooks, increasing durability and long-term wear

Includes

- Coolflex® Standard Calf
- Built-in application sleeve for the independent use and wear
- Pair of non-compressive Basic Liners per leg




Basic Liner (included)

When pairing with a **Basic Liner**, a foot piece is required. See **Coolflex® Foot Piece** on page 59; sold separately.



Transition Liner

When pairing with a **Transition Liner**, a foot piece is not required. You can substitute the Basic Liner for a **Transition Liner** (page 60) by contacting customer care.



COOLFLEX® STANDARD CALF SIZE CHART (A–C1)

	Small	Small Max	Medium / Large	Medium / Large Max	X Large / XX Large
C1	28–39 cm	35–50 cm	34–50 cm	45–64 cm	45–64 cm
C	28–39 cm	35–50 cm	34–50 cm	45–64 cm	45–64 cm
B	23–33 cm	27–38 cm	29–43 cm	35–49 cm	37–54 cm
A	16–26 cm	16–26 cm	23–33 cm	23–33 cm	30–42 cm

G Length:

Regular: 30–35.9 cm

Tall: 36–42 cm

PRODUCT CODE (Available in Black Only)

Left	Regular	2101-NF-BKR-L	2101-NF-BMR-L	2103-NF-BKR-L	2103-NF-BMR-L	2105-NF-BKR-L
	Tall	2101-NF-BKT-L	2101-NF-BMT-L	2103-NF-BKT-L	2103-NF-BMT-L	2105-NF-BKT-L
Right	Regular	2101-NF-BKR-R	2101-NF-BMR-R	2103-NF-BKR-R	2103-NF-BMR-R	2105-NF-BKR-R
	Tall	2101-NF-BKT-R	2101-NF-BMT-R	2103-NF-BKT-R	2103-NF-BMT-R	2105-NF-BKT-R

COMPREFLEX®

20–50 mmHg



INTELLIGENT DESIGN

Garment features a contoured calf, providing a superior fit

PREMIUM MATERIALS

Soft, conforming material provides a comfortable fit

EASY APPLICATION

With a front stretch panel to secure the garment in place and straps that roll back, the garment is easy to don and adjust, improving compliance

SUPERIOR QUALITY

Each garment is double stitched to ensure durability



Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Lymphedema - Active and Maintenance
- Non-specific or post surgical edema
- Patients with fragile and weakened skin
- Opened/Healed wounds
- Non-compliance due to difficulty of donning
- Those whose measurements fall outside sizes for compression socks

Features and benefits

- High quality inelastic wrap composed of durable breathable inelastic material with moisture wicking capacities
- Contoured design helps the wrap conform to the natural shape of the leg
- Front stretch panel secures the garment in place, while the Velcro straps assist for easy application
- Ability to customize with different liner and boot options

Compreflex® Standard

CALF & FOOT

20–50 mmHg



For those who require a calf and foot garment

Includes

- Compreflex® Wrap
- Pair of non-compressive Basic Liners
- Non attached Compreboot® Standard Boot



Accutab® Devices

(See page 61 for details)




Basic Liner

(See page 60 for details)



Compreboot® Standard Boot

(See page 59 for details)



COMPREFLEX® CALF SIZING CHART (A, C)

	Small	Medium	Large	X Large	XX Large
C	29–39 cm	34–44 cm	40–50 cm	46–56 cm	54–64 cm
A	19–26 cm	21–30 cm	26–36 cm	31–41 cm	36–46 cm

G Length:

Regular: 30–33.9 cm

Tall: 34–42 cm

Foot Size

I: Measure circumference of ankle bend and heel

J: Measure circumference across metatarsal heads

K: Measure **medial** length from

*See page 58 for foot sizing chart

PRODUCT CODE

Size Regular	Black	Beige	Size Tall	Black	Beige
Small	1401-BKR	1411-BKR	Small	1401-BKT	1411-BKT
Medium	1402-BKR	1412-BKR	Medium	1402-BKT	1412-BKT
Large	1403-BKR	1413-BKR	Large	1403-BKT	1413-BKT
X Large	1404-BKR	1414-BKR	X Large	1404-BKT	1414-BKT
XX Large	1405-BKR	1415-BKR	XX Large	1405-BKT	1415-BKT

Note:

- 1401-BKR, 1411-BKR, 1401-BKT, and 1411-BKT include a Regular Length Boot
- 1402-BKR, 1412-BKR, 1402-BKT, 1412-BKT, 1403-BKR, 1413-BKR, 1403-BKT, and 1413-BKT include a Medium/Large Regular Length Boot
- 1404-BKR, 1414-BKR, 1404-BKT, 1414-BKT, 1405-BKR, 1415-BKR, 1405-BKT, and 1415-BKT include a X-Large/XX-Large Regular Length Boot

Compreflex® Calf Options

Standard CALF

20–50 mmHg

For those who require an alternative foot garment

Includes:

- Compreflex® Wrap
- Pair of non-compressive Basic Liners

Can be matched with alternative foot pieces* (see page 59).

* When paired with a basic liner, a foot piece is highly recommended



PRODUCT CODE

		Black	Beige
Size Regular	Small	1401-NF-BKR	1411-NF-BKR
	Medium	1402-NF-BKR	1412-NF-BKR
	Large	1403-NF-BKR	1413-NF-BKR
	X Large	1404-NF-BKR	1414-NF-BKR
	XX Large	1405-NF-BKR	1415-NF-BKR
Size Tall	Small	1401-NF-BKT	1411-NF-BKT
	Medium	1402-NF-BKT	1412-NF-BKT
	Large	1403-NF-BKT	1413-NF-BKT
	X Large	1404-NF-BKT	1414-NF-BKT
	XX Large	1405-NF-BKT	1415-NF-BKT

Transition CALF

20–50 mmHg

For those who require a compressive liner for the foot and ankle

Includes:

- Compreflex® Wrap
- Pair of Transition Liners, which provide 10–15 mmHg in the foot and ankle

Ability to wear regular shoe over liner



PRODUCT CODE

		Black	Beige
Size Regular	Small	1401-UC-BKR	1411-UC-BKR
	Medium	1402-UC-BKR	1412-UC-BKR
	Large	1403-UC-BKR	1413-UC-BKR
	X Large	1404-UC-BKR	1414-UC-BKR
	XX Large	1405-UC-BKR	1415-UC-BKR
Size Tall	Small	1401-UC-BKT	1411-UC-BKT
	Medium	1402-UC-BKT	1412-UC-BKT
	Large	1403-UC-BKT	1413-UC-BKT
	X Large	1404-UC-BKT	1414-UC-BKT
	XX Large	1405-UC-BKT	1415-UC-BKT

Complete CALF

20–50 mmHg

For those who require a compressive liner with compression in the foot and ankle, and graduated compression in the calf

Includes:

- Compreflex® Wrap
- Pair of Complete Liners, which provide 10–15 mmHg in the foot and ankle, and graduated compression in the calf

Ability to wear regular shoe over liner




PRODUCT CODE

		Black	Beige
Size Regular	Small	1401-FC-BKR	1411-FC-BKR
	Medium	1402-FC-BKR	1412-FC-BKR
	Large	1403-FC-BKR	1413-FC-BKR
	X Large	1404-FC-BKR	1414-FC-BKR
	XX Large	1405-FC-BKR	1415-FC-BKR
Size Tall	Small	1401-FC-BKT	1411-FC-BKT
	Medium	1402-FC-BKT	1412-FC-BKT
	Large	1403-FC-BKT	1413-FC-BKT
	X Large	1404-FC-BKT	1414-FC-BKT
	XX Large	1405-FC-BKT	1415-FC-BKT

See page 60 for Liner Sizing chart

Foot Sizing Charts



Foot Size:

I: Measure circumference of ankle bend and heel

J: Measure circumference across metatarsal heads

K: Measure **medial** length from heel to 1st metatarsal head

COMPREBOOT® (I, J, K)

	Small / Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long
I	28–36 cm	28–36 cm	≤44 cm	≤44 cm
J	22–26 cm	22–26 cm	≤33 cm	≤33 cm
K	14–18 cm	19–23 cm	20–22 cm	25–27 cm

COMPREBOOT® STANDARD AND PLUS (I, J, K)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–36cm	28–36 cm	≤39 cm	≤39 cm	≤44 cm	<44 cm
J	22–26 cm	22–26 cm	≤30 cm	≤30 cm	≤33 cm	≤33 cm
K	14–18 cm	19–23 cm	18–20 cm	23–25 cm	20–22 cm	25–27 cm

COOLFLEX® STANDARD FOOT (I, J, K)

(Available in Black Only)

	Small		Medium / Large		X Large / XX Large	
	Regular	Long	Regular	Long	Regular	Long
I	28–39 cm	N/A	30–43 cm	N/A	34–49 cm	N/A
J	22–26 cm	N/A	25–30 cm	N/A	28–33 cm	N/A
K	14–18 cm	19–21 cm	18–20 cm	20–22 cm	20–22 cm	23–25 cm
Right	2101-BTR-R	2101-BTL-R	2103-BTR-R	2103-BTL-R	2105-BTR-R	2105-BTL-R
Left	2101-BTR-L	2101-BTL-L	2103-BTR-L	2103-BTL-L	2105-BTR-L	2105-BTL-L

Foot Compression Garments

Foot compression garments below are pictured with the Basic Liner (not included). Basic Liner sold separately (see page 60) or with the purchase of an Inelastic Wrap.

Compreboot® Lite FOOT

Variable compression



- Ultra-thin profile that easily fits into most footwear
- Highly adjustable design may be trimmed to fit
- Provides mild, therapeutic compression



PRODUCT CODE (Available in Black Only)

Size	Length	Black
Small / Medium / Large	Regular	1803-BLR
	Long	1803-BLL
X Large / XX Large	Regular	1805-BLR
	Long	1805-BLL

Compreboot® Standard FOOT

Variable compression



- Four-way stretch for firm compression around ankle
- Ideal for mild to moderate chronic edema
- Low stretch material conforms to foot shape



PRODUCT CODE

Size	Length	Black	Beige
Small	Regular	1801-BTR	1811-BTR
	Long	1801-BTL	1811-BTL
Medium / Large	Regular	1803-BTR	1813-BTR
	Long	1803-BTL	1813-BTL
X Large / XX Large	Regular	1805-BTR	1815-BTR
	Long	1805-BTL	1815-BTL

Compreboot® Plus FOOT

Variable compression



- Chip foam lining creates high / low pressure ridges
- Removable gripper patches provide traction
- Provides moderate, therapeutic compression
- Used for Lymphedema stage I, II, III



PRODUCT CODE (Available in Black Only)

Size	Length	Black
Small	Regular	1801-BPR
	Long	1801-BPL
Medium / Large	Regular	1803-BPR
	Long	1803-BPL
X Large / XX Large	Regular	1805-BPR
	Long	1805-BPL

Coolflex® Standard FOOT

Variable compression



- Low profile design allows boot to fit in most footwear
- Spacer Fabric™ is highly breathable for cool-wearing comfort
- Durable metal clasps, instead of Velcro closure



PRODUCT CODE (Available in Black Only)

Size	Length	Left (black)	Right (black)
Small	Regular	2101-BTR-L	2101-BTR-R
	Long	2101-BTL-L	2101-BTL-R
Medium / Large	Regular	2103-BTR-L	2103-BTR-R
	Long	2103-BTL-L	2103-BTL-R
X Large / XX Large	Regular	2105-BTR-L	2105-BTR-R
	Long	2105-BTL-L	2105-BTL-R

Liner Options

Basic Liner

No compression

- Non-compressive liner for the protection of the skin
- Easy to apply knee-high liners
- drirelease moisture wicking and odour control
- Flat seams for added comfort
- Available in Heathered Silver colour only



PRODUCT CODE

(Included with Compreflex® Standard Calf & Foot, Compreflex® Standard Calf, and Coolflex® Standard Calf)

Size	Pair (gray)	10 pairs (gray)
Universal	1700 - NK1	1700 - NKC

Transition Liner

10–15 mmHg

- Mild compression in foot and ankle only
- drirelease moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Available in Black and Beige



PRODUCT CODE

(Included with Compreflex® Transition)

Size	Pair (black)	Pair (beige)
Small	1601-CK	1611-CK
Medium	1602-CK	1612-CK
Large	1603-CK	1613-CK
X Large	1604-CK	1614-CK
XX Large	1605-CK	1615-CK

* When ordering a Coolflex® standard without a foot piece, please order a transition liner.

Complete Liner

10–15 mmHg

- Mild compression in foot and ankle with graduated compression in calf
- Comfortable knee-high liners
- drirelease moisture wicking and odour control
- Extra thick sole from toe to heel for added comfort
- Available in Black and Beige




PRODUCT CODE

(Included with Compreflex® Complete)

Size	Pair (black)	Pair (beige)
Small	1601-CK-F1	1611-CK-F1
Medium	1602-CK-F1	1612-CK-F1
Large	1603-CK-F1	1613-CK-F1
X Large	1604-CK-F1	1614-CK-F1
XX Large	1605-CK-F1	1615-CK-F1

TRANSITION AND COMPLETE LINERS SIZING CHART

Size	Ankle Circ. (A)	Calf Circ. (C)	Women's Shoe	Men's Shoe
Small	19–26 cm	29–39 cm	<7½	<6½
Medium	21–30 cm	34–44 cm	7.5–9.5	6.5–8.5
Large	26–36 cm	40–50 cm	10–12	9–11
X Large	31–41 cm	46–56 cm	12.5–15	11.5–14
XX Large	36–46 cm	54–64 cm	12.5–15	11.5–14



ACCESSORIES


Compreflex® Strap Extender

- Provides up to 10 cm of additional circumferential size range to Compreflex® calf wraps
- Detaches easily
- Potential solution for patients outside of standard size ranges of Compreflex® calf wraps or needing custom size



PRODUCT CODE		
Qty	Black	Beige
1 unit	1401-SE1	1411-SE1
10 Unit Pack	1401-SEB	1411-SEB

Accutab® Devices




The Accutab® accurately delivers the level of compression prescribed by the health care provider. An easy to use device, the Accutab® ensures the safe delivery of compression in doses of 20–30 mmHg, 30–40 mmHg, or 40–50 mmHg.

PRODUCT CODE	
Size	1 set
Small	1491-TB1
Medium	1492-TB1
Large	1493-TB1
X Large	1494-TB1
XX Large	1495-TB1



COMPRESHORT®

10–15 mmHg



Product available in short, capri and custom lengths

Ideal for

- Lymphedema
- Lipedema

Colour



Black

Features and benefits


- The Compreshort® offers mild, supportive compression for the treatment of abdominal or truncal and upper thigh edema. Patients with lymphedema and lipedema benefit from the therapeutic compression and support the garment provides.
- For custom sizes, please contact customer care on our full product portfolio for custom garments.

COMPRESHORT® SIZING CHART

(Waist, Hip, Midthigh)

PRODUCT CODE

	Waist	Hip	Midthigh	Product SKU (Short)	Product SKU (Capri)
Small	69–76 cm	84–91 cm	41–46 cm	1301-CS	1301-CP
Small Max	76–86 cm	91–102 cm	41–46 cm	1301-CM	1301-CPM
Medium	76–86 cm	91–102 cm	46–53 cm	1302-CS	1302-CP
Medium Max	86–97 cm	102–112 cm	46–53 cm	1302-CM	1302-CPM
Large	86–97 cm	102–112 cm	53–61 cm	1303-CS	1303-CP
Large Max	97–112 cm	112–127 cm	53–61 cm	1303-CM	1303-CPM
X Large	97–112 cm	112–127 cm	61–71 cm	1304-CS	1304-CP
X Large Max	112–132 cm	127–147 cm	61–71 cm	1304-CM	1304-CPM
2X Large	112–132 cm	127–147 cm	71–81 cm	1305-CS	1305-CP
2X Large Max	132–152 cm	147–168 cm	71–81 cm	1305-CM	1305-CPM
3X Large	132–152 cm	147–168 cm	81–91 cm	1306-CS	1306-CP
3X Large Max	152–172 cm	168–188 cm	81–91 cm	1306-CM	1306-CPM
4X Large	152–172 cm	168–188 cm	91–101 cm	1307-CS	1307-CP
4X Large Max	172–193 cm	188–208 cm	91–101 cm	1307-CM	1307-CPM
5X Large	172–193 cm	188–208 cm	101–112 cm	1308-CS	1308-CP



COMPREVEST®



- Straps are easily adjusted
- Anterior hook adjustment for optimal fit

Ideal for

- Post-surgical procedures including mastectomy, lumpectomy or radical reconstruction
- Lymphedema

Colours



Features and benefits

- The Comprevest® offers mild to moderate compression to the chest and back. Designed with features for an improved fit and overall comfort; the Comprevest includes wide straps and adjustable hook closures.
- Offered in multiple band and cup sizes.
- For added stimulation of the affected area, pair the Comprevest® with the **Chip Pad Bra**. Chip Pad Bra is designed to be inserted into the interior side pocket of the Comprevest. With chip foam sewn into channels, the it massages and breaks down fibrotic tissue caused by edema or radiation, while the channeling directs the fluid for drainage.


COMPREVEST® SIZING CHART / PRODUCT CODE

		Small	Medium	Large	X Large	XX Large	
Band Size		28–30 in	32–34 in	36–38 in	40–42 in	44–46 in	
Cup Size	BLACK	A/B	2201-BRR-A	2202-BRR-A	2203-BRR-A	2204-BRR-A	2205-BRR-A
		C/D	2201-BRR-C	2202-BRR-C	2203-BRR-C	2204-BRR-C	2205-BRR-C
		DD/DDD	2201-BRR-D	2202-BRR-D	2203-BRR-D	2204-BRR-D	2205-BRR-D
	WHITE	A/B	2281-BRR-A	2282-BRR-A	2283-BRR-A	2284-BRR-A	2285-BRR-A
		C/D	2281-BRR-C	2282-BRR-C	2283-BRR-C	2284-BRR-C	2285-BRR-C
		DD/DDD	2281-BRR-D	2282-BRR-D	2283-BRR-D	2284-BRR-D	2285-BRR-D

HOW TO MEASURE FOR SIZE

Step 1

Measure (in inches) the circumference at the fullest part of bust



Step 2

Measure (in inches) the circumference directly under the bust



Step 3

Determine the cup by subtracting the band from the bust

2" or less = **A/B**
2" - 4" = **C/D**
4" - 6" = **DD/DDD**

CHIP PAD BRA SIZING CHART

Size	Cup	Extension	Total Length	Colour	Product Code
Half	20cm x 20cm	15cm	35cm	BLACK	2230-BRH
				WHITE	2230-BRF
Full		30cm	50cm	BLACK	2280-BRH
				WHITE	2280-BRF

Chip Pad Bra

Easily inserted into the interior side pocket of the Comprevest®.



Compression socks and hosiery for everyday wellbeing



Care for Yourself

Why choose **Consumer Health**?

- Tired, heavy, restless, aching legs
- Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease, prevention of swelling and varicose veins
- Simplified sizing

CONSUMER HEALTH PRODUCTS

Series	Product	Page
120	Sheer Fashion	66
170	Soft Silhouette Leggings	67
146 / 186	Casual Cotton	68
146/T / 186/T	Traveno Travel Socks	68
151 / 191	Sea Island Cotton	68
152 / 192	All Season Merino Wool	69
153 / 193	Linen	69
143 / 183	Microfiber Shades	69
401	Athletic Recovery Socks	70
412V	Performance Sleeves	70
421	Merino Outdoor Socks	71
160	Eversoft Diabetic Socks	71
–	Standard Sizing Chart	72
–	Unisex Sizing Chart	73

SHEER FASHION 120

for women

15-20 mmHg

Features and Benefits

- Incredible sheer design in a variety of fashionable colours
- Sheer toe and heel area
- Comfortable stretch and soft fabric, excellent for everyday wear

Fiber Content

- 79% Nylon, 21% Spandex (Content can vary slightly by style)
- Latex free



Featured
colour is
Honey

Closed Toe = ●
Open Toe = ○

	CHARCOAL 12	TAUPE 29	HONEY 33	GOLDEN 36	CAFÉ 73	BLACK 99
Calf (C)	●	●●	●●	●●	●	●
Thigh (N)	●	●	●	●	●	●
Pantyhose (P)	●	●	●	●	●	●

Product Code

Example: 120PC33

Series	Style	Size	Colour
120	P	C	33

CALF / THIGH SIZING CHART

Women's Shoe	5-7	7.5-9.5	10-12
Ankle Circumference	18-23 cm (7"-9")	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")
Size	A	B	C

Weight Guidelines For Thigh Style: Size A: <135Lbs ; Size B: <180Lbs ; Size C: <230Lbs.

PANTYHOSE SIZING CHART

cm/ft	Kg/Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
145 cm / 4'9"																												
147 cm / 4'10"																												
150 cm / 4'11"																												
152 cm / 5'0"																												
155 cm / 5'1"																												
157 cm / 5'2"																												
160 cm / 5'3"																												
162 cm / 5'4"																												
165 cm / 5'5"																												
167 cm / 5'6"																												
170 cm / 5'7"																												
173 cm / 5'8"																												
175 cm / 5'9"																												
178 cm / 5'10"																												
180 cm / 5'11"																												
183 cm / 6'0"																												

SIZE A

ANKLE
CIRCUMFERENCE
19-23 cm (7.5"-9")

SIZE B

ANKLE
CIRCUMFERENCE
19-23 cm (7.5"-9")

SIZE C

ANKLE
CIRCUMFERENCE
20.5-25 cm (8"-10")

SIZE D

ANKLE
CIRCUMFERENCE
20.5-25 cm (8"-10")

SIZE E

ANKLE
CIRCUMFERENCE
22-27 cm (8.5"-10.5")

SIZE F


ANKLE
CIRCUMFERENCE
22-27 cm (8.5"-10.5")



15-20 mmHg

Closed Toe Calf Size C
in Charcoal

120CC12



15-20 mmHg

Closed Toe Pantyhose Size A
in Black

120PA99

SOFT SILHOUETTE LEGGINGS 170

for women

15-20 mmHg

Features and Benefits

- Advanced soft stretch microfiber fabric shapes legs and smooths skin
- Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort
- Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts

Fiber Content

- 83% Nylon, 17% Spandex
- Latex free



Non-binding ankle cuff
stays comfortably in place
for all day wear

Open Foot = ☐BLACK
99Leggings (L) ☐

Product Code Example: 170LB99

Series	Style	Size	Colour
170	L	B	99

LEGGING SIZING CHART

cm / ft	Kg / Lbs	45 / 100	48 / 105	50 / 110	52 / 115	54 / 120	57 / 125	59 / 130	61 / 135	64 / 140	66 / 145	68 / 150	70 / 155	73 / 160	75 / 165	77 / 170	80 / 175	82 / 180	84 / 185	86 / 190	89 / 195	91 / 200	93 / 205	95 / 210	98 / 215	100 / 220	102 / 225	105 / 230
147 cm / 4'10"																												
150 cm / 4'11"																												
152 cm / 5'0"																												
155 cm / 5'1"																												
157 cm / 5'2"																												
160 cm / 5'3"																												
162 cm / 5'4"																												
165 cm / 5'5"																												
167 cm / 5'6"																												
170 cm / 5'7"																												
173 cm / 5'8"																												
175 cm / 5'9"																												
178 cm / 5'10"																												
180 cm / 5'11"																												
183 cm / 6'0"																												

Guidelines for Ankle Circumference: **Size A:** 18–23 cm (7"–9"), **Size B & C:** 20–26 cm (8"–10"), **Size D:** 22–29 cm (8.5"–11.5")

CASUAL COTTON 146/186

for women & men


15-20 mmHg

Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation

Fiber Content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free



Calf (C)
Closed toe

Featured
colour
is Black

Product Code Example: 146CB99

Series	Style	Size	Colour
146	C	B	99

TRAVENO TRAVEL SOCKS 146/T | 186/T

for women & men

15-20 mmHg

Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation
- Box cover promotes Traveno as a Travel friendly sock

Fiber Content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free



Calf (C)
Closed toe

Product Code Example: 146CB99/T

Series	Style	Size	Colour	Travel
146	C	B	99	/T

SEA ISLAND COTTON 151/191

for women & men

15-20 mmHg


Also available in 20-30 mmHg, see page 18

Features and Benefits

- Knitted with the world's rarest cotton for a soft, luxurious feel
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free



Calf (C)
Closed toe

Featured colour
is Black

Product Code Example: 151CA99

Series	Style	Size	Colour
151	C	A	99

ALL-SEASON MERINO WOOL 152/192

for women & men

15-20 mmHg


Also available in 20-30 mmHg, see page 20

Features and Benefits

- High quality extra-fine Merino Wool makes this sock soft and breathable
- Moisture wicking capabilities and naturally odour resistant
- Thermoregulating for year-round wear and keeps feet dry
- Smooth, flat toe seam

Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free



Calf (C)
Closed toe

Featured
colour is
Black



Product Code Example: 152CB99

Series	Style	Size	Colour
152	C	B	99

LINEN 153/193

for women & men

15-20 mmHg

Also available in 20-30 mmHg, see page 22

Features and Benefits

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free



Calf (C)
Closed toe

Featured
colour is
Lavender



Product Code Example: 153CB207

Series	Style	Size	Colour
153	C	B	207

MICROFIBER SHADES 143/183

for women & men

15-20 mmHg



Also available in 20-30 mmHg, see page 30

Features and Benefits

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Excellent breathability and moisture management
- Glides on and off with ease
- Smooth, flat toe seam

Fiber Content

- 87% Nylon, 13% Spandex
- Latex free



Calf (C)
Closed toe

Featured colour is
Royal Blue Argyle



Product Code Example: 183CB48

Series	Style	Size	Colour
183	C	B	48

ATHLETIC RECOVERY SOCKS 401

for women & men


15-20 mmHg

Features and Benefits

- drirelease® yarns provide moisture wicking and odour control
- Comfortable and easy to put on and take off
- Helps reduce lactic acid buildup
- Mild cushioning in sole, heel, and toe area for added comfort
- Smooth, flat toe seam

Fiber Content

- 67% drirelease® (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free



Calf (C)
Closed toe



15-20 mmHg
Closed Toe Calf
Size M in Black

401CM99

Featured colour
is White

Product Code Example: 401CM00

Series	Style	Size	Colour
401	C	M	00

PERFORMANCE SLEEVES 412V

for women & men


20-30 mmHg

Features and Benefits

- Stabilizes and supports active muscles
- Ideal for running, triathlons, endurance sports, and gym workouts
- Fabric resists deterioration in salt, chlorine or fresh water
- Quick drying
- High-visibility reflective transfer

Fiber Content

- 75% Nylon, 25% Spandex
- Latex free



Calf Sleeve (V)
Open foot,
healthy legs only

* Purple is not
available in X-Large

Featured colour
is Purple

Product Code Example: 412VL58

Series	Style	Size	Colour
412	V	L	58

MERINO OUTDOOR SOCKS 421

for women & men

15-20 mmHg

Also available in 20-30 mmHg, see page 40

Features and Benefits

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for year-round wear
- Cushioned soles help prevent blisters
- Mesh flex zones for comfort and ventilation
- Smooth, flat toe seam
- Soft and comfortable
- Perfect for outdoor work environments and activities

Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free

CHARCOAL
12

OLIVE
31

Calf (C)
Closed toe



Featured colour
is Charcoal

Product Code Example: 421CS12

Series	Style	Size	Colour
421	C	S	12

EVERSOFT DIABETIC SOCKS 160

for women & men

8-15 mmHg

Features and Benefits

- Padded sole, heel, and toe for comfort and protection
- Smooth, flat toe seam
- Compression-free foot
- drirelease® yarns provide moisture wicking and odour control
- Suitable for individuals with limited dexterity or lower tolerance for compression
- Ideal for individuals with sensitive skin

Fiber Content

- 67% drirelease® (57% Polyester, 10% Cotton), 26% Nylon, 7% Spandex
- Latex free

WHITE
00

Calf (C)
Closed toe



Product Code Example: 160CS00

Series	Style	Size	Colour
160	C	S	00



STANDARD SIZING CHART			
Women's Shoe	5-7	7.5-9.5	10-12
Ankle Circumference	18-23 cm (7"-9")	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")
Men's Shoe	6.5-8.5	9-11	11.5-14
Ankle Circumference	19-25 cm (7.5"-10")	20-27 cm (8"-10.5")	21-29 cm (8.5"-11.5")
Size	A	B	C

UNISEX SIZING CHART



UNISEX SIZING CHART

Women's Shoe	<7.5	7.5–9.5	10–12	12.5–15	≥15.5
Men's Shoe	<6.5	6.5–8.5	9–11	11.5–14	14.5–17
Ankle Circumference	18–24 cm (7"–9.5")	19–26 cm (7.5"–10")	20–27 cm (8"–10.5")	21–29 cm (8.5"–11.5")	21–29 cm (8.5"–11.5")
Size	S	M	L	X	2*

*Size 2 is only available for 401 Athletic Recovery Socks

PERFORMANCE SLEEVE 412V SIZING CHART

Ankle Circumference	18–21.5 cm (7"–8.5")	21.5–25.5 cm (8.5"–10")	25.5–29 cm (10"–11.5")	29–33 cm (11.5"–13")
Calf Circumference	28–39 cm (11"–15.5")	31.5–45 cm (12.5"–17.5")	35.5–51 cm (14"–20")	40.5–58 cm (16"–23")
Size	S	M	L	X

Donning Aids and Product Care Accessories



Make Life Easier

Why choose **Accessories**?

- Aid in putting on and removing stockings
- Stocking cleaning to maximize usage
- Extra liners for added comfort

ACCESSORIES

Product	Page
Gloves	76
Cotton Underliner	77
Washing Solution	77
It Stays® Roll-On Adhesive	77
Sim-Slide®	78
Magnide On/Off®	78
Sigvaris Doff N' Donner® ROLLY®	79
Sigvaris Doff N' Donner® Cone	79
Melany®	79

GLOVES



Textile Gloves*

These high-quality gloves have a latex surface to provide a good grip when putting on and taking off compression stockings.

PACKAGE: 12 PER PACK

SIZE	PRODUCT CODE
X-Small	593R400XS
Small	593R400S
Medium	593R400M
Large	593R400L
X-Large	593R400X
ASST'D 4x XS/S/M	593R400A



Latex-Free Gloves

Gripping donning gloves recommended for use in a latex-free environment.

PACKAGE: 12 PER PACK

SIZE	PRODUCT CODE
Small	592R400S
Medium	592R400M
Large	592R400L
X-Large	592R400X

TEXTILE GLOVE SIZING CHART

Measure Length and Width of the Hand

Hand Length**	20.5–21 cm	21.5–22 cm	22.5–23 cm	23.5–24 cm	24.5–25 cm
Hand Width	7–7.5 cm	7.5–8 cm	8–8.5 cm	9–9.5 cm	10–10.5 cm
	X-Small	Small	Medium	Large	X-Large

LATEX-FREE GLOVE SIZING CHART

Measure Circumference of the Hand at the Palm

Hand Circumference	16.5–18 cm (6.5"–7")	19–20.5 cm (7.5"–8")	21.5–23 cm (8.5"–9")	24–25.5 cm (9.5"–10")
	Small	Medium	Large	X-Large

* **Caution:** Textile Gloves contain natural rubber latex that may cause allergic reaction.

** **Hand Length:** measured from the tip of the longest finger to the crease under the palm.

LINERS



Cotton Underliner

10–15 mmHg

Can be used as an underlayer to protect fragile skin, secure wound dressings, add 10–15 mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings.

For sizing, please refer to the Natural Rubber Calf sizing on page 47.

PACKAGE: 2 PER PACK

SIZE	PRODUCT CODE	NATURAL RUBBER SIZING
Small / Short	588USS	S1/S3
Small / Long	588USL	S2/S4
Medium / Short	588UMS	M1/M3
Medium / Long	588UML	M2/M4
Large-XL / Short	588ULS	L1/L3 - X1/X3
Large-XL / Long	588ULL	L2/L4 - X2/X4

WASHING SOLUTION

Washing Solution

Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!

PACKAGE OPTIONS

CONTAINS	PRODUCT CODE
Six 16 oz. Bottles (473ml)	586W900
One 16 oz. Bottle (473ml)	586W950
Twelve 2 oz. Bottles (60ml)	586W975



IT STAYS!®

Roll-on body adhesive

A gentle roll-on body adhesive that washes off with water and may be used to hold compression garments in place against the skin.

IT STAYS!® is pliable and moves with the body.

SIZE	PRODUCT CODE
One 2 oz Bottle (60ml)	581I350
Twelve 2 oz. Bottles (60ml)	581I300



DONNING AIDS



Sim-Slide®

(For open toe Styles)

SIZE	PRODUCT CODE
Small	587SLIDESM
Medium	587SLIDEMD
Large	587SLIDELG
X-Large	587SLIDEXL

SIZING CHART For Sim-Slide®

Women's Shoe Size	≤4-6	6.5-8	8.5-9.5	≥10
Men's Shoe Size	≤3-4.5	5-6.5	7-8	≥8.5
Size	S	M	L	XL



See how to use
the Sim-Slide®
bit.ly/SIGssyt



Magnide On/Off®

(For closed toe Styles)

SIZE	PRODUCT CODE
Medium	587MAGMD
Large	587MAGLG
X-Large	587MAGXL

SIZING CHART For Magnide On/Off®

Women's Shoe Size	≤8.5	9-12	≥ 12.5
Men's Shoe Size	≤7.5	8-11	≥ 11.5
Size	M	L	XL

* XL suitable for extra-wide feet



See how to use
the Magnide On/Off®
bit.ly/49sxcYX

DONNING AIDS



SIGVARIS **Doff N' Donner®** ROLLY®

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both self-donning and caregiver use. Latex free.

SIZE	PRODUCT CODE
One Size	583D100

Improved design coming in 2025: Product Code **358496**



See how to use Sigvaris
DOFF N' DONNER® ROLLY®
bit.ly/3QSBcKU



SIGVARIS **Doff N' Donner®** Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner® ROLLY®. The bottom of the base firmly attaches to most clean, flat surfaces via suction.

SIZE	PRODUCT CODE
One Size	583D200

New product code in 2025: **358497**



MELANY®

MELANY® from Sigvaris makes both donning and doffing compression garments simple. Its innovative one-piece design with adjustable handles was made especially for users with limited range of motion and finger dexterity.*

SIZE	CALF CIRCUMFERENCE	PRODUCT CODE
X-Small	25–30 cm	579SXS
Small	30–38 cm	579SSM
Medium	38–46 cm	579SMD
Large	46–56 cm	579SLG

* Individual components not sold separately



See how to use
the **MELANY®**
bit.ly/3SCQYLu

Learn more



Help yourself

How can **resources** help?

- Learn about the levels of compression for various indications to better help your clients
- Understand the terms and icons that we use
- Follow specific instructions on how to measure for the perfect fit

RESOURCES

Section	Page
Symptoms & Indications	82
Glossary of Icons	83
Glossary of Terms	84
How to Measure	85
How to Order	86

SYMPTOMS / INDICATIONS	COMPRESSION LEVEL AT THE ANKLE
<ul style="list-style-type: none"> □ Tired, achy legs □ Occupational or evening edema □ Leg discomfort from long hours of standing or sitting □ Prevent varicosities & swelling during pregnancy □ Preventative for legs predisposed to risk □ Reduce swelling during travel □ Prevent DVT during travel 	15–20 mmHg
<ul style="list-style-type: none"> □ Heavy, fatigued, aching legs □ Mild edema in lower extremities □ Mild varicosities with minimal edema □ Mild varicosities during pregnancy □ Lymphedema upper extremities □ Post-procedure of small veins □ Prevention or management of DVT <p><i>Note: A prescription is recommended for stockings with the indications of 20 mmHg or higher</i></p>	20–30 mmHg
<ul style="list-style-type: none"> □ Moderate varicosities with mild to moderate edema (also during pregnancy) □ Post-procedure of larger veins to maintain treatment success □ CEAP C3 (edema), C4 (skin changes without ulceration), C5 (skin changes with healed ulceration), C6 (skin changes with active ulcerations) □ Management and treatment of DVT or Post-Thrombotic Syndrome □ Superficial phlebitis □ Lymphedema after decongestant therapy to maintain reduction □ Orthostatic / postural hypotension □ Post-fracture, post-traumatic edema 	30–40 mmHg
<ul style="list-style-type: none"> □ Lymphedema after decongestant therapy to maintain reduction □ Severe Post-Thrombotic Syndrome □ Severe varicosities and/or edema 	40 mmHg or higher
<ul style="list-style-type: none"> □ Severe arterial insufficiency, intermittent claudication, ischemia □ Uncontrolled congestive heart failure 	Contraindications
<ul style="list-style-type: none"> □ Skin sensitivities or allergies □ Advanced neuropathy with sensory loss in leg or foot □ Diabetes with arterial compromise □ Confinement to bed or non-ambulatory use unless otherwise prescribed by the physician □ Eczema, hypodermatitis, dermatosclerosis, dermatitis 	Cautions
<p>The therapeutic effect of graduated compression stockings is based on the compression level. With the proper CEAP classification, the correct compression level at the ankle can be chosen. For more severe symptoms, consult with a vascular specialist.</p>	
<p><i>No liability accepted for non-observance of contraindications and cautions.</i></p> <p>Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.</p>	



Sensitive Skin



Thermoregulating



Breathable



Moisture Wicking



Soft




Durable



Cotton




Wool



Microfiber



Natural Rubber



Largest Size Range



Sheer



Medium Sheer



Opaque



Stylish



Stylish



Elegant



Colours



Patterns



Sensinno® Grip-top
Band



Anatomically
Contoured



Cushioned Sole



Anti-bacterial



Hypoallergenic
Natural fibers



Moisture
Management

TERM	DEFINITION	SIGVARIS QUALITY STANDARD
CHRONIC VENOUS DISEASE (CVD)	A disease state comprised of a wide range of disorders including tired, achy legs, varicose veins, Chronic Venous Disease (CVD), skin issues, and venous leg ulcers (open and healed).	SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH CARE products to help reduce the risk of, improve symptoms of, and provide treatment for CVD.
DOUBLE-COVERED YARN	Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort.	Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort.
DEEP VEIN THROMBOSIS (DVT)	A blood clot in a deep vein, most commonly in the lower leg or thigh.	SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening.
GRADUATED COMPRESSION VS. ANTI-EMBOLISM	Graduated compression is for ambulatory patients, whereas anti-embolism stockings are for non-ambulatory patients.	SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation.
HYPOALLERGENIC	Designed to minimize the chance of an allergic reaction.	SIGVARIS GROUP offers Latex free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin.
ISO CERTIFIED	SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have the required quality systems in place.	SIGVARIS GROUP follows ISO Standard 13485:2016 for our products manufactured in North and South America.
MICROFIBER	An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural fiber fabric. This fabric has excellent breathability.	Used to create products with a soft and comfortable feel. Find this fiber in Soft Opaque 840, Microfiber Patterns 830, Microfiber Shades 143/183, and Microfiber 820.
mmHg (MILLIMETERS OF MERCURY)	Defined as the extra pressure generated by a column of mercury one millimeter high; the same unit of measurement used for blood pressure.	SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement.
MERINO WOOL	Merino is the softest and finest wool available in the world and is inherently itch-free.	SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in Thermoregulating Wool 422, Merino Outdoor Socks 421, Merino Wool 240, and All-Season Merino Wool 152/192.
SEA ISLAND COTTON	Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness.	SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our Sea Island 220 and Sea Island Cotton 151/191.
SENSINNOV®	A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day.	This top band is featured in Sheer 780, Soft Opaque 840, and Microfiber 820 to help keep thigh-high stockings in place all day long.
SUPIMA COTTON	Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber.	SIGVARIS GROUP uses this premium cotton yarn in the Cotton 230, Sea Island 220, Sea Island Cotton 151/191, Casual Cotton 146/186, and Cushioned Cotton 360, making them comfortable and ideal for those with sensitive skin.

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as early in the morning as possible when swelling is minimal.

CALF STYLE

Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Determine the calf length by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg.



1 Ankle circumference directly above ankle bone



2 Calf circumference at fullest part of the calf


PANTYHOSE AND THIGH-HIGH STYLE

Open/Closed Toe

1	Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM or LARGE. Note: This should be your first measurement.
2	Measure the calf circumference at the fullest part of the calf.
3	Skip Step 3 (but refer to General Guidelines #6 below).
4	Measure the thigh at its widest circumference.
5	Determine the leg length by measuring from the floor up to the gluteal fold.
6	Determine the hip circumference by measuring the widest part of the hips.




3 Calf length from the floor to the right-angle bend of the knee



4 Widest circumference of the thigh

GENERAL GUIDELINES

1	Measurements should always be taken as early in the morning as possible, when swelling is minimal.
2	Measure against bare skin, not over clothing (Steps 1–5).
3	The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to right for examples).
4	Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting.
5	Length measurements should be taken from the ground up.
6	Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles.
7	Record shoe size; it may be helpful in determining the best possible fit.



5 Leg length from the floor up to the gluteal fold



6 Circumference at the hip

HOW TO ORDER MEDICAL AND SPECIALTY PRODUCTS



Product Code Example: 232CMLM99 (page 8)

Series Number	Compression Level	Product Style	Size	Gender	Colour
23 = Cotton <i>(See individual pages, use first 2 digits of series number)</i>	1 = 15–20 mmHg 2 = 20–30 mmHg 3 = 30–40 mmHg 4 = 40–50 mmHg* 5 = 50–60 mmHg*	C = Knee-high or Calf <i>(See individual pages for available styles)</i>	ML = Medium Long <i>(See individual pages for sizing)</i>	M = Men W = Women O = Open toe	99 = Black <i>(See individual pages for available colours)</i>
23	2	C	ML	M	99
* 40–50 mmHg and 50–60 mmHg available in Natural Rubber only					
To order Plus Calf with Grip-top or Plus Size Pantyhose (860P for women only): add “/P” after the colour in the ordering code					
To order Calf Style with Grip-top: add “/S” after the colour in the ordering code <i>(see product pages for availability)</i>					
To order Thigh with Waist Attachment: add “/L” to indicate the Left leg or “/R” to indicate the Right leg after the colour in the ordering code <i>(860W and 500W only)</i>					
To order Motion High Tech, Motion Thermoregulating Wool and Motion Flow Tech®: The letter for Gender is not used in the ordering code. Example: 412CSM59 = Motion High Tech, 20–30 mmHg, Closed toe Calf, Size Small Medium, Colour Red					

HOW TO ORDER CONSUMER HEALTH PRODUCTS



Product Code Examples: 151CA99 (page 68) 191CB99 (page 68)

Series Number	Product Style	Size	Colour
151/191 = SEA ISLAND COTTON First 3 digits indicate gender 151 = woman 191 = man (See individual pages for series number)	C = Calf (See individual pages for available styles)	According to shoe size A= women's shoe sizes 5-7 B = men's shoe sizes 9-11 (See individual pages for sizing)	99 = Black (See individual pages for available colours)
151	C	A	99
191	C	B	99
ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, MERINO OUTDOOR SOCKS, and EVERSOFT DIABETIC SOCKS are UNISEX: the first 3 digits are the same for both genders. Example: 421CL31 = MERINO OUTDOOR SOCK, 15–20 mmHg, Closed toe Calf, Size Large, Colour Olive			

Sigvaris, SIGVARIS GROUP, PFS2, Sensinnov, Compreflex, Coolflex, Flow Tech, Compreboot, Accutab, Doff N' Donner, ROLLY and MELANY are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland, in many countries worldwide. © Copyright 2025 SIGVARIS, Inc.

Spacer Fabric and MeshLink are trademarks of SIGVARIS AG, CH-9014, St. Gallen/Switzerland, in many countries worldwide.

Lycra is a registered trademark of INVISTA. Supima is a registered trademark of the Supima Association of America Corporation. SENSIL is a registered trademark of Nilit. All rights reserved. drirelease is copyright ©2021 drirelease. drirelease is an Optimer Brands company. All rights reserved. Sim-Slide and Magnide are registered trademarks of Arion Holding B.V. IT STAYS is a registered trademark of New Horizons Diversified Interests.

mobilis

by SIGVARIS GROUP

SOFT ORTHOPEDIC SUPPORTS AND BRACES

INDEX

Intro 88

Mode of action 91

How to measure 91

Elbow 92

EPICARE - Elbow Support 93

EPIBAND - Elbow Strap 93

EPIACTIVE - Elbow Brace 94

Wrist and thumb 95

MANUCARE - Wrist Support 96

MANUWRAP - Wrist and Thumb Brace 96

MANUACTIVE - Wrist Brace 97

MANUSUPPORT - Wrist Brace 97

MANUPLUS SUPPORT - Wrist and Thumb Stabilizer 98

RHIZOSUPPORT - Thumb Brace 98

Back 99

LUMBOCARE - Back Support 100

LUMBOSUPPORT - Lumbar Brace 100

LUMBOSUPPORT EVO - Lumbar Brace 101

LUMBOINTENSE - Lumbar Brace 101

Knee 102

GENUCARE - Knee Support 103

GENUWRAP - Knee Wrap 103

GENUBAND - Knee Strap 104

GENUACTIVE STAB - Stabilizing Knee Brace 104

GENUACTIVE PAD - Patellar Knee Brace 105

Ankle 106

MALLEOCARE - Ankle Support 107

MALLEOACTIVE PAD - Ankle Brace 107

MALLEOACTIVE BAND - Ankle Brace 108

MALLEOSUPPORT - Ankle Brace 108

MALLEOINTENSE - Ankle Stabilizer 109

How to order Mobilis 109

MOBILIS: SOFT ORTHOPEDIC SUPPORTS AND BRACES

A collection of high-quality supports and braces from SIGVARIS GROUP designed to provide support, stability, and pain relief while allowing individuals to stay active.




Mobilis product names are made of two parts:

① Body part where the product is used


	LINE
Elbow	→ EPI
Hand	→ MANU
Thumb	→ RHIZO
Back	→ LUMBO
Knee	→ GENU
Ankle	→ MALLEO

② Support level:




Example ① + ② = **MALLEOSUPPORT**

MODE OF ACTION




Compression effect

- Decreases swelling and pressure on joints
- Provides pain relief
- Improves mobility



Massage effect with padding

- Targeted swelling reduction
- Improves blood circulation
- Aids healing



Joint stabilization

- Rigid elements stabilize the joint and improve muscle activity



Proprioception


- Aids in perception of movements and body position



Warming effect


- Increases comfort and may provide pain relief

HOW TO MEASURE



Elbow

Measure the arm circumference at 3 cm below the elbow (two fingers)




Measure the elbow circumference with 30 degrees of inclination

Wrist and thumb




Measure wrist circumference

Back




Measure waist circumference at belly button

Knee



Measure knee circumference with 30 degrees of inclination



Measure thigh circumference at 15 cm above the knee

Ankle



Measure ankle circumference above ankle bone

Elbow



	EPI CARE	EPI BAND	EPI ACTIVE
INDICATIONS			
Arthritis	○		○
Joint fluid and swelling	○		○
Chronic irritation/ Tendomyopathy	○		○
Post trauma, Post surgery	○		○
Tennis/golfer's elbow (Epicondylopathy)	○	○	○
SUPPORT LEVEL	 Support level mild	 Support level medium	 Support level medium

EPICARE

Elbow Support




Contains natural rubber latex that may cause allergic reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, tennis/golfer's elbow (epicondylopathy)

ACTION MODE




FEATURES

- Mild support
- Maximum comfort
- Discreet look
- Can be worn on either elbow
- Optimal adhesion due to silicone band

COMPOSITION

64% Polyamide, 27% Elastodiene, 5% Elastane, 4% Silicone



Size	Circumference	SKU
Size XS	20-23 cm	MECE-XS-66
Size S	23-26 cm	MECE-SM-66
Size M	26-29 cm	MECE-MD-66
Size L	29-32 cm	MECE-LG-66
Size XL	32-35 cm	MECE-XL-66

EPIBAND

Elbow Strap



INDICATIONS

Tennis/golfer's elbow (epicondylopathy)

ACTION MODE



FEATURES

- Adjustable strap to modify support level
- Silicone pad provides relief for tendon
- Can be worn on either elbow

COMPOSITION

89% Polyamide, 11% Elastane



Measure the arm circumference at 3 cm below the elbow (two fingers)



Size	Circumference	SKU
One size	20 - 35 cm	MEBE-91

EPIACTIVE

Elbow Brace



Support level
medium



Contains natural
rubber latex that
may cause allergic
reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, tennis/golfer's elbow (epicondylopathy)

ACTION MODE



Compression



Massaging effect



Proprioception

FEATURES

- Elbow brace with silicone inserts
- Removable elastic strap to adjust support
- Anatomical shape and 3D knitting for maximum comfort
- Can be worn on either elbow

COMPOSITION

56% Thermoplastic elastomer, 32% Polyamide, 7% Elastodiene, 3% Elastane, 2% Polyester



Measure the elbow
circumference with
30 degrees of inclination



Size	Circumference	SKU
Size XS	20 - 23 cm	MEAE-XS-91
Size S	23 - 26 cm	MEAE-SM-91
Size M	26 - 29 cm	MEAE-MD-91
Size L	29 - 32 cm	MEAE-LG-91
Size XL	32 - 35 cm	MEAE-XL-91



Watch a
product video



Wrist & Thumb



	MANU CARE	MANU WRAP	MANU ACTIVE	MANU SUPPORT	MANU PLUS SUPPORT	RHIZO SUPPORT
INDICATIONS						
Arthritis	○	○	○	○	○	○
Chronic irritation/ Tendomyopathy	○	○	○	○	○	○
Carpal tunnel syndrome			○	○	○	
Instability			○	○	○	○
Post trauma, Post surgery		○	○	○	○	
Rhizarthrosis					○	○
Skier's thumb/Torn thumb ligament (Lateral metacarpo- phalangeal ligament lesion)					○	○
SUPPORT LEVEL	 Support level mild	 Support level moderate	 Support level medium	 Support level strong	 Support level strong	 Support level strong


MANUCARE

Wrist Support

INDICATIONS

Arthritis, chronic irritation/tendomyopathy

ACTION MODE



FEATURES

- Mild support
- Soft compression
- Ultimate comfort
- Discreet look
- Can be worn on either wrist

COMPOSITION

66% Polyamide, 29% Elastodiene, 5% Elastane



Contains natural rubber latex that may cause allergic reaction



Size	Circumference	SKU
Size XS	13 - 15 cm	MMCW-XS-66
Size S	15 - 17 cm	MMCW-SM-66
Size M	17 - 19 cm	MMCW-MD-66
Size L	19 - 21 cm	MMCW-LG-66
Size XL	21 - 23 cm	MMCW-XL-66


MANUWRAP

Wrist and Thumb Brace

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery

ACTION MODE



FEATURES

- Moderate adjustable support level
- Breathable fabric
- Can be worn on either hand

COMPOSITION

96% Polyamide, 4% Elastane



Size	Circumference	SKU
Sizes S/M	13 - 17 cm	MMWW-SM/MD-91
Sizes L/XL	17 - 23 cm	MMWW-LG/XL-91


MANUACTIVE

Wrist Brace

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery

ACTION MODE



FEATURES

- Brace with silicone inserts around cubitus bone
- Metal rod to limit wrist movement
- Elastic strap allows support level adjustment
- Anatomically shaped 3D knit for comfort and fit
- Specific for left or right hand

COMPOSITION

56% Polyamide, 25% Thermoplastic Elastomer, 9% Elastodiene, 5% Polyester, 5% Elastane



Measure wrist circumference



Contains natural rubber latex that may cause allergic reaction



Watch a product video



Size	Circumference	SKU Left Hand	SKU Right Hand
Size XS	13 - 15 cm	MMAW-XS-L-91	MMAW-XS-R-91
Size S	15 - 17 cm	MMAW-SM-L-91	MMAW-SM-R-91
Size M	17 - 19 cm	MMAW-MD-L-91	MMAW-MD-R-91
Size L	19 - 21 cm	MMAW-LG-L-91	MMAW-LG-R-91
Size XL	21 - 23 cm	MMAW-XL-L-91	MMAW-XL-R-91


MANUSUPPORT

Wrist Brace

INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery

ACTION MODE



FEATURES


- Immobilization/stabilization brace: rigid palm structure, external and internal rods to restrain wrist movements
- Three hook-and-loop straps to adjust the fit
- Can be worn on either wrist

COMPOSITION

92% Polyamide, 5% Elastane, 3% Cotton



Watch a product video



Measure wrist circumference

Size	Circumference	SKU
Sizes S/M	13 - 17 cm	MMSWS-SM/MD-91
Sizes L/XL	17 - 23 cm	MMSWS-LG/XL-91

MANUPLUS SUPPORT


Wrist & Thumb Stabilizer



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, carpal tunnel syndrome, instability, post trauma, post surgery, rhizarthrosis, skier's thumb/torn thumb ligament (lateral metacarpophalangeal ligament lesion)

ACTION MODE




FEATURES

- Immobilization/stabilization brace for wrist and thumb
- The two hook-and-loop straps can be clamped by laces to adjust the fit
- Specific for right or left hand

COMPOSITION

50% Polyester, 46% Polyamide, 4% Polyurethane



Measure wrist circumference



Size	Circumference	SKU Left Hand	SKU Right Hand
Size XS	13 - 15 cm	MMPSW-XS-L-91	MMPSW-XS-R-91
Size S	15 - 17 cm	MMPSW-SM-L-91	MMPSW-SM-R-91
Size M	17 - 19 cm	MMPSW-MD-L-91	MMPSW-MD-R-91
Size L	19 - 21 cm	MMPSW-LG-L-91	MMPSW-LG-R-91
Size XL	21 - 23 cm	MMPSW-XL-L-91	MMPSW-XL-R-91

RHIZOSUPPORT


Thumb Brace



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, instability, rhizarthrosis, skier's thumb/torn thumb ligament (lateral metacarpophalangeal ligament lesion)

ACTION MODE



FEATURES

- Immobilization/stabilization brace for thumb due to external and internal rods
- Two hook-and-loop straps to adjust the fit
- Can be worn on either hand

COMPOSITION

92% Polyamide, 8% Elastane




Watch a product video




Measure wrist circumference

Size	Circumference	SKU
Sizes S/M	13 - 17 cm	MRST-SM/MD-91
Sizes L/XL	17 - 23 cm	MRST-LG/XL-91



Back



	LUMBO CARE	LUMBO SUPPORT	LUMBO SUPPORT EVO	LUMBO INTENSE
INDICATIONS				
Muscular tension	○	○	○	○
Lower back pain (Lumbar syndrome)	○	○	○	○
Muscular weakness and loss (degeneration)		○	○	○
Post surgery		○	○	○
Sciatica		○	○	○
SUPPORT LEVEL	Support level mild	Support level strong	Support level strong	Support level extra strong

LUMBOCARE


Back Support



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome)

ACTION MODE




FEATURES

- Adjustable belt provides mild support
- Comfortable material that provides a warming effect
- Can be worn against the skin

COMPOSITION

65% Cotton, 31% Polyester, 4% Elastodiene



Size	Circumference	SKU
Size S	70 - 95 cm	MLCB-SM-66
Size M	95 - 115 cm	MLCB-MD-66
Size L	115 - 135 cm	MLCB-LG-66

LUMBOSUPPORT


Lumbar Brace



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



FEATURES


- Four stabilizing rods to support and secure
- Adjustable straps for a precise fit
- To be worn over clothing or underwear
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

54% Polyamide, 21% Elastane, 20% Polypropylene, 5% Polyester



Watch a product video



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size S	70 - 95 cm	MLSB-SM-S-91	MLSB-SM-T-91
Size M	95 - 115 cm	MLSB-MD-S-91	MLSB-MD-T-91
Size L	115 - 135 cm	MLSB-LG-S-91	MLSB-LG-T-91

LUMBOSUPPORT EVO


Lumbar Brace



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



Watch a product video

FEATURES


- Lumbar spine stabilizing belt with four posterior rods
- The double elastic hook-and-loop harness is removable and permits individual adjustment
- To be worn over the clothes
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less;
Tall (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

44% Polyamide, 23% Polyester, 8% Viscose, 7% Elastodiene, 6% Ethylene vinyl acetate, 5% Polyvinyl chloride, 4% Polyurethane, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size S	56 - 71 cm	MLSEB-SM-S-99	MLSEB-SM-T-99
Size M	68 - 84 cm	MLSEB-MD-S-99	MLSEB-MD-T-99
Size L	81 - 96 cm	MLSEB-LG-S-99	MLSEB-LG-T-99
Size XL	94 - 109 cm	MLSEB-XL-S-99	MLSEB-XL-T-99
Size XXL	107 - 122 cm	MLSEB-2X-S-99	MLSEB-2X-T-99

LUMBOINTENSE


Lumbar Brace



INDICATIONS

Muscular tension, lower back pain (lumbar syndrome), muscular weakness and loss (degeneration), post surgery, sciatica

ACTION MODE



Watch a product video

FEATURES


- Stabilizing belt with four posterior rods provides stronger pressure on the lumbar spine
- The double elastic hook-and-loop harness is removable and permits individual adjustment
- To be worn over the clothes
- Two different heights: **Short** (21 cm): Recommended for individuals with a height of 5'3" (160 cm) or less; **Tall** (26 cm): Recommended for individuals with a height of 5'4" (162 cm) or more.

COMPOSITION

47% Polyamide, 27% Polyester, 2% Viscose, 4% Elastodiene, 6% Ethylene-acetate, 5% Polyvinyl chloride, 4% Polyurethane, 2% Polypropylene, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



Size	Circumference	SKU Short (S)	SKU Tall (T)
Size S	56 - 71 cm	MLIB-SM-S-99	MLIB-SM-T-99
Size M	68 - 84 cm	MLIB-MD-S-99	MLIB-MD-T-99
Size L	81 - 96 cm	MLIB-LG-S-99	MLIB-LG-T-99
Size XL	94 - 109 cm	MLIB-XL-S-99	MLIB-XL-T-99
Size XXL	107 - 122 cm	MLIB-2X-S-99	MLIB-2X-T-99

Knee



	GENU CARE	GENU WRAP	GENU BAND	GENU ACTIVE STAB	GENU ACTIVE PAD
INDICATIONS					
Arthritis	○	○		○	○
Joint fluid and swelling	○			○	○
Chronic irritation/ Tendomyopathy	○	○	○	○	○
Post trauma, Post surgery	○	○		○	○
Kneecap (patellar) pain		○	○	○	○
Deep knee pain (Femoro-patellar syndrome)		○	○	○	○
Instability, Hyperlaxity, Meniscus lesion, Gonarthrosis				○	
SUPPORT LEVEL	 Support level mild	 Support level moderate	 Support level medium	 Support level medium	 Support level medium

GENUCARE

Knee Support



Contains natural rubber latex that may cause allergic reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery

ACTION MODE




FEATURES

- Mild support
- Optimal adhesion due to silicone band
- Discreet look
- Can be worn on either knee

COMPOSITION

64% Polyamide, 28% Elastodiene
5% Elastane, 3% Silicone



Size	Circumference	SKU
Size XS	30 - 34 cm	MGCK-XS-66
Size S	34 - 38 cm	MGCK-SM-66
Size M	38 - 42 cm	MGCK-MD-66
Size L	42 - 46 cm	MGCK-LG-66
Size XL	46 - 50 cm	MGCK-XL-66

GENUWRAP

Knee Wrap



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE




FEATURES

- Moderate support knee wrap with three straps for optimum fit
- Adjustable size
- Breathe-O-Prene material allows skin to breathe
- Can be worn on either knee

COMPOSITION

94% Polyamide, 6% Elastane



Size	Circumference	SKU
Size S/M	30 - 42 cm	MGWK-SM/MD-91
Size L/XL	42 - 55 cm	MGWK-LG/XL-91

GENUBAND

Knee Strap



INDICATIONS

chronic irritation/tendomyopathy, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE




FEATURES

- Adjustable knee strap to modify support level
- Silicone insert is located on kneecap tendon
- Can be worn on either knee
- Anatomical shape for extra comfort

COMPOSITION

89% Polyamide, 11% Elastane



Size	Circumference	SKU
One Size	30 - 55 cm	MGBK-91

GENUACTIVE STAB


Stabilizing Knee Brace



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome), instability, hyperlaxity, meniscus lesion, gonarthrosis

ACTION MODE




FEATURES

- Stabilizing knee brace with lateral supports
- Silicone insert around the kneecap with three Velcro straps
- The anatomically shaped fabric is compressive
- Wearable on both knees

COMPOSITION

8% Polyamide, 23% Elastomer, 15% Elastodiene, 6% Polyester, 5% Elastane, 2% Polyvinyl chloride, 1% Silicone



Watch a product video



Size	Knee Circumference	Thigh Circumference	SKU
Size XS	30 - 34 cm	40 - 44 cm	MGSK-XS-91
Size S	34 - 38 cm	44 - 48 cm	MGSK-SM-91
Size M	38 - 42 cm	48 - 52 cm	MGSK-MD-91
Size L	42 - 46 cm	52 - 56 cm	MGSK-LG-91
Size XL	46 - 50 cm	56 - 60 cm	MGSK-XL-91
Size 2X	50 - 54 cm	60 - 64 cm	MGSK-2X-91

GENUACTIVE PAD

Patellar Knee Brace



Support level
medium




Contains natural
rubber latex that
may cause allergic
reaction

INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/ tendomyopathy, post trauma, post surgery, kneecap (patellar) pain, deep knee pain (femoro-patellar syndrome)

ACTION MODE




FEATURES

- Stabilizing knee brace with lateral supports
- Silicone insert around the kneecap
- Anatomically shaped 3D knit for comfort and fit

COMPOSITION

52% Polyamide, 26% Thermoplastic Urethane, 14% Elastodiene, 5% Polyester, 2% Elastane, 1% Silicone



Measure knee circumference below the kneecap with 30 degrees of inclination, and thigh circumference 15 cm above the knee



Watch a
product video




Size	Knee Circumference	Thigh Circumference	SKU
Size XS	30 - 34 cm	40 - 44 cm	MGAK-XS-91
Size S	34 - 38 cm	44 - 48 cm	MGAK-SM-91
Size M	38 - 42 cm	48 - 52 cm	MGAK-MD-91
Size L	42 - 46 cm	52 - 56 cm	MGAK-LG-91
Size XL	46 - 50 cm	56 - 60 cm	MGAK-XL-91
Size 2X	50 - 54 cm	60 - 64 cm	MGAK-2X-91

3D Knit
Provides breathability and comfort for both everyday and physical activities.

Lateral Supports
Provides high stability

Silicone Insert
The anatomical shape helps the kneecap to remain stable

Bend Zone
Knitted fabric is less compressive for more comfort and flexibility





Ankle

	MALLEO CARE	MALLEO ACTIVE PAD	MALLEO ACTIVE BAND	MALLEO SUPPORT	MALLEO INTENSE
INDICATIONS					
Arthritis	○	○	○	○	○
Joint fluid and swelling	○	○	○		
Chronic irritation/ Tendomyopathy	○	○	○	○	○
Post trauma, Post surgery, Soft tissue lesion	○			○	○
Hyperlaxity				○	○
Chronic instability				○	○
Post trauma, Post surgery, Ligament injury				○	○
Moderate sprains				○	○
SUPPORT LEVEL	 Support level mild	 Support level medium	 Support level medium	 Support level strong	 Support level extra strong

MALLEOCARE

Ankle Support




Contains natural rubber latex that may cause allergic reaction



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/ tendomyopathy, post trauma, post surgery, soft tissue lesion

ACTION MODE



FEATURES

- Mild support
- Adjustable level of stability
- Can be worn on either ankle

COMPOSITION

70% Polyamide, 12% Elastane,
9% Elastodiene, 8% Polyester,
1% Polypropylene



Size	Circumference	SKU
Size S	18 - 22 cm	MMCA-SM-66
Size M	22 - 25 cm	MMCA-MD-66
Size L	25 - 28 cm	MMCA-LG-66

MALLEOACTIVE PAD

Ankle Brace




Contains natural rubber latex that may cause allergic reaction



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/ tendomyopathy

ACTION MODE




FEATURES

- 2 lateral pads stabilize and support the ankle
- Gentle massaging effect
- Breathable 3D knit for optimal fit and comfort
- Left and right foot option
- Includes a foot slip for easier donning

COMPOSITION

45% Polyamide, 40% Thermoplastic Urethane, 10% Elastodiene, 3% Elastane, 2% Polyester



Watch a product video

Size	Circumference	SKU Left Ankle	SKU Right Ankle
Size S	16 - 21 cm	MMAA-SM-L-91	MMAA-SM-R-91
Size M	21 - 26 cm	MMAA-MD-L-91	MMAA-MD-R-91
Size L	26 - 31 cm	MMAA-LG-L-91	MMAA-LG-R-91
Size XL	31 - 37 cm	MMAA-XL-L-91	MMAA-XL-R-91

MALLEOACTIVE BAND


Ankle Brace



INDICATIONS

Arthritis, joint fluid and swelling, chronic irritation/tendomyopathy

ACTION MODE



FEATURES

- The elastic band allows for adjustable support
- The anatomically shaped breathable knit fabric is compressive
- Can be worn in a shoe
- Can be worn on either ankle

COMPOSITION


62% Polyamide, 18% Elastodiene, 12% Viscose, 5% Polyester, 3% Elastane



Contains natural rubber latex that may cause allergic reaction



Watch a product video



Size	Circumference	SKU
Size XS	18 - 20 cm	MMBA-XS-91
Size S	20 - 23 cm	MMBA-SM-91
Size M	23 - 25 cm	MMBA-MD-91
Size L	25 - 27 cm	MMBA-LG-91
Size XL	27 - 30 cm	MMBA-XL-91

MALLEOSUPPORT


Ankle Brace



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion, hyperlaxity, chronic instability, ligament injury, moderate sprains

ACTION MODE




FEATURES

- Stabilizing ankle brace with lateral rigid rods
- Two hook-and-loop straps for added support
- The lace closure system allows for an adjustable fit
- Can be worn on either ankle

COMPOSITION

91% Polyamide, 5% Elastane, 2% Polyester, 2% Cotton



Size	Circumference	SKU
Size S	19 - 22 cm	MMSA-SM-91
Size M	22 - 26 cm	MMSA-MD-91
Size L	26 - 30 cm	MMSA-LG-91



MALLEOINTENSE

Ankle Stabilizer



Contains natural rubber latex that may cause allergic reaction




Watch a product video



INDICATIONS

Arthritis, chronic irritation/tendomyopathy, post trauma, post surgery, soft tissue lesion, hyperlaxity, chronic instability, ligament injury, moderate to severe sprains

ACTION MODE




FEATURES

- Immobilizing ankle stabilizer based on two rigid shells
- Adaptable fit with two gel cushions to provide comfort and cooling effect
- Two Velcro straps as well as heel support strap provide ability to adjust tightness and support
- Can be worn on either ankle

COMPOSITION

50% Gel, 22% Polyvinyl chloride, 20% Polypropylene, 7% Polyamide, 1% Neoprene



Size	Height	SKU
One Size	≥ 150 cm	MMIA-91

HOW TO ORDER MOBILIS PRODUCTS

M	E	A	E	-	MD	-	-	-	91
Mobilis	Letter matches first part of the product name: E = Epi M = Manu MP = Manu Plus R = Rhizo L = Lumbo G = Genu M = Malleo	Letter matches line of product: C = Care W = Wrap A = Active S = Support or Stab B = Band I = Intense SE = Support Evo	Letter matches body part where the product is used: E = Elbow W = Wrist T = Thumb B = Back K = Knee A = Ankle		Size				Available colours: 66 = Beige 91 = Graphite 99 = Black

Add **L** to indicate the Left side or **R** to indicate the Right side after the Size in the ordering code. **L** or **R** is not specified in products that fit both left and right. **Example:** MMAW-MD-L-91.

For back braces, add **S** to indicate the Short height or **T** to indicate the Tall height after the Size in the ordering code.

You can find products without a size. **Example:** MEBE-91

You can find products with combined sizes. **Example:** MRST-SM/MD-91



SIGVARIS
GROUP

We are dedicated to **helping people feel their best.** Every day.



Get more done faster:

Chat with us live on the
Online Dealer Portal:
sigvaris-online.com



Download
Center

To order:

sigvaris-online.com
canadaorders@sigvaris.com
1.800.363.4999



sigvarisgroup.canada



sigvarisgroup_canada



[SIGVARIS GROUP Canada](#)



[SIGVARIS GROUP Canada](#)