
Compression Solutions

SIGVARIS

GROUP



Graduated Compression Socks, Hosiery and Inelastic Garments

SIGVARIS GROUP

Established in 1864, SIGVARIS GROUP is, and has always been, a family-owned company. Founded and headquartered in Switzerland, our passion for helping people feel their best is why we've been a global market leader for 158 years and counting.

Our compression garment portfolio includes Medical, Consumer Health Care and products for Medically Complex Edema that are specifically designed to restore, assist and enhance people's lives. We are renowned for making exceptional products that combine advanced technology and quality craftsmanship with diverse styles that fit our modern customers' lives.

When heritage, competence, and dedication come together, there's no limit to where we'll go.

Live your life. Be yourself.

For more information, please visit: sigvaris.com

The World-Wide Leader in Compression Therapy

WHY SIGVARIS GROUP?

We produce the best compression garments to enhance our customers' health and quality of life. We are dedicated to helping legs look and feel their best. Using premium yarns, we've spun ourselves into being the world-wide leader in graduated compression socks and hosiery by embracing our core beliefs of QUALITY, COMFORT, FIT and INNOVATION. Staying true to these values, we are continuing our efforts by expanding our business to help customers with medically-complex needs.

COMFORT

We offer the most comfortable fabric options More than 20 people handle our products from available, including the finest cottons, wools and microfibers. We double-cover the yarns for ease of highest quality. donning and doffing.

PRECISE FIT

Our products are designed using the most recent human data measurement studies, fitting 99.49% of the population.

VARIETY

Constructed with body-shaping contouring and a We offer you more choices of materials, colours, cotton lined gusset for breathability and comfort. patterns, and styles than any other medical compression brand.

GRIP TOPS

High stretch, high memory components to deliver comfortable "stay up" power.

INNOVATION

We use the finest textiles and materials to create Canada-wide salesforce and a live customer care modern styles in our circular-knit products, while team combine to provide dedicated support maintaining medical efficacy. We continue to & expertise. invest in medical research to help people feel their best, every day.

TO ORDER: 800.363.4999 (phone) 800.263.8736 (fax) canadaorders@sigvaris.com

Online Dealer Portal: sigvaris-online.com

QUALITY

start to finish, ensuring every garment is of the

NATURAL FIBERS

We have the highest use of natural materials, including Sea Island Cotton, Supima Cotton, Linen and Merino Wool.

PANTYHOSE

RECIPROCATED TOE

Many of our products are constructed with a smooth, flat toe seam and a contoured toe box for added comfort.

CUSTOMER SERVICE

| | PRODUCT | | STYLES | COMPRESSION LEVEL | PAGE |
|-----------|-------------------------|-----------------|---|--|--------|
| IAL | COTTON 230 | | calf, calf w/grip-top, thigh w/grip-top calf, thigh w/grip-top | 20–30 mmHg 30–40 mmHg | 10 |
| ESSENTIAL | OPAQUE 860 | | calf, petite calf, calf plus, calf w/grip-top, thigh w/grip-top, petite thigh w/grip-top, extra-large thigh w/grip-top, pantyhose, petite pantyhose, plus-sized pantyhose calf, calf plus, calf w/grip-top, thigh w/grip-top, extra-large thigh w/grip-top, thigh w/waist attachment | 20–30 mmHg 30–40 mmHg | 12, 14 |
| | SHEER 780 (for v | women only) | calf, thigh w/grip-top, pantyhose, maternity pantyhose calf, thigh w/grip-top, pantyhose | 15–20 mmHg 20–30 mmHg 30–40 mmHg | 18 |
| | MEDIUM SHEER 750 (for v | women only) | ● calf, thigh w/grip-top, pantyhose | 20–30 mmHg | 20 |
| STYLE | SOFT OPAQUE 840 (for v | women only) | calf, calf w/grip-top, thigh w/grip-top, pantyhose, maternity pantyhose calf, thigh w/grip-top, pantyhose | 15–20 mmHg 20–30 mmHg 30–40 mmHg | 22 |
| | PATTERNS 710 (for v | women only) | ● thigh w/grip-top, pantyhose | 20–30 mmHg | 24 |
| | SEA ISLAND COTTON 22 | 0 | ● calf | 20–30 mmHg | 26 |
| | MERINO WOOL 240 | | ● calf | 20–30 mmHg | 28 |
| | LINEN 250 | | • calf | 20–30 mmHg | 30 |
| | MICROFIBER 820 (fo | or men only) | calf, calf w/grip-top, thigh w/grip-top, extended foot sizes | 15–20 mmHg 20–30 mmHg 30–40 mmHg | 32 |
| | MICROFIBER PATTERNS | 830 | ● calf | 20–30 mmHg | 34 |
| | CUSHIONED COTTON 36 | 50 | ● calf | 20–30 mmHg | 38 |
| MOTION | HIGH TECH 412C | | ● calf | 20–30 mmHg | 40 |
| MOT | THERMOREGULATING W | VOOL 422 | ● calf | 20–30 mmHg | 42 |
| | COMFORT 850 | | ● calf | 20–30 mmHg | 44 |
| SPECIALTY | NATURAL RUBBER 500 | | ● calf, thigh, thigh w/grip-top, thigh w/waist attachment, pantyhose | 30–40 mmHg 40–50 mmHg 50–60 mmHg | 48 |
| SPEC | ADVANCE ARMSLEEVE | 910 | • armsleeve w/ grip-top w/gaunlet and without gauntlet | 15–20 mmHg 20–30 mmHg 30–40 mmHg | 50 |
| | | | | | |

PRODUCT **CASUAL COTTON** 146/186 TRAVENO TRAVEL SOCKS 14 CUSHIONED COTTON 142/182 SEA ISLAND COTTON 151/191 ALL-SEASON MERINO WOOL CONSUMER HEALTH LINEN 153/193 MICROFIBER SHADES 143/183 SHEER FASHION 120 SOFT SILHOUETTE LEGGINGS ATHLETIC RECOVERY SOCKS **PERFORMANCE SLEEVES** 412 MERINO OUTDOOR SOCKS **EVERSOFT DIABETIC SOCKS** MEDICALLY COMPLEX EDEM & INELASTIC WRAPS ACCESSORIES SYMPTOMS & INDICATIONS GLOSSARY

Closed Toe = • Open Toe = •

| PRODUCT | STYLES | COMPRESSION LEVEL | PAGE |
|--|---|-------------------|------|
| CASUAL COTTON 146/186 | • calf | 15–20 mmHg | 54 |
| TRAVENO TRAVEL SOCKS 146/T - 186/T | ● calf | 15–20 mmHg | 56 |
| CUSHIONED COTTON 142/182 | ● calf | 15–20 mmHg | 58 |
| SEA ISLAND COTTON 151/191 | • calf | 15–20 mmHg | 60 |
| ALL-SEASON MERINO WOOL 152/192 | ● calf | 15–20 mmHg | 62 |
| LINEN 153/193 | ● calf | 15–20 mmHg | 64 |
| MICROFIBER SHADES 143/183 | • calf | 15–20 mmHg | 66 |
| SHEER FASHION 120 | calf, thigh w/grip-top, pantyhose calf | 15–20 mmHg | 68 |
| SOFT SILHOUETTE LEGGINGS 170 | • open foot | 15–20 mmHg | 70 |
| ATHLETIC RECOVERY SOCKS 401 | ● calf | 15–20 mmHg | 72 |
| PERFORMANCE SLEEVES 412V | • open foot | 20–30 mmHg | 74 |
| MERINO OUTDOOR SOCKS 421 | ● calf | 15–20 mmHg | 76 |
| EVERSOFT DIABETIC SOCKS 160 | ● calf | 8–15 mmHg | 78 |
| MEDICALLY COMPLEX EDEMA & INELASTIC WRAPS | | | 80 |
| ACCESSORIES | | | 98 |
| SYMPTOMS & INDICATIONS | | | 103 |
| GLOSSARY | | | 104 |
| HOW TO MEASURE & ORDER | | | 106 |

Natural fibers



sigvaris

Be yourself.

ESSENTIAL

Support Yourself

Why choose **Essential**?



STYLE

Express Yourself

Why choose **Style**?

MOTION

Push Yourself

Why choose **Motion**?







Support Yourself

- Why choose **Essential**?
- Long-distance travel

| ESSENTIAL PR |
|--------------|
| Series |
| 230 |
| 860 |
| 860 |
| |

BEST SUITED FOR

ESSENTIAL

• Reduce swelling and discomfort

• Everyday comfort and reliability

RODUCTS

| Product | Page |
|----------------|------|
| Cotton | 10 |
| Opaque (Women) | 12 |
| Opaque (Men) | 14 |

• Containment of various types of edema, including venous edema,

phlebolymphedema, and early stages of lymphedema





Compression Level



Features and Benefits



- 100% Supima Cotton against the skin
- Ideal for sensitive skin and allergies
- Breathable and durable
- Excellent containment of various types of edema
- \cdot Classic ribbed design for everyday wear

Fiber Content

- 50% Nylon, 25% Spandex, 25% Supima Cotton (content can vary slightly by style)
- Latex free



20-30 mmHg Closed Toe Calf w/ Grip-top Medium Long in Black

232CMLM99/S



Silicone dotted grip-top available for calf and thigh



Scan here to see product images online siavaris.com

SIZING CHART

| SIZING CHART | | | | | | | | | | |
|---------------------------------|--------------------------|--------------------------------|------------------------------------|--------------------------------|----------------------------------|--------------------------------|--------------------------------|-------------------------|--|--|
| | Sm | nall | Med | lium | Lar | rge X L | | arge | | |
| Ankle Circumference | 18–21 (7"– | .5 cm 8.5") | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | | | |
| Calf Circumference | | 9 cm 15.5") | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | | | |
| | Short | Long | Short | Long | Short | Long | Short | Long | | |
| Calf Length | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | > 41 cm (16") | | |
| Thigh Circumference | | 45.5–66 cm (18"–26") | | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | | | | |
| | Short | Long | Short | Long | Short | Long | | | | |
| Leg Length: Women & Open Toe | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | NA | | | |
| | Short | Long | Short | Long | Short | Long | | | | |
| Leg Length: Men | ≤76 cm (30") | > 76 cm (30") | ≤76 cm (30") | > 76 cm (30") | ≤76 cm (30") | > 76 cm (30") | | | | |
| Size | SS | SL | MS | ML | LS | LL | XS | XL | | |
| | | | | | | | | | | |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

| Closed Toe = ● Open Toe = ○ | WHITE 00 | NAVY 10 | BLACK MIST 14 | LIGHT BEIGE 66 | BLACK 99 |
|--|-------------|------------|------------------|-------------------|-------------|
| Calf (C)* | • | • | • | •0 | • |
| Calf w/Grip-Top (C/S) Silicone dotted grip-top band | | | | •** | • |
| Thigh w/Grip-Top (N) Silicone dotted grip-top band | | | | •0 | • |

* 30–40 mmHg extra large calf only available in black with grip-top ** 20–30 mmHg only

Product Code Example: 232NMLW66

| Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour | Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour | Grip-top |
|--------|---------------|-------|------|---------------------|--------|--------|---------------|-------|------|---------------------|--------|----------|
| 23 | 2 | Ν | ML | W | 66 | 23 | 2 | с | LS | М | 66 | /s |

Product Code Example: 232CLSM66/S



OPAQUE 860



Features and Benefits



- Largest size range, from petite to plus
- Durable construction combined with a smooth opaque look
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema
- Plus size pantyhose are an alternative to maternity pantyhose

Fiber Content

- 64% Nylon, 36% Spandex (content can vary slightly by style)
- Latex free



20-30 mmHg Closed Toe Calf Medium Long in Black



Silicone dotted grip-top available for calf and thigh



Scan here to see product images online sigvaris com

| SIZING CHART | | | | | | | | | | | |
|---|---------------------------|-------------------------------------|---------------------------|--------------------------------|---------------------------------|-----------------------------|-------------------------------------|-----------------------------|-----------------------------|-----------------------------|--|
| | | Small | | | Medium | | Large | | X Large | | |
| Ankle Circumference | 18- | - 21.5 cm (7"–8 | .5") | 21.5–25.5 cm (8.5"–10") | | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | | |
| Calf Circumference | 28 | -39 cm (11"–15 | .5") | 31.5- | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cr | n (14"–20") | 41–61 cm | ı (16"–24") | |
| Calf Circumference: Calf Plus w/Grip-Top | | NA | | NA | 39.5–53 cm (15.5"–21") | | 45.5–61 cr | n (18"–24") | 51–71 cm (20"–28") | | |
| | Petite | Short | Long | Petite | Short | Long | Short | Long | Short | Long | |
| Calf Length | ≤ 35.5 cm (14") | ≤41 cm (16") | > 41 cm (16") | ≤35.5 cm (14") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | |
| Thigh Circumference | 45. | 5–66 cm (18"–: | 26") | 48- | 48-73 cm (19"-28.5") | | | 52–84 cm (20.5"–33") | | 64–95 cm (25"–37.5") | |
| | Petite | Short | Long | Petite | Short | Long | Short | Long | Short | Long | |
| Leg Length | ≤65 cm (25.5") | ≤72 cm (28.5") | > 72 cm (28.5") | ≤65 cm (25.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | >72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | |
| Hip Circumference | | 83–112.5 cm (32.5"–44.5") | | | 87–127 cm (34.5"–50") | | 98–150.5 cm (38.5"–59.5") | | NA | | |
| Plus or Maternity Size Hip Circumference | | ≤142 cm (56") | | : | ≤ 157.5 cm (62' |) | ≤178 cm (70") | | NA | | |
| Size | SP | SS | SL | MP | MS | ML | LS | LL | XS | XL | |
| Calf Plus w/Grip-Top | | | | | MS/P | ML/P | LS/P | LL/P | XS/P | XL/P | |
| | - | | | - | | | | | | | |

Women's Shoe Size : Petite Leg 4–8, Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

| STILLS AND COLOORS | | | | |
|---|-------------|--------------|----------------------|-------------|
| Closed Toe = ● Open Toe = ○ | HONEY 33 | GOLDEN 36 | LIGHT BEIGE 66 | BLACK 99 |
| Calf (C) | • | • | Ο | •0 |
| Calf Petite Sizes (C) Only available in 20–30 mmHg, not available with dotted grip-top band | • | • | | • |
| Calf w/Grip-Top (C/S) Silicone dotted grip-top band | | | •0 | • |
| Calf Plus w/Grip-Top (C/P) Only available with silicone dotted grip-top band | | | •0 | |
| Thigh w/Grip-Top (N)** Silicone dotted grip-top band | • | • | 0 | • |
| Thigh Petite Sizes (N) Only available in 20–30 mmHg | • | • | | • |
| Pantyhose (P) | • | • | | • |
| Pantyhose Petite (P) Only available in 20–30 mmHg | • | • | | • |
| Pantyhose Plus and Maternity (P/P) | | | • | • |
| Thigh w/Waist Attachment (W)* Indicate leg LEFT (/L) or RIGHT (/R) | | | 0 | |
| * Not available in XL with waist attachment | | | | |

** Thigh w/Grip-Top (N) and Open Toe Thigh w/Grip-Top (N) sizes XS and XL are only available in Light Beige 66

Product Code Example: 863PMSW33

| Series | mmHg Level | Style | Size | Gender, Open To |
|--------|---------------|-------|------|--------------------|
| 86 | 3 | P | MS | W |

862CMLW99

OPAQUE 860

Product Code Example: 862CLLO66/P

| olour | Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour | Grip-top with Plus size |
|-------|--------|---------------|-------|------|---------------------|--------|----------------------------|
| 33 | 86 | 2 | С | LL | 0 | 66 | /P |



OPAQUE 860





Features and Benefits



- Durable construction with a ribbed design on closed toe calf products
- Smooth opaque design on open-toe styles.
- Comfortable, stretchable material easily glides on and off
- Containment of various types of edema

Fiber Content

- 64% Nylon, 36% Spandex
- (content can vary slightly by style)
- Latex free



20–30 mmHg Open Toe Thigh Medium Long in Light Beige



Silicone dotted grip-top available for calf and thigh



Scan here to see product images online sigvaris com

SIZING CHART

| | _ | | | | | | | | |
|---|--------------------------|--------------------------------|-----------------------------------|------------------------------------|-----------------------------------|--------------------------------|--------------------------------|------------------------------|--|
| | Sm | all | Med | lium | La | rge | X La | arge | |
| Ankle Circumference | 18–21 (7"– | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | | |
| Calf Circumference | | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | il cm -24") | |
| Calf Circumference: Calf Plus w/Grip-Top | N | NA | | 39.5–53 cm (15.5"–21") | | 45.5–61 cm (18"–24") | | 51–71 cm (20"–28") | |
| | Short | Long | Short | Long | Short | Long | Short | Long | |
| Calf Length | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | |
| Thigh Circumference | 45.5–66 cr | n (18"–26") | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | | 64–95 cm (25"–37.5") | | |
| | Short | Long | Short | Long | Short | Long | Short | Long | |
| Leg Length | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | |
| Size | SS | SL | MS | ML | LS | LL | xs | XL | |
| Calf Plus w/Grip-Top | | | MS/P | ML/P | LS/P | LL/P | XS/P | XL/P | |

Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = • Open Toe = **o**

Calf (C)

Calf w/Grip-Top (C/S) Silicone dotted grip-top band

Calf Plus w/Grip-Top (C/P)

Only available with silicone dotted grip-top band Open Toe Thigh w/Grip-Top (N)

Silicone dotted grip-top band Thigh w/Waist Attachment (W)** Indicate leg LEFT (/L) or RIGHT (/R)

* 20-30mmHg only

** Not available in XL with waist attachment

Product Code Example: 862CMLM99

| Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour | Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour | Grip-top with Plus size |
|--------|---------------|-------|------|---------------------|--------|--------|---------------|-------|------|---------------------|--------|----------------------------|
| 86 | 2 | С | ML | М | 99 | 86 | 2 | С | LL | 0 | 66 | /P |

862NMLO66

OPAQUE 860

RIBBED BLACK 99 SMOOTH BLACK 99 RIBBED LIGHT BEIGE 66 SMOOTH LICHT BEICE 66 •* 0 0 •* 0 0 0 0

Product Code Example: 862CLLO66/P





Express Yourself

- Why choose Style?
- \cdot Special occasion and formal wear
- Our most luxurious feel

STYLE PRODUCTS

| Series | Product | Page |
|--------|---------------------|------|
| 780 | Sheer | 18 |
| 750 | Medium Sheer | 20 |
| 840 | Soft Opaque | 22 |
| 710 | Patterns | 24 |
| 220 | Sea Island Cotton | 26 |
| 240 | Merino Wool | 28 |
| 250 | Linen | 30 |
| 820 | Microfiber | 32 |
| 830 | Microfiber Patterns | 34 |

BEST SUITED FOR


- Management of Chronic Venous Disease (CVD)
- \cdot Pre- and post-vein treatment with minimal edema
- Pregnancy
- People who prefer fashionable options
- Ideal for business or dress



• Fashionable, modern patterns and designs



SHEER 780



Elegant **Breathable**

- Fashionable for professional and dress wear
- Consistent sheerness throughout the leg
- Excellent breathability and moisture management
- Maternity pantyhose options available

Fiber Content

Sheer

- 15–20 and 20–30 mmHg: 69% Nylon, 31% Spandex (content can vary slightly by style)
- 30–40 mmHg: 60% Nylon, 40% Spandex (content can vary slightly by style)
- Latex free



15–20 mmHg Closed Toe Thigh Medium Long in Dark Navy



Sensinnov grip-top on thigh



Scan here to see product images online siavaris.com

SIZING CHART

| | Sm | nall | Med | lium | Lar | ge |
|-----------------------------|--------------------------|---------------------------|---|---------------------------|----------------------------------|---------------------------|
| Ankle Circumference | 18–21 (7"– | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | |
| Calf Circumference | 28–3 (11"–1 | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| | Short | Long | Short | Long | Short | Long |
| Calf Length | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") |
| Thigh Circumference | 45.5-6 (18"- | | 48-73 cm (19"–28.5") 52-84 cm (20.5"–33") | | | |
| | Short | Long | Short | Long | Short | Long |
| Leg Length | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") |
| Hip Circumference | 83–112 (32.5"- | | | 2 7 cm "–50") | 98–150 (38.5"- | |
| Maternity hip circumference | ≤142 (50 | | ≤ 157.5 cm (62") ≤178 c (70" | | | |
| Size | SS | SL | MS | ML | LS | LL |
| | | | | | | |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

| STILES AND COLOU | | | | | | | | | |
|---|--------------------|-------------|---------------------|--------------|--------------------------|-------------|--------------|-----------------------|-------------|
| Closed Toe = ● Open Toe = ○ | DARK NAVY 08 | HONEY 33 | WARM SAND 34* | GOLDEN 36 | TOASTED ALMOND 37* | CAFÉ 73* | MOCHA 85* | NIGHT- SHADE 94 | BLACK 99 |
| Calf (C) | • | •0 | •0 | •0 | •0 | •0 | •0 | • | •0 |
| Thigh w/Grip-Top (N) Sensinnov grip-top band | • | •0 | •0 | •0 | •0 | •0 | •0 | • | •0 |
| Pantyhose (P) | • | •0 | •0 | •0 | •0 | •0 | •0 | • | •0 |
| Maternity Pantyhose (M)* | | • | | | | | | | • |

* 20–30 mmHg only

Product Code Example: 781NMLW36

| Series | mmHg Level | Style | Size | Gender/ Open To |
|--------|---------------|-------|------|--------------------|
| 78 | 1 | N | ML | W |

781NMLW08

36

Product Code Example: 783PSSO33

| Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour |
|--------|---------------|-------|------|---------------------|--------|
| 78 | 3 | Р | SS | 0 | 33 |





Compression Level

20–30 mmHg

Features and Benefits



- \cdot Fashionable, durable, medium sheer hosiery for everyday wear
- \cdot Excellent breathability and moisture management
- Comfortable 2" top band on calf style

Fiber Content

- 68% Nylon, 32% Spandex (content can vary slightly by style)
- Latex free



20-30 mmHg Closed Toe Thigh Medium Long in Black

Silicone dotted grip-top on thigh



Scan here to see product images online sigvaris.com

SIZING CHART

| SIEING CHAN | | | | | | |
|---------------------|--------------------------|---------------------------|------------------------------------|---------------------------|--|---------------------------|
| | Sm | nall | Med | lium | Lar | ge |
| Ankle Circumference | | .5 cm 8.5") | | 5.5 cm '–10") | 25.5–2 (10"– | |
| Calf Circumference | | '9 cm 15.5") | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| | Short | Long | Short | Long | Short | Long |
| Calf Length | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") |
| Thigh Circumference | | 66 cm -26") | | '3 cm 28.5") | (16") (16") 52-84 cm (20.5"-33") | |
| | Short | Long | Short | Long | Short | Long |
| Leg Length | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") |
| Hip Circumference | | 2.5 cm -44.5") | | 2 7 cm "–50") | 98–150 (38.5"- | |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

| STYLES AND COLOURS | | | |
|---|-------------|-------------------|-------------|
| Closed Toe = ● | HONEY 33 | NIGHT-SHADE 94 | BLACK 99 |
| Calf (C) | ٠ | ٠ | • |
| Thigh w/Grip-Top (N) Silicone dotted grip-top band | • | • | • |
| Pantyhose (P) | ٠ | • | • |

Product Code Example: 752NMLW33

| Series | mmHg Level | Style | Size |
|--------|---------------|-------|------|
| 75 | 2 | N | ML |

752NMLW99

| Gender | Colour |
|--------|--------|
| W | 33 |



OPAQUE 840



Compression Level



Features and Benefits



- Fashionable opaque hosiery for professional and everyday wear
- \cdot Long-lasting microfiber softness provides all day comfort and allows product to easily glide on and off
- Durable and breathable
- Maternity pantyhose option available

Fiber Content

- 79% Nylon, 21% Spandex
- (content can vary slightly by style)
- Latex free

Scan here to see product images online sigvaris.com

SIZING CHART

| | Sm | nall | Mec | lium | Lai | ge | |
|-----------------------------|--------------------------|---------------------------|-----------------------------------|---------------------------|----------------------------------|---------------------------|--|
| Ankle Circumference | 18–21 (7"– | .5 cm 8.5") | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | | |
| Calf Circumference | 28–3 (۱۱"–1 | | | 45 cm –17.5") | 35.5–51 cm (14"–20") | | |
| | Short | Long | Short | Long | Short | Long | |
| Calf Length | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | >41 cm (16") | |
| Thigh Circumference | | 56 cm -26") | | | | 34 cm 5"–33") | |
| | Short | Long | Short | Long | Short | Long | |
| Leg Length | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | |
| Hip Circumference | 83–112 (32.5"- | 2.5 cm -44.5") | | 2 7 cm "–50") | 98–15 ((38.5"- | | |
| Maternity hip circumference | ≤142 (50 | ? cm 6") | | | ≤17 8 (70 | | |
| Size | SS | SL | MS ML | | LS | LL | |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS

| STILES AND COLOURS | | | | | | | |
|---|------------------------|------------|-------------|-----------------|-----------------|----------------|-------------|
| Closed Toe = ● Open Toe = ● | MIDNIGHT BLUE 09 | CHAI 35 | PECAN 40 | MULBERRY 60* | ESPRESSO 89* | GRAPHITE 91 | BLACK 99 |
| Calf (C) | • | •0 | • | • | • | • | •0 |
| Calf w/Grip-Top (C/S)* Silicone dotted grip-top | | | | | | | • |
| Thigh w/Grip-Top (N) Sensinnov [®] grip-top band | • | •0 | • | • | • | • | •0 |
| Pantyhose (P) | • | •0 | • | • | • | • | •0 |
| Maternity Pantyhose (M)* | | | | | | | • |

* 20–30 mmHg only

Product Code Example: 842CMSW09

| Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour |
|--------|---------------|-------|------|---------------------|--------|
| 84 | 2 | С | MS | W | 09 |



20-30 mmHg Closed Toe Pantyhose Small Long in Black





Sensinnov grip-top on thigh and silicone dotted grip-top available on calf



SOFT OPAQUE 840



PATTERNS 710

Compression Level

20–30 mmHg

Features and Benefits



- Stylish diamond-dot and stripe patterns for professional and dress wear
- \cdot Sheer hosiery with excellent breathability and moisture management
- Unique and elegant for those looking for a fashionable option

Fiber Content

- 64% Nylon, 36% Spandex (content can vary slightly by style)
- Latex free



20-30 mmHg Closed Toe Thigh Medium Long in Graphite





Sensinnov Scallopped for Onyx Stripes and Sensinnov band on thigh



Scan here to see product images online sigvaris.com

SIZING CHADT

| SIZING CHART | | | | | | |
|---------------------|-------------------------------------|--------------------------------|------------------------------------|-----------------------------------|-------------------------------------|---------------------------|
| | Small | | Medium | | Large | |
| Ankle Circumference | | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"−JO") | | !9 cm 11.5") |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| Thigh Circumference | 45.5–66 cm (18"–26") | | 48–73 cm (19"–28.5") | | 52–84 cm (20.5"–33") | |
| | Short | Long | Short | Long | Short | Long |
| Leg Length | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") | ≤72 cm (28.5") | > 72 cm (28.5") |
| Hip Circumference | 83-112.5 cm (32.5°-44.5°) | | 87–127 cm (34.5"–50") | | 98–150.5 cm (38.5"–59.5") | |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12

STYLES AND COLOURS Closed Toe = • Thigh w/Grip-Top (N) Sensinnov grip-top band Pantyhose (P)

Product Code Example: 712NLSW19

| Series | mmHg Level | Style | Size | |
|--------|---------------|-------|------|--|
| 71 | 2 | Ν | LS | |

712NMLW91

| | XX | |
|----------------------|----------------|-------------|
| ONYX STRIPE 19 | GRAPHITE 91 | BLACK 99 |
| | | |
| • | | |
| • | • | • |

| Gender | Colour | | | |
|--------|--------|--|--|--|
| W | 19 | | | |



SEA ISLAND COTTON 220



Compression Level

20–30 mmHg

Also available in 15–20mmHg, see page 60

Features and Benefits



- Knitted with the world's rarest cotton for a soft,
- Ideal for sensitive skin
- Breathable and comfortable
- Suitable for those who prefer natural fibers for everyday wear
- Smooth, flat toe seam

Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex, 7% Supima Cotton
- Latex free



20-30 mmHg Closed Toe Calf Small Long in Brown





20–30 mmHg Closed Toe Calf Medium Long in Black

222CMLM99



Scan here to see product images online siavaris.com

SIZING CHART

| SIZING CHARI | | | | | | | | |
|---------------------|--------------------------------|---------------------------|------------------------------------|-------------------------|----------------------------------|-------------------------|-------------------------------|---------------------------|
| | Small | | Medium | | Large | | X Large | |
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | | | 25.5–29 cm (10"–11.5") | | 6 cm '-14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| | Short | Long | Short | Long | Short | Long | Short | Long |
| Calf Length | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Size | SS | SL | MS | ML | LS | LL | XS | XL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14


| STYLES AND COLOURS | |
|--------------------|--|
| Closed Toe = ● | |
| Calf (C) | |

Product Code Example: 222CLLM10

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|---------------|-------|------|--------|--------|
| 22 | 2 | С | LL | М | 10 |



SEA ISLAND COTTON 220







Compression Level



Also available in 15–20 mmHg, see page 62

Features and Benefits



- High quality extra-fine Merino wool makes this sock soft and breathable
- Moisture wicking capabilities keep feet dry
- Naturally odour resistant
- Thermoregulating for year-round wear
- Smooth, flat toe seam

Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free



20-30 mmHg Closed Toe Calf Medium Long in Brown





20–30 mmHg Closed Toe Calf Medium Long in Navy

242CMLW10



Scan here to see product images online siavaris.com

SIZING CHADT

| SIZINU CHARI | | | | | | | | | |
|---------------------|--------------------------------|---------------------------|------------------------|-------------------------|------------------------|----------------------------------|------------------------|-------------------------------|--|
| | Small | | Medium | | Large | | X Large | | |
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | | | | 25.5–29 cm (10"–11.5") | | 36 cm "–14") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | | | | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| | Short | Long | Short | Long | Short | Long | Short | Long | |
| Calf Length | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | > 41 cm (16") | |
| Size | SS | SL | MS | ML | LS | LL | XS | XL | |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7-11, Long Leg 9.5-14

| STYLES AND COLOURS | |
|--------------------|--|
| Closed Toe = ● | |
| Calf (C) | |

Product Code Example: 242CMSW99

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|---------------|-------|------|--------|--------|
| 24 | 2 | С | MS | W | 99 |

MERINO WOOL 240





LINEN 250



Compression Level



Also available in 15–20 mmHg, see page 64

Features and Benefits



Hypoallergenic Natural fibers

Moisture Management

Strong and Durable

• Lightweight and breathable for all day comfort

- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free



20-30 mmHg Closed Toe Calf Medium Long in Brown





20–30 mmHg Closed Toe Calf Medium Long in Denim

252CMLW201



Scan here to see product images online siavaris.com

SIZING CHART


| | Small | | Medium | | Large | |
|---------------------|--------------------------------|---------------------------|------------------------------------|---------------------------|----------------------------------|---------------------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | |
| | Short | Long | Short | Long | Short | Long |
| Calf Length | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

| STYLES AND COLOURS | |
|--------------------|---|
| Closed Toe = ● | B |
| Calf (C) | |

Product Code Example: 252CMLW201

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|---------------|-------|------|--------|--------|
| 25 | 2 | С | ML | W | 201 |





MICROFIBER 820



Compression Level



Features and Benefits



- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking.
- Glides on and off with ease
- Comfortable and durable sock for professional and everyday wear
- Available in extended foot sizes
- Smooth, flat toe seam

Fiber Content

- 87% Nylon , 13% Spandex
- (content can vary slightly by style)
- Latex free



30-40 mmHg Closed Toe Thigh Medium Long in Black





Silicone dotted grip-top option available for calf and designed with Sensinnov grip-top for thigh



Scan here to see product images online siavaris.com

SIZING CHART

| SIZINU CHARI | | | | | | | | |
|---------------------|------------------------|--------------------------------|------------------------|-------------------------|----------------------------------|-------------------------|--------------------------------|---------------------------|
| | Small | | Medium | | Large | | X Large | |
| Ankle Circumference | 18–21 (7"– | .5 cm 8.5") | | 5.5 cm '–10") | 25.5–29 cm (10"–11.5") | | 29–36 cm (11.5"–14") | |
| Calf Circumference | | 28–39 cm (11"–15.5") | | 45 cm –17.5") | 35.5–51 cm (14"–20") | | 41–61 cm (16"–24") | |
| | Short | Long | Short | Long | Short | Long | Short | Long |
| Calf Length | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Thigh Circumference | 45.5-6 (18"- | | | '3 cm 28.5") | | 4 cm "–33") | | |
| | Short | Long | Short | Long | Short | Long | NA | |
| Leg Length | ≤76 cm (30") | > 76 cm (30") | ≤76 cm (30") | > 76 cm (30") | ≤76 cm (30") | > 76 cm (30") | | |
| Size | SS 7-11 | SL 9.5–14 | MS 7-11 | ML 9.5–14 | LS 7–11 | LL 9.5–14 | XS 7–11 | XL 9.5–14 |
| Shoe Size | | | | MX* 12.5–17.5 | | LX* 12.5–17.5 | | XX* 12.5–17.5 |

Men's shoe sizes: Short Leg 7-11, Long Leg 9.5-14 (For extended foot sizes: long leg 12.5-17.5) *MX, LX, and XX extended foot sizes only available in calf style in Black & Tan Khaki in 15–20 & 20–30 mmHg only. Extended foot sizes not available with grip top

| STYLES AND COLOURS | |
|--|----|
| Closed Toe = ● | ТА |
| Calf (C) | |
| Calf w/Grip-Top (C/S)** Silicone dotted grip-top band | |
| Thigh w/Grip-Top (N) Sensinnov grip-top band | |
| | |

** Calf w/ Grip-Top not available in 15–20 mmHg *** Steel Grey only available in 15–20 mmHg and 20–30 mmHg

Product Code Example: 822NLSM99

| Series | mmHg Level | Style | Size | Gender | Colour | Series | mmHg Level | Style | Size | Gender | Colour | Grip-top |
|--------|---------------|-------|------|--------|--------|--------|---------------|-------|------|--------|--------|----------|
| 82 | 2 | N | LS | М | 99 | 82 | 2 | С | ML | М | 32 | /s |

823NMLM99



| Product Code | Example: | 822CMLM32/ | S |
|--------------|----------|------------|---|
|--------------|----------|------------|---|



MICROFIBER PATTERNS 830




Compression Level



Also available in 15–20 mmHg, see page 66

Features and Benefits



- 100% of microfiber on inside of sock to improve comfort, softness, and facilitate moisture wicking.
- · Fashionable, durable sock for everyday wear
- Glides on and off with ease
- Smooth, flat toe seam

Fiber Content

- 87% Nylon, 13% Spandex
- (content can vary slightly by style)
- Latex free



20–30 mmHg Closed Toe Calf Small Long in Royal Blue Argyle

832CSLM48



20–30 mmHg Closed Toe Calf Large Short in Pink Stripe

832CLSW15



Scan here to see product images online siavaris.com

SIZING CHART

| SIZING CHARI | | | | | | | |
|---------------------|--------------------------------|---------------------------|------------------------------------|-------------------------|---------------------------------|-------------------------|--|
| | Small | | Мес | dium | Large | | |
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | | 29 cm -11.5") | |
| Calf Circumference | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | 35.5–51 cm (14"–20") | | |
| | Short | Long | Short | Long | Short | Long | |
| Calf Length | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | |
| Size | SS | SL | MS | ML | LS | LL | |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

STYLES AND COLOURS

Closed Toe = •

Calf (C)

Calf (C)

* Available for women only ** Available for men only


Product Code Example: 832CMLM17

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|---------------|-------|------|--------|--------|
| 83 | 2 | С | ML | М | 17 |

MICROFIBER PATTERNS 830

| MARINER | GRAPHITE | DARK NAVY | PINK | GRAPHITE | ONYX | GRAPHITE |
|---------|-------------|-----------|--------|----------|------------|----------|
| STRIPE | MINI-STRIPE | STRIPE | STRIPE | STRIPE | STRIPE | HEATHER |
| 02* | 07 | 13 | 15* | 17 | 19 | 27 |
| • | • | • | • | • | • | • |
| GREY | BLACK | PINK | PURPLE | CRAPHITE | ROYAL BLUE | ONYX |
| PLAID | PLAID | ARGYLE | ARGYLE | ARCYLE | ARGYLE | ARGYLE |
| 22* | 42** | 45* | 46 | 47 | 48 | 49 |
| • | • | • | • | • | • | • |





Push Yourself

- Why choose **Motion**?

MOTION PRODUCTS

| ç | Series |
|---|--------|
| | 360 |
| | 412 |
| | 422 |
| | 850 |
| | |

BEST SUITED FOR

| Product | Page |
|-----------------------|------|
| Cushioned Cotton | 38 |
| High Tech | 40 |
| Thermoregulating Wool | 42 |
| Comfort | 44 |



CUSHIONED COTTON 360



Compression Level



Also available in 15–20 mmHg, see page 58

Features and Benefits



- $\cdot\,$ Supima cotton for highest level of comfort and breathability
- Extra-padded sole, heel, and toe to protect feet
- Ideal for sensitive skin and allergies
- Flat, low friction toe seam
- Diabetic friendly
- Ideal for people who stand for prolonged periods of time



Closed Toe Calf Medium Long in



Fiber Content

- 64% Supima Cotton, 28% Nylon, 8% Spandex
- Latex free



Scan here to see product images online sigvaris com

SIZING CHART

| | Sm | nall | Med | lium | La | rge | X La | arge | | |
|--|------------------------|--------------------------------|------------------------|---------------------------|------------------------|---------------------------|------------------------|-------------------------|--|--|
| Ankle Circumference | | 18–21.5 cm (7"–8.5") | | 5.5 cm '–10'') | | 29 cm -11.5") | | 6 cm –14") | | |
| Calf Circumference | 28–3 (11"–1 | 9 cm 15.5") | | 45 cm –17.5") | | 51 cm -20") | 41–6 (16"- | 1 cm 24") | | |
| | Short | Long | Short | Long | Short | Long | Short | Long | | |
| Calf Length | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | >41 cm (16") | ≤41 cm (16") | > 41 cm (16") | | |
| Size | SS | SL | MS | ML | LS | LL | XS | XL | | |
| Nomon's Shoo Sizos: Short Log 55 95 Long | 100.8.12 | | | | | | | | | |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe Sizes: Short Leg 7–11, Long Leg 9.5–14



Product Code Example: 362CMLM00

| Series | mmHg Level | Style | Size |
|--------|---------------|-------|------|
| 36 | 2 | С | ML |



20–30 mmHg

Black

CUSHIONED COTTON 360

| Gender | Colour |
|--------|--------|
| М | 00 |



HIGH TECH 412



Compression Level

20–30 mmHg

Features and Benefits



Moisture Wicking

Anatomically Contoured

- High-Tech fibers for moisture wicking and improved thermoregulation
- Cushioned footbed, toe, and Achilles tendon zones
 for optimal comfort
- Anatomically contoured to support the natural shape of the foot and calf
- $\cdot\,$ Smooth, flat toe seam
- $\cdot\,$ Perfect for everyday activities and active individuals



20–30 mmHg Closed Toe Calf Small Medium in Limeade



20–30 mmHg Closed Toe Calf Medium Large in Pink

412CML56

Fiber Content

- 64% Nylon, 17% Polyester, 19% Spandex
- Latex free



Scan here to see product images online sigvaris.com

SIZING CHART

| SIZING CHARI | | | | | | | | | | | | |
|---------------------|--------------------------------|--|---------|---------|--------------------------------|----------------------------------|---------|--------------------------------|-------|---------|---------|-----|
| | | Small | | I | Mediun | ۱ | | Large | | | X Large | • |
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | : | | | 25.5–29 cm (10"–11.5") | | 29–33 cm (11.5"–13") | | | | |
| Calf Circumference | | 28–39 cm 31.5–45 cm 35.5–51 cm (11"–15.5") (12.5"–17.5") (14"–20") | | | 40.5–58 cm (16"–23") | | | | | | | |
| Women's Shoe Size | 4.5–7.5 | 8–10.5 | 11–13.5 | 4.5–7.5 | 8–10.5 | 11–13.5 | 8–10.5 | 11–13.5 | >14 | 8–10.5 | 11–13.5 | >14 |
| Men's Shoe Size | ≤6.5 7–9.5 10–12.5 | | ≤6.5 | 7–9.5 | 10–12.5 | 7–9.5 | 10–12.5 | 13–15.5 | 7–9.5 | 10–12.5 | 13–15.5 | |
| Size | SS | SM | SL | MS | ММ | ML | LS | LM | LL | XS | ХМ | XL |

| STYLES AND CO | STYLES AND COLOURS | | | | | | | | | | | |
|----------------|--------------------|------------|------------|------------|-----------|---------------------|-----------------------|---------------|-------------|--|--|--|
| Closed Toe = ● | WHITE 00 | BLUE 50 | LIME 54 | РІМК 56 | RED 59 | STEEL BLUE 62 | BLOOD ORANGE 63 | LIMEADE 64 | BLACK 99 | | | |
| Calf (C) | • | • | • | • | • | • | • | • | • | | | |

Product Code Example: 412CML62

| Series | mmHg Level | Style | Size | Colour |
|--------|---------------|-------|------|--------|
| 41 | 2 | С | ML | 62 |

2



THERMO-REGULATING WOOL 422



Compression Level



Also available in 15–20 mmHg, see page 76

Features and Benefits



- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for yearround wear
- Cushioned, padded soles help prevent blisters
- Soft and comfortable
- Smooth, flat toe seam
- Perfect for outdoor work environments and activities

Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- \cdot Latex free



Scan here to see product images online sigvaris.com

SIZING CHART

| SIZING CHART | | | | | | | | | | |
|---------------------|--------------------------------|--------------------------------|-----------------------------------|------------------------------------|--------|----------------------------------|---------|---------------------------------|-----|--|
| | | Small | | | Medium | | | Large | | |
| Ankle Circumference | 18–21.5 cm (7"–8.5") | | 21.5–25.5 cm (8.5"–10") | | | 25.5–29 cm (10"–11.5") | | | | |
| Calf Circumference | | 28–39 cm (11"–15.5") | | 31.5–45 cm (12.5"–17.5") | | | | 35.5–51 cm (14"–20") | | |
| Women's Shoe Size | 4.5–7.5 | 8–10.5 | 11–13.5 | 4.5–7.5 | 8–10.5 | 11–13.5 | 8–10.5 | 11–13.5 | >14 | |
| Men's Shoe Size | ≤6.5 | 7–9.5 | 10–12.5 | ≤6.5 7-9.5 10-12.5 | | 7–9.5 | 10–12.5 | 13–15.5 | | |
| Size | SS | SM | SL | MS | ММ | ML | LS | LM | LL | |

STYLES AND COLOURS

| Calf (C) | |
|----------------|--|
| Closed Toe = ● | |

Product Code Example: 422CML12

| Series | mmHg Level | Style | Size | | |
|--------|---------------|-------|------|--|--|
| 42 | 2 | С | ML | | |



20-30 mmHg Closed Toe Calf Large Long in Charcoal




20–30 mmHg Closed Toe Calf Medium Short in Olive

422CMS31

42

THERMOREGULATING WOOL 422







COMFORT 850



Compression Level

20–30 mmHg

Features and Benefits



- NILIT Aquarius fibers provide moisture wicking to keep feet dry and comfortable
- $\cdot\,$ Cushioned soles to absorb impact and help prevent blisters
- Smooth, flat toe seam
- Suitable for uniforms requiring solid black color

Fiber Content

- 86% Nylon, 14% Spandex
- Latex free



20–30 mmHg Closed Toe Calf Large Long in Black





852CMLW99

Scan here to see product images online siavaris.com

SIZING CHART

| | Sm | nall | Med | lium | La | rge |
|---------------------|--|-------------------------|-----------------------------------|-------------------------|----------------------------------|---------------------------|
| Ankle Circumference | 18–21.5 cm (7"–8.5") 28–39 cm (11"–15.5") | | 21.5–25.5 cm (8.5"–10") | | 25.5–29 cm (10"–11.5") | |
| Calf Circumference | | | | 45 cm –17.5") | 35.5–51 cm (14"–20") | |
| | Short | Long | Short | Long | Short | Long |
| Calf Length | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | > 41 cm (16") | ≤41 cm (16") | >41 cm (16") |
| Size | SS | SL | MS | ML | LS | LL |

Women's Shoe Sizes: Short Leg 5.5–9.5, Long Leg 8–12 Men's Shoe sizes: Short Leg 7–11, Long Leg 9.5–14

| STYLES AND COLOURS | |
|--------------------|--|
| | |

| | | |
|----------|------|--|
| Calf (C) | | |


Closed Toe = •

Product Code Example: 852CMLM99

| Series | mmHg Level | Style | Size | Gender | Colour |
|--------|---------------|-------|------|--------|--------|
| 85 | 2 | С | ML | М | 99 |







Take care of Yourself

- Why choose **Specialty**?
- Specific medical needs
- Need low fatigue fibers
- Natural rubber (500)

SPECIALTY PRODUCTS

| Series |
|--------|
| 500 |
| 910 |

BEST SUITED FOR

- Chronic edema post DVT
- Severe venous stasis
- \cdot Early stages of lymphedema
- Active venous ulcers and prevention of ulcer recurrence

SPECIALTY

| Product | Page |
|-------------------|------|
| Natural Rubber | 48 |
| Advance Armsleeve | 50 |

- People with early stages of arm lymphedema



NATURAL **RUBBER** 500 for women & men

Compression Level



Features and Benefits



Natural Rubber

Largest Size Range

- Durable, long-term containment solution with minimal fiber fatigue throughout the day
- Extended sizing for a custom-like fit
- Ideal for patients with more serious indications, including:
 - Severe varicosities
 - Chronic edema post DVT
 - Prevention and treatment of Active Venous Ulcers
 - Severe venous stasis
 - Lymphedema



30–40 mmHg Open Toe Thigh Small Average Long in Beige

Silicone dotted grip-top available on thigh

Fiber Content

- 70% Nylon, 30% Natural Latex Rubber (Calf & Thigh-high styles)
- 75% Nylon, 20% Natural Latex Rubber, 5% other fibers (Pantyhose style)
- **Caution:** contains natural rubber latex that may cause allergic reaction

SIZING CHART

| SIZING CHART | _ | | | | | | | | _ | | | | | | | |
|---|---------------------------|----------------------------|---------------------------|----------------------------|--------------------------------|----------------------------|---------------------------|----------------------------|---|--------------------------------|---------------------------|-------------------------|--------------------------------|----------------------------|---------------------------|----------------------------|
| | | Sm | nall | | | Med | lium | | | La | rge | | | X La | arge | |
| Ankle Circumference | | 19–23 cm | ר (7.5"–9") | | 23–30 cm (9"–12") | | | | 30–40 cm (12"–15.5") | | | | 40–46 cm (15.5"–18") | | | |
| | Ave | rage | Full | | Ave | Average | | Full | | Average Full | | III | Average | | Fu | ull |
| Calf Circumference | 28–3 (11"– | 8 cm -15") | ≤43 cm (17") | | 32–43 cm (12.5"–17") | | | ≤48 cm (19") | | 36–48 cm (14.5"–19") | | cm 5") | 45–57 cm (18"–22.5") | | ≤68 cm (27") | |
| Calf Circumference for 505C (Full only) | | : full only | 32–4 (12.5' | | | C: full only | 36–4 (14.5 | 8 cm '–19") | 505C: full 45–57 cm (18"–22.5") | | 505C: not available in XL | | XL | | | |
| | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long |
| Calf Length | ≤40 cm (≤15.5") | > 40 cm (>15.5") | ≤40 cm (≤15.5") | > 40 cm (>15.5") | ≤40 cm (≤15.5") | > 40 cm (>15.5") | ≤40 cm (≤15.5") | > 40 cm (>15.5") | ≤40 cm (≤15.5") | > 40 cm (>15.5") | ≤40 cm (≤15,5") | >40 cm (>15.5") | ≤40 cm (≤15.5") | > 40 cm (>15.5") | ≤40 cm (≤15.5") | > 40 cm (>15.5") |
| | Ave | rage | F | III | Ave | rage | F | ıll | Ave | rage | Fu | ıll | | | | |
| Thigh Circumference | 41–5 4 (16"- | | <u><</u> 61 (2- | | | 51 cm -24") | ≤71 (2) | cm B") | | '1 cm -28") | <u><</u> 81 (32 | | | | | |
| Thigh Circumference for 505T (Full only) | 5051 thigh | | <u>≤</u> 61 (24 | | | F: full n only | ≤71 (2 | cm B") | | f: full n only | <u><</u> 81 (32 | | | N | /^ | |
| | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | Short | Long | | IN, | A | |
| Leg Length | ≤71 cm (28") | > 71 cm (28") | ≤71 cm (28") | > 71 cm (28") | ≤71 cm (28") | > 71 cm (28") | ≤71 cm (28") | > 71 cm (28") | ≤71 cm (28") | > 71 cm (28") | ≤71 cm (28") | > 71 cm (28") | | | | |
| Hip Circumference | | | 2 cm 8") | | | | 7 cm 4") | | | | 5 cm 51") | | | | | |
| Size | S1 | S2 | S 3 | S4 | M1 | M2 | M3 | M4 | LI | L2 | L3 | L4 | XI | X2 | Х3 | X4 |

STYLES AND COLOURS

| Open Toe = o | BEIGE 77 | Caution: this product contains natural rubber latex which may cause allergic reactions. For all indications, cautions, and contraindications, please see page 101 for details. Sigvaris accepts no liability for non-observance of contraindications or cautions. |
|---|-------------|--|
| Calf (C) | o | Calf only available in S, M, and L full circumference for 50–60 mmHg |
| Thigh w/out Grip-Top (T) | o | |
| Thigh w/Grip-Top (N) Silicone dotted grip-top band | ο | Not available in 50-60 mmHg or XL |
| Pantyhose (P) | o | Pantyhose only available in 30-40 mmHg and only in full thigh circumference |
| Thigh w/Waist Attachment (W) Indicate Leg LEFT (/L) or RIGHT (/R) | ο | Thigh w/ waist attachment is available in average and full thigh circumferences in 30–40 mmHg and in full thigh circumferences only in 40–50 mmHg |

Product Code Example: 504CS2O77

| Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour | Series | mmHg Level | Style | Size | Gender/ Open Toe | Colour | Leg |
|--------|---------------|-------|------|---------------------|--------|--------|---------------|-------|------|---------------------|--------|-----|
| 50 | 4 | С | S2 | 0 | 77 | 50 | 3 | W | LI | 0 | 77 | /L |

503NS2077

Product Code Example: 503WL1077/L





20–30 mmHg Armsleeve With Gauntlet

ADVANCE ARMSLEEVE 910



Compression Level



Features and Benefits



Grip-top Band

Breathable

- Ideal for patients with arm lymphedema or postmastectomy care
- Provides moisture wicking properties to keep skin dry and comfortable
- Breathable and soft

Wicking

- Sensinnov grip-top band ensures armsleeve stays in place all day long
- Available with or without gauntlet

Fiber Content

- 50% Nylon, 28% Rayon, 22% Spandex
- Latex free

SIZING CHART

| | _ | | | | | | | | | | | | | |
|-------------------|--|-------------------------|------------------------|------------------------------|--------------------------------|------------------------------|--|-------------------------|------------------------------|----------------------------------|----------------------------|--------------------------------|------------------------|-------------------------|
| | Sm | nall | | Med | lium | | | Lai | rge | | | X La | arge | |
| Palm (A) | Palm (A) 17-21 cm (6.5"-8") 18-22 cm (7"-8.5") | | | | | 19–23 cm (7.5"–9") | | | 20–24 cm (8"–9.5") | | | | | |
| Wrist (C) | | 6 cm -6.5") | | 16-18 cm (6.5"-7") | | | 18–20 cm (7"–8") | | | | 20–23 cm (8"–9") | | | |
| | Reg | Jular | Regular | | Plus | | Reg | gular Plus | | Regular | | Plus | | |
| Upper Mid-Arm (F) | 26–3 (10"– | 4 cm 13.5") | 29–3 (11.5"– | | 37–43 cm (14.5"–17") | | 32–40 cm (12.5"–15.5") 40–46 c (15.5"–18 | | | 34–42 cm (13.5"–16.5") | | 42–48 cm (16.5"–19") | | |
| | Regular | Long | Regular | Long | Regular | Long | Regular | Long | Regular | Long | Regular | Long | Regular | Long |
| Arm Length (C-G)* | ≤43 cm (17") | > 43 cm (17") | ≤43 cm (17") | > 43 cm (17") | ≤43 cm (17") | > 43 cm (17") | ≤43 cm (17") | > 43 cm (17") | ≤43 cm (17") | > 43 cm (17") | ≤43 cm (17") | > 43 cm (17") | ≤43 cm (17") | > 43 cm (17") |
| Size | SR | SL | MR | ML | MR/P | ML/P | LR | LL | LR/P | LL/P | XR | XL | XR/P | XL/P |

*Arm should have slight bend in arm when measuring length

| STYLES AND COLOURS | |
|---|---|
| Open Sleeve = O | |
| Armsleeve w/Gauntlet (G) Sensinnov grip-top band | |
| Armsleeve w/out Gauntlet (O) Sensinnov grip-top band | |
| Note: all armsleeves come with Grip-top (/S |) |

15–20 mmHg not available with gauntlet or plus sizing

ARMSLEEVE WITHOUT GAUNTLET


Product Code Example: 912ASLO77/S

| Series | mmHg Level | Style | Size | W/out Gauntlet |
|--------|---------------|-------|------|-------------------|
| 91 | 2 | А | SL | 0 |

Small Long in Beige 912ASLG77/S

ADVANCE ARMSLEEVE 910





ARMSLEEVE WITH GAUNTLET

Product Code Example: 912ALRG77/PS

| Colour | Grip-Top | Series | mmHg Level | Style | Size | W/ Gauntlet | Colour | Plus- Size | Grip-Top |
|--------|----------|--------|---------------|-------|------|----------------|--------|---------------|----------|
| 77 | /s | 91 | 2 | А | LR | G | 77 | /P | S |

SIGVARIS GROUP

MARINER STRIPE

Care for Yourself

Why choose **Consumer Health**?

- \cdot Tired, heavy, restless, aching legs
- \cdot Travelers, occupational leg health risks, expecting mothers, family history of Chronic Venous Disease, prevention of swelling and varicose veins
- Simplified sizing

CONSUMER HEALTH PRODUCTS

| Series |
|---------------|
| 146 / 186 |
| 146/T / 186/1 |
| 142 / 182 |
| 151 / 191 |
| 152 / 192 |
| 153 / 193 |
| 143 / 183 |
| 120 |
| 170 |
| 401 |
| 412 |
| 421 |
| 160 |

CONSUMER HEALTH

| Product | Page |
|--------------------------|------|
| Casual Cotton | 54 |
| Traveno Travel Socks | 56 |
| Cushioned Cotton | 58 |
| Sea Island Cotton | 60 |
| All-Season Merino Wool | 62 |
| Linen | 64 |
| Microfiber Shades | 66 |
| Sheer Fashion | 68 |
| Soft Silhouette Leggings | 70 |
| Athletic Recovery Socks | 72 |
| Performance Sleeves | 74 |
| Merino Outdoor Socks | 76 |
| Eversoft Diabetic Socks | 78 |


SIGVARIS GROUP



CASUAL COTTON 146/186



Compression Level



Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- \cdot Subtle ribbed design
- \cdot Flat, low-friction toe seam to prevent irritation

Fiber Content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- \cdot Latex free



STYLES AND COLOURS

Closed Toe = •

Calf (C)

Product Code Example: 146CB10

| Series | Style | Size |
|--------|-------|------|
| 146 | С | В |



15–20 mmHg Closed Toe Calf Size B in Khaki





186 Men's Calf

Style


15–20 mmHg Closed Toe Calf Size C in Black

186CC99

54

CASUAL COTTON 146/186

| 5-7 | 7.5–9.5 | 10–12 |
|-------------------------------|-------------------------------|---------------------------------|
| 18–23 cm (7"–9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| 6.5–8.5 | 9–11 | 11.5–14 |
| 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| А | В | с |



| Colour |
|--------|
| 10 |

CONSUMER HEALTH


SIGVARIS GROUP



TRAVENO TRAVEL SOCKS 146/T - 186/T



Compression Level



Features and Benefits

- Made with premium Supima Cotton for optimum softness and durability
- Subtle ribbed design
- Flat, low-friction toe seam to prevent irritation

Fiber Content

- 66% Supima Cotton, 26% Nylon, 8% Spandex
- Latex free

STYLES AND COLOURS



SIZING CHART

| 6/Т | Women's Shoe | 5–7 | 7.5-9.5 | 10–12 |
|-----|---------------------|----------------------------|----------------------------|------------------------------|
| 14(| Ankle Circumference | 18–23 cm (7"–9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| ٦/۲ | Men's Shoe | 6.5–8.5 | 9–11 | 11.5-14 |
| 186 | Ankle Circumference | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| | Size | А | В | с |

IMPROVE YOUR CIRCULATION IN THE AIR OR ON THE GROUND

Travens

TRAVEL SOCKS



TRAVENO TRAVEL SOCKS 146/T - 186/T







Travel Socks for Women & Men Available in 15–20 mmHg in black

Countertop displays available. Contact your Territory Manager or call **1.800.363.4999** for more information



CUSHIONED COTTON 142/182



Compression Level



Also available in 20–30 mmHg, see page 38

Features and Benefits

- Supima cotton for highest level of comfort and breathability
- $\cdot\,$ Extra-padded sole, heel, and toe to protect feet
- Ideal for sensitive skin and allergies
- Flat, low friction toe seam
- \cdot Diabetic friendly
- Ideal for people who stand for prolonged periods of time

Fiber Content

- 68% Supima Cotton, 25% Nylon, 7% Spandex
- $\cdot\,$ Latex free

SIZING CHART




STYLES AND COLOURS

Closed Toe = •

Calf (C)

Product Code Example: 142CB00

| Series | Style | Size |
|--------|-------|------|
| 142 | С | В |



15–20 mmHg Closed Toe Calf Size B in White

142CB00



15–20 mmHg Closed Toe Calf Size C in Black

182CC99

CUSHIONED COTTON 142/182

| 5-7 | 7.5–9.5 | 10-12 |
|----------------------------|----------------------------|------------------------------|
| 18–23 cm (7"–9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| 6.5-8.5 | 9–11 | 11.5–14 |
| 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| A | В | С |



Colour

00



SEA ISLAND COTTON 151/191



Compression Level



Also available in 20–30 mmHg, see page 26

Features and Benefits

- Knitted with the world's rarest cotton for a soft,
- luxurious feel
- \cdot Ideal for sensitive skin
- \cdot Breathable and comfortable
- Suitable for those who prefer natural fibers
 for everyday wear
- Smooth, flat toe seam

Fiber Content

- 43% Sea Island Cotton, 39% Nylon, 11% Spandex,
 7% Supima Cotton
- Latex free

SIZING CHART Image: Size Women's Shoe Image: Ankle Circumference Men's Shoe Image: Ankle Circumference Ankle Circumference Image: Size Size

STYLES AND COLOURS

Closed Toe = •

Calf (C)

Product Code Example: 151CA10

| Series | Style | Size |
|--------|-------|------|
| 151 | С | А |



15–20 mmHg Closed Toe Calf Size A in Navy

151CA10



15–20 mmHg Closed Toe Calf Size B in Brown

191CB11

60

SEA ISLAND COTTON 151/191

| 5-7 | 7.5-9.5 | 10–12 |
|----------------------------|----------------------------|------------------------------|
| 18–23 cm (7"–9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| 6.5–8.5 | 9–11 | 11.5-14 |
| 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| A | В | С |



Colour 10

CONSUMER HEALTH


SIGVARIS GROUP



ALL-SEASON MERINO WOOL 152/192



Compression Level



Also available in 20–30 mmHg, see page 28

Features and Benefits

- High quality extra-fine Merino Wool makes this sock soft and breathable
- Moisture wicking capabilities and naturally odour resistant
- Thermoregulating for year-round wear and keeps feet dry
- \cdot Smooth, flat toe seam

Fiber Content

- 53% Extra-fine Merino Wool, 38% Nylon, 9% Spandex
- Latex free



STYLES AND COLOURS

Closed Toe = •

Calf (C)

Product Code Example: 192CB11

| Series | Style | Size |
|--------|-------|------|
| 192 | С | В |



15–20 mmHg Closed Toe Calf Size A in Black

152CA99




15–20 mmHg Closed Toe Calf Size C in Navy

192CC10

62

ALL-SEASON MERINO WOOL 152/192

| 5-7 | 7.5-9.5 | 10-12 |
|----------------------------|----------------------------|------------------------------|
| 18-23 cm (7"-9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| 6.5–8.5 | 9–11 | 11.5–14 |
| 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| A | В | С |



| Colour | |
|--------|--|
| 11 | |

SIGVARIS GROUP



LINEN 153/193

Compression Level



Also available in 20–30 mmHg, see page 30

Features and Benefits

- Lightweight and breathable for all day comfort
- Ideal for sensitive skin with hypoallergenic properties
- Strong, stylish and sustainable in four unique colours

Fiber Content

- 48% Cotton, 23% Linen, 22% Nylon and 7% Spandex
- Latex free



STYLES AND COLOURS

Closed Toe = •

Calf (C)

Product Code Example: 153CA207

| Series | Style | Size |
|--------|-------|------|
| 153 | С | А |



15–20 mmHg Closed Toe Calf Size A in Brown

153CA111




15–20 mmHg Closed Toe Calf Size C in Denim

193CC201

64

| 57 | 7.5–9.5 | 10–12 |
|----------------------------|----------------------------|------------------------------|
| 18–23 cm (7"–9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| 6.5-8.5 | 9–11 | 11.5-14 |
| 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| Α | В | С |



| Colour | |
|--------|--|
| 207 | |





MICROFIBER SHADES 143/183



Compression Level



Also available in 20–30 mmHg, see page 34


Features and Benefits

- 100% of microfiber sits next to skin for a soft feel
- Fashionable, durable socks for everyday wear
- Excellent breathability and moisture management
- Glides on and off with ease
- \cdot Smooth, flat toe seam

Fiber Content

- 87% Nylon, 13% Spandex
- Latex free

SIZING CHART



STYLES AND COLOURS

Closed Toe = •

Actual colour may vary

Calf (C)

Calf (C)

* Available for women only ** Available for men only

Product Code Example: 143CB02

| Series | Style | Size | | | |
|--------|-------|------|--|--|--|
| 143 | С | В | | | |



15–20 mmHg Closed Toe Calf Size B in Pink Argyle

143CB45



15–20 mmHg Closed Toe Calf Size C in Black Plaid

183CC42

66

MICROFIBER SHADES 143/183

| 5-7 | 7.5-9.5 | 10-12 |
|----------------------------|----------------------------|-------------------------------|
| 18–23 cm (7"–9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| 6.5-8.5 | 9-11 | 11.5-14 |
| | 5-11 | ب ار <u>-</u> ۱۱ |
| 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5""–11.5") |
| Α | В | С |



Colour 02

SIGVARIS GROUP

PANTYHOSE SIZING CHART



SHEER FASHION 120

for women

Compression Level




Features and Benefits

- Incredible sheer design in a variety of fashionable colours
- Sheer toe and heel area
- · Comfortable stretch and soft fabric, excellent for everyday wear

Fiber Content

- 79% Nylon, 21% Spandex (Content can vary slightly by style)
- Latex free



CALF / THIGH SIZING CHART

180 cm /5'11"

183 cm / 6'0"

| 120 | Women's Shoe | 57 | 7.5-9.5 | 10–12 |
|-----|---------------------|-------------------------|----------------------------|----------------------------|
| | Ankle Circumference | 18–23 cm (7"–9") | 19–25 cm (7.5"–10") | 20–27 cm (8"–10.5") |
| | Size | Α | В | С |

Weight Guidelines For Thigh Style: Size A: <135Lbs; Size B: <180Lbs; Size C: <230Lbs.

| STYLES AND COLOURS | | | | | | |
|--------------------------------|----------------|-------------|-------------|--------------|------------|-------------|
| Closed Toe = ● Open Toe = ● | CHARCOAL 12 | TAUPE 29 | HONEY 33 | GOLDEN 36 | CAFÉ 73 | BLACK 99 |
| Calf (C) | • | •0 | •0 | •0 | • | • |
| Thigh (N) | • | • | • | • | ٠ | • |
| Pantyhose (P) | • | • | • | • | • | • |

Product Code Example: 120CCO36

| Series | Style | Size | Open Toe | Colour | |
|--------|-------|------|----------|--------|--|
| 120 | С | С | 0 | 36 | |



15-20 mmHg Closed Toe Calf Size C in Charcoal

120CC12



Closed Toe Pantyhose Size C in Black

120PC99

SHEER FASHION 120

| | 30 | 35 | 40 | 45 | 50 | 55 | 60 | 65 | 70 | 175 | 80 | 185 | 06 | 95 | 00 | 205 | 210 | 215 | 220 | 225 | 230 |
|----|--------|--------|--------|--------|----------|--------|----------|--------|----------|--------|--------|--------|--------|--------|--------|--------|--------|-------|---------|-----------|---------|
| | 59/130 | 61/135 | 64/140 | 66/145 | 68 / 150 | 70/155 | 73 / 160 | 75/165 | 77 / 170 | 80/175 | 82/180 | 84/185 | 86/190 | 89/195 | 91/200 | 93/205 | 95/210 | 98/3 | 100/220 | 102 / 225 | 105/230 |
| | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | |
| E | | | | | | S | ZE | С | | | | | | | | | | | | | |
| ') | | | | | | | NKL | | | | | | | | | | | | | | |
| | | | | | CI | RCU | MFE | REN | CE | | | | | | | SIZ | ΕE | | | | |
| | | | | | 20 |).5–25 | 5 cm | (8"–1 | O") | | | | | | | | KLE | | | | |
| | | | | | | | | | | | | | | | CIRC | | | NCE | | | |
| | | | | | | | | | | | | | | 2 | 22–27 | cm | (8.5" | -10.5 | ") | | |
| E | | | | | | SI | ZE | D | | | | | | | | | | | | | |
| ') | | | | | | | NKL | | | | | | | | | | | | | | |
| | | | | | CI | RCU | | | CE | | | | | | | SIZ | ΕF | | | | |
| | | | | | |).5–25 | | | | | | | | | | | KLE | | | | |
| | | | | | | | | | | | | | | | CIRC | | | NCE | | | |
| | | | | | | | | | | | | | | | 2–27 | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | |

Product Code Example: 120PE73

| Series | Style | Size | Colour |
|--------|-------|------|--------|
| 120 | Ρ | E | 73 |

CONSUMER HEALTH

SIGVARIS GROUP



SOFT SILHOUETTE LEGGINGS 170

Compression Level



Features and Benefits

- Advanced soft stretch microfiber fabric shapes legs
 and smooths skin
- Foldable, comfort-stretch waist band provides a seamless look for all-day wearing comfort
- Stylish, subtle ribbed design looks chic and pairs perfectly with dresses, blouses, tunics, and shirts

Fiber Content

- 83% Nylon, 17% Spandex
- Latex free

SIZING CHART

| <g lbs<="" th=""><th>+5 / 100</th><th>8/105</th><th>0110</th><th>2/115</th><th>54 / 120</th></g> | +5 / 100 | 8/105 | 0110 | 2/115 | 54 / 120 |
|--|----------|--------|--|--|---|
| | V | V | <u>(</u>) | L() | LU) |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | Kg / Lbs | Kg/Lbs | Kg/Lbs 1 45/100 48/105 | Kg/Lbs Kg/lbs 45/100 48/105 50/110 | Kg/Lbs Kg/Lbs |

Guidelines for Ankle Circumference Size A: 18–23 cm (7"–9"), Size B & C: 20–26 cm (8"–10"), Size D: 22–29 cm (8.5"–11.5")

STYLES AND COLOURS

Open Foot = **O**

Leggings (L)



| Series | Style | Size | | | |
|--------|-------|------|--|--|--|
| 170 | L | В | | | |




15–20 mmHg Legging Size B in Black

170LB99



Non-binding ankle cuff stays comfortably in place for all day wear







CONSUMER HEALTH

SIGVARIS GROUP





Compression Level



Features and Benefits

- Drirelease yarns provide moisture wicking and odour control
- \cdot Comfortable and easy to put on and take off
- Helps reduce lactic acid buildup
- Mild cushioning in sole, heel, and toe area for added comfort

Fiber Content

- 67% drirelease (57% Polyester, 10% Cotton), 26% Nylon,
 7% Spandex
- Latex free

| SIZING CHART | | |
|--------------|---------------------|--|
| | Women's Shoe | |
| 401 | Men's Shoe | |
| | Ankle Circumference | |
| | Size | |

STYLES AND COLOURS

Closed Toe = •

Calf (C)

Product Code Example: 401C200

| Series | Style | Size |
|--------|-------|------|
| 401 | С | 2 |



15–20 mmHg Closed Toe Calf Size X-Large in Black

401CX99



15–20 mmHg Closed Toe Calf Size Medium in White

401CM00

72

ATHLETIC RECOVERY SOCKS 401

| <7.5 | 7.5–9.5 | 10–12 | 12.5–15 | ≥15.5 |
|-----------------------|------------------------|-----------------|--------------------------|--------------------------|
| <6.5 | 6.5–8.5 | 9–11 | 11.5–14 | 14.5–17 |
| 18–24 cm | 19–26 cm | 20–27 cm | 21–29 cm | 21–29 cm |
| (7"–9.5") S | (7.5"–10") M | (8"–10.5") L | (8.5"–11.5") X | (8.5"–11.5") 2 |




Colour



PERFORMANCE **SLEEVES** 412

- Stabilizes and supports active muscles
- Ideal for running, triathlons, endurance sports, and
- Fabric resists deterioration in salt, chlorine,

SIZING CHART



STYLES AND COLOURS

Calf Sleeve = **O**

Calf Sleeve (V)

* Pink and purple not available in X-Large

Product Code Example: 412VL56

| Series | Style | Size | |
|--------|-------|------|--|
| 412 | V | L | |



20-30 mmHg Calf Sleeve Size Small in Purple

412VS58



20–30 mmHg Calf Sleeve Size Medium in Blue

412VM50

74

PERFORMANCE SLEEVES 412

| 18–21.5 cm (7"–8.5") | 21.5–25.5 cm (8.5"–10") | 25.5–29 cm (10"–11.5") | 29–33 cm (11.5"–13") |
|--------------------------------|-----------------------------------|----------------------------------|--------------------------------|
| 28–39 cm | 31.5–45 cm | 35.5–51 cm | 40.5–58 cm |
| (11"–15.5") | (12.5"–17.5") | (14"–20") | (16"–23") |
| S | м | L | x |

| WHITE 00 | BLUE 50 | PINK 56* | PURPLE 58* | BLACK 99 |
|-------------|------------|-------------|---------------|-------------|
| 0 | ο | ο | ο | ο |

Colour 56

CONSUMER HEALTH

SIGVARIS GROUP



MERINO OUTDOOR SOCKS 421

for women & men

Compression Level

15–20 mmHg

Also available in 20–30 mmHg, see page 42

Features and Benefits

- Natural properties of wool make these socks odour resistant
- Thermoregulating and moisture-wicking for yearround wear
- Cushioned, padded soles help prevent blisters
- Mesh flex zones for comfort and ventilation
- Smooth, flat toe seam
- Soft and comfortable
- Perfect for outdoor work environments and activities

Fiber Content

- 67% Fine Merino Wool, 26% Nylon, 7% Spandex
- Latex free



STYLES AND COLOURS

Closed Toe = •

Calf (C)

Product Code Example: 421CS31

| Series | Style | Size |
|--------|-------|------|
| 421 | С | S |



15–20 mmHg Closed Toe Calf Size Large in Olive

421CL31




15–20 mmHg Closed Toe Calf Size Medium in Charcoal

421CM12

76

MERINO OUTDOOR SOCKS 421

| <7.5 | 7.5–9.5 | 10–12 | 12.5–15 |
|------------------------------|-------------------------------|-------------------------------|---------------------------------|
| <6.5 | 6.5–8.5 | 9–11 | 11.5–14 |
| 18–24 cm (7"–9.5") | 19–26 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21–29 cm (8.5"–11.5") |
| S | М | L | X |



Colour 31

CONSUMER HEALTH

SIGVARIS GROUP



EVERSOFT DIABETIC SOCKS 160



Compression Level

8–15 mmHg

Features and Benefits

- Padded sole, heel, and toe for comfort and protection
- $\cdot\,$ Flat, low friction toe seam
- Compression-free foot
- Drirelease yarns provide moisture wicking and odour control
- Suitable for individuals with limited dexterity or lower tolerance for compression
- Ideal for individuals with sensitive skin

Fiber Content

- 67% drirelease (57% Polyester, 10% Cotton), 26% Nylon,
 7% Spandex
- Latex free

SIZING CHART Women's Shoe Men's Shoe Ankle Circumference Size

STYLES AND COLOURS

Closed Toe = •

Calf (C)

Product Code Example: 160CS00

| Series | Style | Size |
|--------|-------|------|
| 160 | С | S |



8-15 mmHg Closed Toe Calf Size Small in White **8–15 mmHg** Closed Toe Calf Size Large in White

160CL00

78

160CS00

EVERSOFT DIABETIC SOCKS 160

| <7.5 | 7.5–9.5 | 10–12 | 12.5–15 |
|------------------------------|-------------------------------|-------------------------------|---------------------------------|
| <6.5 | 6.5-8.5 | 9–11 | 11.5–14 |
| 18–24 cm (7"–9.5") | 19–26 cm (7.5"–10") | 20–27 cm (8"–10.5") | 21-29 cm (8.5"-11.5") |
| S | м | L | x |



| Colour | |
|--------|---|
| 00 |] |





Why Choose Inelastic Wraps?

- Adjustable pressure
- Large variety of solutions

Product

- Compreflex W
- Compreflex St
- Compreflex L
- Compreboot
- Compreflex St
- Accutab
- Compreflex T
- Compreflex K
- Compreknee
- Coolflex
- Compreshort

The SIGVARIS GROUP line of products for medically complex edema (MCE) are for the treatment and management of conditions impacting the lymphatic and venous circulatory systems. The product portfolio contains ready-to-wear garments for the active treatment and maintenance of lymphedema and for patients managing chronic venous disorders in the lower extremities. Sigvaris is looking towards a brighter future with new therapies for complex venous and lymphatic diseases. Our inelastic wraps deliver maximum therapeutic results with increased patient compliance.

MEDICALLY COMPLEX EDEMA & INELASTIC WRAPS

MEDICALLY COMPLEX EDEMA PRODUCTS

| Page |
|------|
| 82 |
| 84 |
| 86 |
| 87 |
| 88 |
| 88 |
| 89 |
| 90 |
| 91 |
| 92 |
| 96 |
| |

COMPREFLEX

Wrap Details

ACCURATE, MEASURED THERAPEUTIC COMPRESSION FOR SELF-DIRECTED TREATMENT AND SELF-MANAGED CARE.

The Compreflex range of products are low profile compression devices designed to deliver moderate to firm compression to the upper and lower extremities. Each Compreflex device is made of a soft conforming inelastic material, offering low to minimal stretch for proper fit and the gradient application of compression. The Compreflex garments are designed with key features that support the patient, including our patented Accutab® Devices which ensure the accurate application of the prescribed dosage of therapeutic compression, an interior stretch panel for easy donning, a contoured calf for limb shaping, and opposing straps that are rounded for improved fit and patient comfort.





Features & Benefits


- High quality inelastic wrap composed of durable breathable inelastic material with moisture wicking capacities
- Contoured design helps the wrap conform to the natural shape of the leg
- Front stretch panel secures the garment in place, while the Velcro straps assist for easy application
- · Ability to customize with different liner and boot options

Conditions

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- · Lymphedema Active and Maintenance
- · Non-specific or post surgical edema

How it works

PROVEN ACCURACY FOR RESTING PRESSURE



Accutab has been rigorously tested according to established guidelines in the International Compression Club Consensus documents. Testing was conducted using PicoPress, which is recognized for its high accuracy and reproducibility

Effective working pressure

Inelastic compression reduces leg edema and improves venous blood flow

· Inelastic compression can increase arterial perfusion in patients with mild to moderate peripheral arterial disease

Resting Resting Pressure Pressure Donnina Laying Down Compreflex Sitting Still

Working pressure

Walking, stair climbing, sit to stand, toe raises or other therapeutic exercises to actuate the calf muscle

Safe, precise resting pressure

- Under the supervision of a healthcare provider, inelastic compression can be used safely for patients with mild to moderate peripheral arterial disease
- Accutab delivers precise dosages for controlled initial resting pressure when donning a Sigvaris Compreflex garment



Accutab[®] Devices (see page 88 for details)











Non-Compressive One pair of Liners included with each leg (see page 86 for details)

Note: • 1401-BKR, 1411-BKR, 1401-BKT, and 1411-BKT includes a Small Regular Boot

• 1402-BKR, 1412-BKR, 1402-BKT, 1412-BKT, 1403-BKR, 1413-BKR, 1403-BKT, and 1413-BKT includes a Medium /Large Regular Boot

• 1404-BKR, 1414-BKR, 1404-BKT, 1414-BKT, 1405-BKR, 1415-BKR, 1405-BKT, and 1415-BKT includes a X-Large / XX-Large Regular Boot



COMPREFLEX

Standard

For those who require a calf and foot garment

CALF & FOOT

20–50 mmHg

• Compreflex Wrap

Non-compressive Basic Liner

Black

1401-BKR

1402-BKR

1403-BKR

1404-BKR

1405-BKR

Beige

1411-BKR

1412-BKR

1413-BKR

1414-BKR

1415-BKR

• Non attached foot piece

PRODUCT CODE

Size

Small

Regular

Medium

Large

X Large

XX Large

Includes:

| | CALI | F SIZING | CHART (| A-C) | | |
|---|------|----------|----------|----------|----------|----------|
| | | Small | Medium | Large | X Large | XX Large |
| - | с | 29–39 cm | 34–44 cm | 40–50 cm | 46–56 cm | 54–64 cm |
| | А | 19–26 cm | 21–30 cm | 26–36 cm | 31–41 cm | 36–46 cm |

Size Tall

Small

Large

X Large

XX Large

Medium

Black

1401-BKT

1402-BKT

1403-BKT

1404-BKT

1405-BKT

G Length Regular: 30-33.9cm Tall: 34–42cm

Foot Size:

I: Measure circumference of ankle bend and heel J: Measure circumference across metatarsal heads

K: Measure medial length from heel to 1st metatarsal head



Beige

1411-BKT

1412-BKT

1413-BKT

1414-BKT

1415-BKT

Compreflex Calf options



Standard CALF 20–50 mmHg


For those who require an alternative foot garment or do not require a foot piece Includes

Compreflex Wrap

Non-compressive Basic Liner

Can be matched with alternative foot pieces (see page 87)

| PRODUCT CODE | | | | |
|--------------|-------------|-------------|--|--|
| Size Regular | Black | Beige | | |
| Small | 1401-NF-BKR | 1411-NF-BKR | | |
| Medium | 1402-NF-BKR | 1412-NF-BKR | | |
| Large | 1403-NF-BKR | 1413-NF-BKR | | |
| X Large | 1404-NF-BKR | 1414-NF-BKR | | |
| XX Large | 1405-NF-BKR | 1415-NF-BKR | | |
| Size Tall | | | | |
| Small | 1401-NF-BKT | 1411-NF-BKT | | |
| Medium | 1402-NF-BKT | 1412-NF-BKT | | |
| Large | 1403-NF-BKT | 1413-NF-BKT | | |
| X Large | 1404-NF-BKT | 1414-NF-BKT | | |
| XX Large | 1405-NF-BKT | 1415-NF-BKT | | |





Transition CALF 20–50 mmHg

For those who require a compressive liner for the foot and ankle

Includes:

- Compreflex Wrap
- Transition Liner, which provides 10-15 mmHg in the foot and ankle

Ability to wear regular shoe over liner

| | Complete CALF |
|-----------------------|-------------------------|
| and the second second | 20–50 mmHg |
| | |

For those who require a compressive liner with compression in the foot and ankle, and graduated compression in the calf

Includes:

- \cdot Compreflex Wrap
- Complete Liner, which provides 10-15mmHg in the foot and ankle, and graduated compression in the calf

Ability to wear regular shoe over liner

| PRODUCT CODE | | | |
|--------------|-------------|-------------|--|
| Size Regular | Black | Beige | |
| Small | 1401-FC-BKR | 1411-FC-BKR | |
| Medium | 1402-FC-BKR | 1412-FC-BKR | |
| Large | 1403-FC-BKR | 1413-FC-BKR | |
| X Large | 1404-FC-BKR | 1414-FC-BKR | |
| XX Large | 1405-FC-BKR | 1415-FC-BKR | |
| Size Tall | | | |
| Small | 1401-FC-BKT | 1411-FC-BKT | |
| Medium | 1402-FC-BKT | 1412-FC-BKT | |
| Large | 1403-FC-BKT | 1413-FC-BKT | |
| X Large | 1404-FC-BKT | 1414-FC-BKT | |
| XX Large | 1405-FC-BKT | 1415-FC-BKT | |

See page 84 for Sizing Chart

| PRODUCT CODE | | | |
|--------------|-------------|-------------|--|
| Size Regular | Black | Beige | |
| Small | 1401-UC-BKR | 1411-UC-BKR | |
| Medium | 1402-UC-BKR | 1412-UC-BKR | |
| Large | 1403-UC-BKR | 1413-UC-BKR | |
| X Large | 1404-UC-BKR | 1414-UC-BKR | |
| XX Large | 1405-UC-BKR | 1415-UC-BKR | |
| Size Tall | | | |
| Small | 1401-UC-BKT | 1411-UC-BKT | |
| Medium | 1402-UC-BKT | 1412-UC-BKT | |
| Large | 1403-UC-BKT | 1413-UC-BKT | |
| X Large | 1404-UC-BKT | 1414-UC-BKT | |
| XX Large | 1405-UC-BKT | 1415-UC-BKT | |

Liner Options



- Non-compressive liner for the protection of the skin
- \cdot Easy to apply knee-high liners
- drirelease moisture wicking and odour control
- \cdot Flat seams for added comfort
- Silver ionized treatment for bacteriostatic properties
- \cdot Available in Heathered Silver colour only

| Size: One size fits all | | |
|--|--|--|
| Pair | | |
| 1700 - NK1 | | |
| Included With Compreflex Standard Calf & Foot Compreflex Standard Calf | | |



Transition Liner

- Mild compression in foot and ankle only
- drirelease moisture wicking and odour control
 Extra thick sole from toe to heel for added comfort
- Silver ionized treatment for bacteriostatic properties

Pair (Black)

1601-CK

1602-CK

1603-CK

1604-CK

1605-CK

Pair (Beige)

1611-CK

1612-CK

1613-CK

1614-CK

1615-CK

 \cdot Available in Black and Beige

PRODUCT CODE

Small

Large

X Large

XX Large

Medium

Size



- Mild compression in foot and ankle with graduated compression in calf
- \cdot Comfortable knee-high liners
- Drirelease moisture wicking and odour control
 Extra thick sole from toe to heel for added
- comfort • Silver ionized treatment for bacteriostatic
- properties
- Available in Black and Beige


| | PRODUCT CODE (Included with Compreflex Complete) | | |
|----------|--|--------------|--|
| Size | Pair (Black) | Pair (Beige) | |
| Small | 1601-CK-F1 | 1611-CK-F1 | |
| Medium | 1602-CK-F1 | 1612-CK-F1 | |
| Large | 1603-CK-F1 | 1613-CK-F1 | |
| X Large | 1604-CK-F1 | 1614-CK-F1 | |
| XX Large | 1605-CK-F1 | 1615-CK-F1 | |



Compreboot Lite FOOT

- Ultra-thin profile that easily fits into most footwear
- \cdot Highly adjustable design may be trimmed to fit
- \cdot Provides mild, therapeutic compression

| PRODUCT CODE (Available in Black Only) | | | | |
|---|----------|--|--|--|
| Size | Black | | | |
| Small / Medium / Large | | | | |
| Regular | 1803-BLR | | | |
| Long | 1803-BLL | | | |
| X Large / XX Large | | | | |
| Regular | 1805-BLR | | | |
| Long | 1805-BLL | | | |



Foot Size:

- I: Measure circumference of ankle bend and heel
- **J:** Measure circumference across metatarsal heads
- K: Measure medial length from heel to \mathbf{l}^{st} metatarsal head

LINER

TRANSITION AND COMPLETE LINERS SIZING CHART

| Size | Ankle Circ. (A) | Calf Circ. (C) | Women's Shoe | Men's Shoe |
|----------|-----------------|----------------|-----------------|------------|
| Small | 19–26 cm | 29–39 cm | <71/2 | <61⁄2 |
| Medium | 21–30 cm | 34–44 cm | 7.5–9.5 | 6.5–8.5 |
| Large | 26–36 cm | 40–50 cm | 10–12 | 9–11 |
| X Large | 31–41 cm | 46–56 cm | 12.5–15 | 11.5–14 |
| XX Large | 36–46 cm | 54–64 cm | 12.5–15 | 11.5–14 |

COMPREBOOT LITE

FOOT SIZING CHART (I, J, K)

| | Small / Medium / Large | | X Large |
|---|------------------------|----------|----------|
| | Regular | Long | Regular |
| I | 28–36 cm | 28–36 cm | ≤44cm |
| J | 22–26 cm | 22–26 cm | ≤33 cm |
| к | 14–18 cm | 19–23 cm | 20–22 cm |

Compreboot Options



- Four-way stretch for firm compression around ankle
- \cdot Ideal for mild to moderate chronic edema
- \cdot Low stretch material conforms to foot shape

| PRODUCT CODE | | | | | |
|--------------------|----------|----------|--|--|--|
| Size | Black | Beige | | | |
| | Small | | | | |
| Regular | 1801-BTR | 1811-BTR | | | |
| Long | 1801-BTL | 1811-BTL | | | |
| Medium / Large | | | | | |
| Regular | 1803-BTR | 1813-BTR | | | |
| Long | 1803-BTL | 1813-BTL | | | |
| X Large / XX Large | | | | | |
| Regular | 1805-BTR | 1815-BTR | | | |
| Long | 1805-BTL | 1815-BTL | | | |





Variable Compression

- Chip foam lining creates high / low pressure ridges
- Removable gripper patches provide traction
- Provides moderate, therapeutic compression

| PRODUCT CODE (Available in Black Only) | | | |
|---|----------|--|--|
| Size | Black | | |
| Small | | | |
| Regular | 1801-BPR | | |
| Long | 1801-BPL | | |
| Medium / Large | | | |
| Regular | 1803-BPR | | |
| Long | 1803-BPL | | |
| X Large / XX Large | | | |
| Regular | 1805-BPR | | |
| Long | 1805-BPL | | |

COMPREBOOT STANDARD AND PLUS

| | FOOT SIZING CHART (I, J, K) | | | | | | | |
|------|-----------------------------|---|----------|----------|----------|----------|-----------|----------|
| ge / | XX Large | | Sm | nall | Medium | / Large | X Large / | XX Large |
| | Long | | Regular | Long | Regular | Long | Regular | Long |
| | ≤44 cm | 1 | 28–36cm | 28–36 cm | <39 cm | <39 cm | <44 cm | <44 cm |
| | ≤33 cm | J | 22–26 cm | 22–26 cm | <30 cm | <30 cm | <33 cm | <33 cm |
| ı | 25–27 cm | к | 14–18 cm | 19–23 cm | 18–20 cm | 23–25 cm | 20–22 cm | 25–27 cm |



Compreflex Strap Extender

Provides up to 10 cm of additional circumferential size range to Compreflex calf wraps

- Detaches easily
- Potential solution for patients outside of standard size ranges of Compreflex calf wraps or needing custom size

| | 1 Unit | 10 Unit Pack |
|-------|----------|--------------|
| Black | 1401-SE1 | 1401-SEB |
| Beige | 1411-SE1 | 1411-SEB |




Accutab Devices

The Accutab accurately delivers the level of compression prescribed by the health care provider. An easy to use device, the Accutab ensures the safe delivery of compression in doses of 20–30mmHg, 30–40mmHg, or 40–50mmHg.





G1 Length Regular: 20–24.9cm | Tall: 25–30cm





Standard Thigh

20–50 mmHg



Features

- The Compreflex Standard Thigh is an inelastic compression device for the thigh. The product includes a Belt and Hip Attachment to support the product throughout daily use.
- Pair the Compreflex Standard Thigh with the Compreflex Standard Knee (sold separately) to create a complete compression system for the upper leg.

Ideal for

- \cdot Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- \cdot Lymphedema Active and Maintenance
- \cdot Non-specific or post surgical edema

Includes

- Compreflex Standard Thigh
- \cdot 1 Pair of non-compressive Cotton Liners per leg
- \cdot Belt and Hip Attachment
- Accutab Devices

Colours



| | COMPREFLEX STANDARD THIGH SIZE CHART (D-F) & PRODUCT CODE | | | | | |
|----------|--|------------|------------|------------|------------|--|
| | | Small | Medium | Large | X Large | |
| | F | 48–58cm | 56–66cm | 64–74cm | 74–84cm | |
| | E | 43–53cm | 51–61cm | 58–68cm | 68–78cm | |
| | D | 38–48cm | 46–56cm | 53–63cm | 63–73cm | |
| BLACK | | | | | | |
| Left | Regular | 1401-TCR-L | 1402-TCR-L | 1403-TCR-L | 1404-TCR-L | |
| Ľ | Tall | 1401-TCT-L | 1402-TCT-L | 1403-TCT-L | 1404-TCT-L | |
| Right | Regular | 1401-TCR-R | 1402-TCR-R | 1403-TCR-R | 1404-TCR-R | |
| , Big | Tall | 1401-TCT-R | 1402-TCT-R | 1403-TCT-R | 1404-TCT-R | |
| | | | BEIGE | | | |
| Left | Regular | 1411-TCR-L | 1412-TCR-L | 1413-TCR-L | 1414-TCR-L | |
| Ľ | Tall | 1411-TCT-L | 1412-TCT-L | 1413-TCT-L | 1414-TCT-L | |
| Right | Regular | 1411-TCR-R | 1412-TCR-R | 1413-TCR-R | 1414-TCR-R | |
| , Big | Tall | 1411-TCT-R | 1412-TCT-R | 1413-TCT-R | 1414-TCT-R | |



COMPREFLEX STANDARD KNEE SIZE CHART (C1-D)

| | Small | Medium | Large | X Large |
|-------|----------|----------|----------|----------|
| D | 38–48 cm | 46–56 cm | 53-63 cm | 58–73 cm |
| C1 | 29–39 cm | 34-44 cm | 39-49 cm | 44–65 cm |
| Black | 1401-KP | 1402-KP | 1403-KP | 1404-KP |
| Beige | 1411-KP | 1412-KP | 1413-KP | 1414-KP |


Standard Knee

20–50 mmHg

Features

• Provides accurate, inelastic compression levels of 20–50 mmHg

- Includes inner donning sleeve for easy application
- Product includes a Butterfly Strap for hard to compress areas



*Arrows indicate reference points of D and C1. Length measurement is not required

COMPREKNEETM

Standard Knee **CLOSED PATELLA**



| COMPREKNEE SIZE CHART (C1-D) | | | | |
|------------------------------|----------|----------|-----------------|--|
| | Small | Medium | Large / X Large | |
| D | 38–48 cm | 46–56 cm | 53-63 cm | |
| C1 | 29–39 cm | 34-44 cm | 39-49 cm | |
| Black | 1101-KP | 1102-KP | 1103-KP | |
| Beige | 1111-КР | 1112-КР | 1113-КР | |


(-)•(-

Variable compression

Features

- Garment covers the knee, providing full coverage
- Straps provide a full-range of adjustment to create a comfortable fit
- Conforming material for a firm secure fit

Colours



*Arrows indicate reference points of D and C1. Length measurement is not required

MEDICALLY COMPLEX EDEMA & INELASTIC WRAPS

COOLFLEXTM

Wrap Details

LIGHTWEIGHT, BREATHABLE COMPRESSION WEAR

Coolflex is a lightweight compression device, optimal for individuals with fragile or broken skin. Made of a non-stretch mesh material, the Coolflex transfers heat and wicks away moisture, leaving the skin cool and dry.

The features of the closure system, along with the built-in guidelines for application, promote independence and encourage improved patient outcomes.

This product is patent pending.

comfort

20–50 mmHg

Ideal for

- Chronic Venous Disease (CVD) including, but not limited to, ulcer due to venous stasis
- Well-managed Lymphedema
- Non-specific or post-surgical edema

Features

- Lightweight and breathable for cool-wearing comfort
- Easy to put on for patients with limited dexterity
- Durable metal clasps, instead of Velcro closure
- · Low profile design allows boot to fit in most footwear
- Spacer Fabric TR is highly breathable for cool-wearing comfort



ONE-HANDED DONNING

EZ-Close Pressure System of eight tabs makes one-handed, lateral donning simple, even for patients with dexterity issues

MOST BREATHABLE GARMENT

The Spacer Fabric[®] TR is designed to be highly breathable for cool-wearing

INNOVATIVE DESIGN

Durable, metal hook closures for maximum strength and durability



ACCURATE COMPRESSION IS EASY WITH THE COOLFLEX

Each Coolflex calf garment includes the unique, EZ-Close™ pressure system. The EZ-Close pressure system is positioned laterally, allowing the patient to apply the product without bending forward. The bands of the closure system are made to an optimal length, making them easy to grip and close using minimal hand strength and little dexterity.



The integrated EZ-Close Pressure System has visual indicators for the accurate delivery of 20–30 mmHg, 30–40 mmHg or 40–50 mmHg on each tab



Standard Calf



Ideal for both venous and lymphatic limb disorders, the garments are available in standard and max sizing and can accommodate limbs that feature disproportionately small ankles with a larger calf.

(-<u>)</u>(()

Features

- Lightweight and breathable for cool-wearing comfort
- Easy to put on for patients with limited dexterity
- Durable metal clasps, instead of velcro closure

Includes

- Coolflex Standard Calf
- 1 Pair of non-compressive basic Liners per leg


Colour





| COOLFLEX STANDARD CALF SIZE CHART (A- | (1) |
|---------------------------------------|-----|
| COOLI LEA STANDARD CALL SIZE CHART (A | |
| (Available in Black Only) | |

| | | Small | Small Max | Medium / Large | Medium / Large Max | X Large / XX Large |
|-------|---------|---------------|---------------|-------------------|-----------------------|-----------------------|
| | С1 | 28–39cm | 35–50cm | 34–50cm | 45–64cm | 45–64cm |
| | с | 28–39cm | 35–50cm | 34–50cm | 45–64cm | 45–64cm |
| | в | 23–33cm | 27–38cm | 29–43cm | 35–49cm | 37–54cm |
| | А | 16–26cm | 16–26cm | 23–33cm | 23–33cm | 30–42cm |
| Left | Regular | 2101-NF-BKR-L | 2101-NF-BMR-L | 2103-NF-BKR-L | 2103-NF-BMR-L | 2105-NF-BKR-L |
| | Tall | 2101-NF-BKT-L | 2101-NF-BMT-L | 2103-NF-BKT-L | 2103-NF-BMT-L | 2105-NF-BKT-L |
| Right | Regular | 2101-NF-BKR-R | 2101-NF-BMR-R | 2103-NF-BKR-R | 2103-NF-BMR-R | 2105-NF-BKR-R |
| А | Tall | 2101-NF-BKT-R | 2101-NF-BMT-R | 2103-NF-BKT-R | 2103-NF-BMT-R | 2105-NF-BKT-R |



COOLFLEX STANDARD FOOT SIZE CHART (I, J, K) (Available in Black Only)

| | Small | | Medium / Large | | X Large / XX Large | |
|-------|------------|------------|----------------|------------|--------------------|------------|
| | Regular | Long | Regular | Long | Regular | Long |
| I | 28–39cm | N/A | 30–43cm | N/A | 34–49cm | N/A |
| J | 22–26cm | N/A | 25–30cm | N/A | 28–33cm | N/A |
| к | 14–18cm | 19–21cm | 18–20cm | 20–22cm | 20–22cm | 23–25cm |
| Right | 2101-BTR-R | 2101-BTL-R | 2103-BTR-R | 2103-BTL-R | 2105-BTR-R | 2105-BTL-R |
| Left | 2101-BTR-L | 2101-BTL-L | 2103-BTR-L | 2103-BTL-L | 2105-BTR-L | 2105-BTL-L |

COOLFLEX

Standard Foot


Variable compression



Foot compression garments are pictured with the Basic Liner. Basic liner is included with the wrap but not with the footpiece

Features

- Low profile design allows boot to fit in most footwear
- Spacer Fabric TR is highly breathable for cool-wearing comfort
- Durable metal clasps, instead of velcro closure



Measuring Instructions Foot Size:

I: Measure circumference of ankle bend and heel J: Measure circumference across metatarsal heads

- K: Measure medial length from heel to 1st
- metatarsal head

COMPRESHORT®



IDEAL FOR PATIENTS WITH ABDOMINAL OR TRUNCAL EDEMA WHO NEED MILD, SUPPORTIVE COMPRESSION

Product available in short and capri.



Features

The Compreshort offers mild, supportive compression for the treatment of abdominal or truncal and upper thigh edema. Available in short and capri lengths, patients with lymphedema and lipedema benefit from the therapeutic compression and support the garment provides.


Ideal for

- Lymphedema
- Lipedema

Colour

Black

| COMPRESHO | ORT SIZING CHAP | PRODUCT CODE | | | |
|--------------|-----------------|--------------|-----------|---------------------|---------------------|
| | Waist | Нір | Midthigh | Product SKU (Short) | Product SKU (Capri) |
| Small | 69–76cm | 84–91cm | 41–46cm | 1301-CS | 1301-CP |
| Small Max | 76–86cm | 91–102cm | 41–46cm | 1301-CM | 1301-CPM |
| Medium | 76–86cm | 91–102cm | 46–53cm | 1302-CS | 1302-CP |
| Medium Max | 86–97cm | 102–112cm | 46–53cm | 1302-CM | 1302-CPM |
| Large | 86–97cm | 102–112cm | 53–61cm | 1303-CS | 1303-CP |
| Large Max | 97–112cm | 112–127cm | 53–61cm | 1303-CM | 1303-CPM |
| X Large | 97–112cm | 112–127cm | 61–71cm | 1304-CS | 1304-CP |
| X Large Max | 112–132cm | 127–147cm | 61–71cm | 1304-CM | 1304-CPM |
| 2X Large | 112–132cm | 127–147cm | 71–81cm | 1305-CS | 1305-CP |
| 2X Large Max | 132–152cm | 147–168cm | 71–81cm | 1305-CM | 1305-CPM |
| 3X Large | 132–152cm | 147–168cm | 81–91cm | 1306-CS | 1306-CP |
| 3X Large Max | 152–172cm | 168–188cm | 81–91cm | 1306-CM | 1306-CPM |
| 4X Large | 152–172cm | 168–188cm | 91–101cm | 1307-CS | 1307-CP |
| 4X Large Max | 172–193cm | 188–208cm | 91–101cm | 1307-CM | 1307-CPM |
| 5X Large | 172–193cm | 188–208cm | 101–112cm | 1308-CS | 1308-CP |



Note: Available in Black Only

Short Length

Capri Length



Help Yourself

- Aid in putting on and removing stockings
- Stocking cleaning to maximize usage

ACCESSORIES

Gloves

Cotton Underl

Sim-Slide

Magnide On/0

Washing Solu

Sigvaris Doff I

Sigvaris Doff

Melany

It Stays

ACCESSORIES

Why choose **Accessories**?

| | Page |
|----------------|------|
| | 100 |
| liner | 100 |
| | 101 |
| Off | 101 |
| ition | 102 |
| N' Donner | 102 |
| N' Donner Cone | 102 |
| | 102 |
| | 102 |
| | |

GLOVES



Rubber Gloves*

SIZE

Medium

X-Large

ASST'D 4X

Small

Large

S/M/L

PACKAGE: 12 PER PACK

Ridged - The practical and protective way to apply and remove your compression stockings.

PRODUCT CODE

591R400S

591R400M

591R400L

591R400X

591R400A



Latex-Free Gloves

PACKAGE: 12 PER PACK

SIZE

Small

Large

X-Large

Medium

Gripping donning gloves recommended for use in a latex-free environment.

PRODUCT CODE

592R400S

592R400M

592R400L

592R400X



Textile Gloves*

These high-quality gloves have a latex surface to provide a good grip when putting on and taking off compression stockings.

| PACKAGE: 12 PER PACK | | | |
|----------------------|--------------|--|--|
| SIZE | PRODUCT CODE | | |
| X-Small | 593R400XS | | |
| Small | 593R400S | | |
| Medium | 593R400M | | |
| Large | 593R400L | | |
| X-Large | 593R400X | | |
| ASST'D 4x XS/S/M | 593R400A | | |

| RUBBER GLOVE AND LATEX-FREE GLOVE SIZING CHART Measure Circumference of the Hand at the Palm | | | |
|--|-----------------------------|-----------------------------|------------------------------|
| 16.5–18 cm (6.5"–7") | 19–20.5 cm (7.5"–8") | 21.5–23 cm (8.5"–9") | 24–25.5 cm (9.5"–10") |
| Small | Medium | Large | X-Large |

| | GLOVE SI Length an | | | | |
|--------|------------------------------|------------|------------|------------|------------|
| Length | 20.5–21 cm | 21.5–22 cm | 22.5–23 cm | 23.5–24 cm | 24.5–25 cm |
| Width | 7–7.5 cm | 7.5–8 cm | 8–8.5 cm | 9–9.5 cm | 10–10.5 cm |
| | X-Small | Small | Medium | Large | X-Large |

LINERS



Cotton Underliner 10–15 mmHg

Can be used as an underlayer to protect fragile skin, secure wound dressings, add 10–15 mmHg of compression, or create an easy glide surface to aid in putting on and taking off compression stockings.

For sizing, please refer to the Natural Rubber Calf sizing on page 49.

* Caution: Rubber Gloves and Textile Gloves contain natural rubber latex that may cause allergic reaction

DONNING AIDS





Scan here to see how to use the Sim-Slide® bit.ly/SIGssyt





Scan here to see how to use the Magnide On/Off® bit.ly/SIGmofyt

| PACKAGE: 2 PER PACK | | |
|---------------------|-----------------|--------------------------|
| SIZE | PRODUCT CODE | NATURAL RUBBER SIZING |
| Small / Short | 588USS | S1/S3 |
| Small / Long | 588USL | S2/S4 |
| Medium / Short | 588UMS | M1/M3 |
| Medium / Long | 588UML | M2/M4 |
| Large-XL / Short | 588ULS | L1/L3 - X1/X3 |
| Large-XL / Long | 588ULL | L2/L4 - X2/X4 |

GLOVES & LINERS / DONNING AIDS

Sim-Slide (For open toe Styles)

| SIZE | PRODUC | T CODE | | | |
|--------------------------------------|------------|--------|---------|------|--|
| Small | 587SLIDESM | | | | |
| Medium | 587SLIDEMD | | | | |
| Large | 587SLIE | DELG | | | |
| X-Large | 587SLIDEXL | | | | |
| Sizing Chart For Sim-Slide | | | | | |
| Women's Shoe Size | ≤4–6 | 6.5–8 | 8.5–9.5 | ≥10 | |
| Men's Shoe Size | ≤3–4.5 | 5–6.5 | 7–8 | ≥8.5 | |
| Size | S | м | L | XL | |
| | | | | | |

Magnide On/Off (For closed toe Styles)

| SIZE | PRODUCT CODE |
|---------|--------------|
| Medium | 587MAGMD |
| Large | 587MAGLG |
| X-Large | 587MAGXL |

| SIZING CHART For Magnide On/Off | | | |
|---|------|------|--------|
| Women's Shoe Size | ≤8.5 | 9-12 | ≥ 12.5 |
| Men's Shoe Size | ≤7.5 | 8-11 | ≥ 11.5 |
| Size | М | L | XL |

* XL suitable for extra-wide feet

DONNING AIDS



SIGVARIS Doff N' Donner

A soft device that allows you to quickly guide a compression sock or stocking up and down your leg or arm, making it simple to put on and remove all styles of graduated compression. It is perfect for both selfdonning and caregiver use. Latex free.



SIGVARIS Doff N' Donner Cone

This loading device makes it easy to roll stockings onto the SIGVARIS Doff N' Donner. The bottom of the base firmly attaches to most clean, flat surfaces via s



MELANY

MELANY from Sigvaris makes both donning and doffing compression garments simple. Its innovative one-piece design with adjustable handles was made especially for users with limited range of motion and finger dexterity.

| SIZE | PRODUCT COD |
|----------|-------------|
| One Size | 583D100 |
| | |



Scan here to see how to use Sigvaris DOFF N' DONNER bit.ly/sigdnd



Scan here to see how to use the MELANY

| SIZE | CALF CIRCUMFERENCE | PRODUCT CODE |
|---------|-----------------------|-----------------|
| X-Small | 25–30 cm | 579SXS |
| Small | 30–38 cm | 579SSM |
| Medium | 38–46 cm | 579SMD |
| Large | 46–56 cm | 579SLG |

WASHING SOLUTION



Gentle washing solution specifically formulated to help maintain the quality of SIGVARIS GROUP compression socks and stockings. Also great for washing delicate garments!

| PACKAGE OPTIONS | | |
|-----------------------------|-----------------|--|
| CONTAINS | PRODUCT CODE | |
| Six 16 oz. Bottles (473ml) | 586W900 | |
| One 16 oz. Bottle (473ml) | 586W950 | |
| Twelve 2 oz. Bottles (60ml) | 586W975 | |

IT STAYS





Roll-on body adhesive

A gentle roll-on body adhesive that washes off with water and may be used to hold compression garments in place against the skin. IT STAYS! is pliable and moves with the body.

| SIZE | PRODUCT CODE |
|--------------------------------|-----------------|
| One 2 oz Bottle (60ml) | 5811300 |
| Twelve 2 oz. Bottles (60ml) | 5811350 |

SYMPTOMS / INDICATIONS

Tired, achy legs

- Occupational or evening edema
- Leg discomfort from long hours of standing
- Prevent varicosities & swelling during preg
- Preventative for legs predisposed to risk
- Reduce swelling during travel
- Prevent DVT during travel
- Heavy, fatigued, aching legs
- Mild edema in lower extremities
- Mild varicosities with minimal edema
- Mild varicosities during pregnancy
- Lymphedema upper extremities
- Post-procedure of small veins
- Prevention or management of DVT

Note: A prescription is recommended for sto the indications of 20 mmHg or higher

- Moderate varicosities with mild to moderate
- Post-procedure of larger veins to maintain
- CEAP C3 (edema), C4 (skin changes without a change wit healed ulceration), C6 (skin changes with
- Management and treatment of DVT or Po
- Superficial phlebitis
- Lymphedema after decongestant therapy
- Orthostatic / postural hypotension
- Post-fracture, post-traumatic edema
- Lymphedema after decongestant therapy
- Severe Post-Thrombotic Syndrome
- Severe varicosities and/or edema
- Severe arterial insufficiency, intermittent
- Uncontrolled congestive heart failure
- Skin sensitivities or allergies
- Advanced neuropathy with sensory loss in
- Diabetes with arterial compromise
- Confinement to bed or non-ambulatory u prescribed by the physician
- Eczema, hypodermatitis, dermatosclerosis

The therapeutic effect of graduate With the proper CEAP classificatio For more severe symptoms, consul

No liability accepted for non-observance of contraindications and cautions. Fronek HS, Bergan JJ, et al. The Fundamentals of Phlebology: Venous Disease for Clinicians. 2004. pg 151. Blaettner, et al. Leg Symptoms of Healthy People and Their Treatment with Compression Hosiery. Phlebology 2008. Compression therapy of the extremities. Partsch, Rabe, Stemmer. ISBN: 2.85480.770.7.

| la suction. | |
|-------------|--------------|
| SIZE | PRODUCT CODE |
| One Size | 583D200 |
| a wara | |

goto.sigvaris.com/bbiMdz

| ion. | |
|------|--------------|
| | PRODUCT CODE |
| Sizo | 5830200 |

DONNING AIDS & PRODUCT CARE / SYMPTOMS & INDICATIONS

| | COMPRESSION LEVEL AT THE ANKLE |
|---|-----------------------------------|
| ng or sitting gnancy | 15–20 mmHg |
| ockings with | 20–30 mmHg |
| ate edema (also during pregnancy) n treatment success out ulceration), C5 (skin changes with active ulcerations) ost-Thrombotic Syndrome y to maintain reduction | 30–40 mmHg |
| y to maintain reduction | 40 mmHg or higher |
| claudication, ischemia | Contraindications |
| n leg or foot use unless otherwise s, dermatitis | Cautions |
| d compression stockings is bas n, the correct compression leve It with a vascular specialist. | |

GLOSSARY

| | Term | Definition | Sigvaris Usage |
|--|---|--|---|
| | CHRONIC VENOUS DISEASE <i>(CVD)</i> | A disease state comprised of a wide range of disorders including tired, achy legs, varicose veins, Chronic Venous Disease (CVD), skin issues, and venous leg ulcers (open and healed). | SIGVARIS GROUP manufactures MEDICAL and CONSUMER HEALTH CARE products to help reduce the risk of, improve symptoms of, and provide treatment for CVD. |
| Sensitive Skin Thermoregulating Breathable Moisture Wicking Soft Durable | DOUBLE-COVERED YARN | Spandex fibers are covered twice with yarn to prevent bare spandex from touching the skin for improved wearing comfort. | Used in nearly all SIGVARIS GROUP graduated compression socks and hosiery to create products that are easy to put on and take off, while also providing greater comfort. |
| | DEEP VEIN THROMBOSIS <i>(DVT)</i> | A blood clot in a deep vein, most commonly in the lower leg or thigh. | SIGVARIS GROUP cares about DVT prevention. Clinical evidence supports the use of graduated compression to help prevent DVTs and the more serious complication of PE (Pulmonary Embolism), which can be life-threatening. |
| Cotton Wool Microfiber Natural Rubber Largest Size Range Strong and Durable | GRADUATED COMPRESSION VS. ANTI-EMBOLISM | Graduated compression is for ambulatory patients, whereas anti- embolism stockings are for non- ambulatory patients. | SIGVARIS GROUP offers true graduated compression, meaning products are tightest at the ankle and decrease in pressure going up the leg to help improve circulation. |
| | HYPOALLERGENIC | Designed to minimize the chance of an allergic reaction. | SIGVARIS GROUP offers Latex free products and cotton fabrics for those with sensitive skin and allergies. Also, many SIGVARIS GROUP accessories are designed for consumers with sensitive skin. |
| Sheer Medium Sheer Opaque Stylish Stylish Elegant | ISO CERTIFIED | SIGVARIS GROUP is ISO (International Organization for Standardization) certified, indicating we have the required quality systems in place. | SIGVARIS GROUP follows ISO Standard 13485:2016 for our products manufactured in North and South America. |
| | MICROFIBER | An extremely fine synthetic fiber knitted into textiles with the texture and drape of natural-fiber fabric. This fabric has excellent breathability. | Used to create products with a soft and comfortable feel. Find this fiber in 840 Soft Opaque, 830, 143 & 183 Microfiber Patterns, and 820 Microfiber. |
| Colours Patterns Sensinnov Grip-top Anatomically Cushioned Sole Anti-bacterial | mmHg (MILLIMETERS OF MERCURY) | Defined as the extra pressure generated by a column of mercury one millimeter high; the same unit of measurement used for blood pressure. | SIGVARIS GROUP products are offered in several ranges of compression levels based on this measurement. |
| | MERINO WOOL | Merino is the softest and finest wool available in the world and is inherently itch-free. | SIGVARIS GROUP uses this premium wool yarn to create soft products that are thermoregulating for use in all seasons. Used in 422 Thermoregulating Wool, 421 MERINO OUTDOOR SOCKS, 240 Merino Wool, and 152 & 192 ALL-SEASON MERINO WOOL. |
| Hypoallergenic Moisture Natural fibers Management | SEA ISLAND COTTON | Sea Island Cotton is the rarest cotton in the world and represents the finest in cotton softness. | SIGVARIS GROUP uses the highest-quality yarns in the industry to create products that are best in class. Look for this fine quality cotton in our 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS. |
| | SENSINNOV | A patented SIGVARIS silicone top band that gently grips the skin and keeps the band in place all day. | This top band is featured in 780 Sheer, 840 Soft Opaque, 710 Patterns, and 820 Microfiber to help keep thigh-high stockings in place all day long. |
| | SUPIMA COTTON | Supima accounts for only about three percent of annual cotton production in the United States. Its fineness and longer staple length make it a premium cotton fiber. | SIGVARIS GROUP uses this premium cotton yarn in the 230 Cotton, 220 Sea Island and 151 & 191 SEA ISLAND COTTON SOCKS, 146 & 186 CASUAL COTTON, and 360, 142 & 182 Cushioned Cotton products, making them comfortable and ideal for those with sensitive skin. |

Please follow these instructions carefully to determine the correct size for different product styles.

Measurements should always be taken in the following order, and as early in the morning as possible when swelling is minimal.

CALF STYLE

Open/Closed Toe

Open/Closed Toe

your first measurement.

1

6

| 1 | Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to determine if ankle is SMALL, MEDIUM, LARGE or X-LARGE. Note: This should be your first measurement. |
|---|--|
| 2 | Measure the calf circumference at the fullest part of the calf. |
| 3 | Determine the calf length by measuring from the floor up to the bend of the knee. This measurement is taken on the outside of the leg. |

Measure the ankle circumference at the narrowest part of the ankle, above the ankle bone. Use the appropriate series sizing chart to

2 Measure the calf circumference at the fullest part of the calf.

5 Determine **the leg length** by measuring from the floor up to

Determine **the hip circumference** by measuring the

3 Skip Step 3 (but refer to General Guidelines #6 below).


4 Measure the thigh at its widest circumference.

determine if ankle is SMALL, MEDIUM or LARGE. Note: This should be

PANTYHOSE AND THIGH-HIGH STYLE



Ankle circumference directly above ankle bone







widest part of the hips.

the gluteal fold.

| 1 | Measurements should always be taken as early in the morning as possible, when swelling is minimal. |
|---|--|
| 2 | Measure against bare skin, not over clothing (Steps 1–5). |
| 3 | The measuring system is based on having the foot flat against the floor and the calf or leg at a right angle to the foot (see images to right for examples). |
| 4 | Measurements should be taken by pulling the tape around the leg so that it is snug but not constricting. |
| 5 | Length measurements should be taken from the ground up. |
| 6 | Take all measurements, even if the patient only requires calf-length products, as the patient might decide at a later date to purchase additional styles. |



Calf length from the floor to the right-angle bend of the knee

3

5 Leg length from the floor up to the gluteal fold



2 Calf circumference at fullest

Widest circumference of the thigh



6 Circumference at the hip

| Series Number | Compression Level | Product Style | Size | Gender | Colour |
|---|--|---|--|--------------------------------------|--|
| 23 = Cotton See individual pages for series number) | 1 = 15–20 mmHg 2 = 20–30 mmHg 3 = 30–40 mmHg 4 = 40–50 mmHg* 5 = 50–60 mmHg* | C = Knee-high or Calf (See individual pages for available styles) | ML = Medium Long (See individual pages for sizing) | M = Men W = Women O = Open toe | 99 = Black (See individual pages for available colours, |
| 23 | 2 | С | ML | М | 99 |

* 40–50 mmHg and 50–60 mmHg available in Natural Rubber only

To order Plus Calf with Grip-top or Plus Size Pantyhose (860P for women only): add "/P" after the colour in the ordering code

To order Calf Style with Grip-top: add "/S" after the colour in the ordering code (see product pages for availability)

To order Thigh with Waist Attachment:

To order Motion High Tech and Motion Thermoregulating Wool: The letter for Gender is not used in the ordering code.

| OW TO ORDER CONSUMER HEALTH PRODUCTS | | | | |
|--|--|---|---|--|
| Series Number | Product Style | Size | Colour | |
| 120 = SHEER FASHION (See individual pages for series number) | N = Thigh-high (See individual pages for available styles) | A = Small (See individual pages for sizing) | 33 = Honey (See individual pages for available colours) | |
| 120 | Ν | Α | 33 | |

and EVERSOFT DIABETIC SOCKS:

Size is only one letter and the letter for Gender is not used in the ordering code. Example: 421CL31 = MERINO OUTDOOR SOCK, 15–20 mmHg, Closed toe Calf, Size Large, Colour Olive

in many countries worldwide. © Copyright 2021 SIGVARIS, Inc.

add "/L" to indicate the Left leg or "/R" to indicate the Right leg after the colour in the ordering code (860W and 500W only)

Example: 412CSM59 = Motion High Tech, 20–30 mmHg, Closed toe Calf, Size Small Medium, Colour Red

To order ATHLETIC RECOVERY SOCKS, PERFORMANCE SLEEVES, MERINO OUTDOOR SOCKS,

Sigvaris, SIGVARIS GROUP, PFS², Sensinnov, Compreflex, Compreboot, and Accutab are registered trademarks of SIGVARIS AG, CH-9014 St.Gallen/Switzerland,

Coolflex, Spacer Fabric TR, and CoverOne are trademarks of SIGVARIS AG, CH-9014, St. Gallen/Switzerland, in many countries worldwide.

Lycra is a registered trademark of INVISTA. Supima is a registered trademark of the Supima Association of America Corporation. NILIT is a registered trademark of Nilit. All rights reserved. dri release is copyright ©2021 dri release. dri release is an Optimer Brands company. All rights reserved. Sim-Slide and Magnide are registered trademarks of Arion Holding B.V.

Scan to see digital catalog



To Order :

1.800.363.4999 Fax: 1.800.263.8736 canadaorders@sigvaris.com **Online Dealer Portal:** sigvaris-online.com







(in)



